

Product datasheet for **SC336825**

CRMP1 (NM_001288661) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CRMP1 (NM_001288661) Human Untagged Clone
Tag:	Tag Free
Symbol:	CRMP1
Synonyms:	CRMP-1; DPYSL1; DRP-1; DRP1; ULIP-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336825 representing NM_001288661.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGAAAATGAACAACGAAGTGGTGGACAAGAGTGACCGACTCCTCATCAAAGGTGGACGGATCATC
AACGATGACCAATCCCTTTATGCTGACGCTACCTGGAGGATGGACTTATCAAACAAATAGGAGAGAAC
TTAATCGTTCCTGGTGGAGTGAAGACCATTGAAGCCAACGGGCGGATGGTTATCCCGGAGGTATTGAT
GTCAACACGTACCTGCAGAAGCCCTCCCAGGGGATGACTGCGGCTGATGACTTCTTCCAAGGACCAGG
GCGGCACTGGTGGGCGGACCACGATGATCATTGACCATGTTGTTCTGAACCTGGGTCCAGCCTACTG
ACCTCTTTCGAGAAGTGGCACGAAGCAGCTGACACCAAATCCTGCTGTGATTACTCCCTCCACGTGGAC
ATCACAAGCTGGTACGATGGCGTTCGGGAGGAGCTGGAGGTGCTGGTGCAGGACAAAGGCGTCAATTCC
TTCCAAGTCTACATGGCCTATAAGGATGTCTACCAAATGTCCGACAGCCAGCTCTATGAAGCCTTTACC
TTCCCTAAGGGCCTGGGAGCTGTGATCTTGGTCCATGCAGAAAATGGAGATTTGATAGCTCAGGAACAA
AAGCGGATCCTGGAGATGGGCATCACGGTCCCGAGGGCCATGCCCTGAGCAGACCTGAAGAGCTGGAG
GCCGAGGCGGTGTTCCGGGCCATCACCATTGCGGGCCGGATCAACTGCCCTGTGTACATACCAAGGTC
ATGAGCAAGAGTGCAGCCGACATCATCGCTCTGGCCAGGAAGAAAGGGCCCTAGTTTTTGGAGAGCCC
ATTGCCGCCAGCCTGGGACCGATGGCACCATTACTGGAGCAAGAACTGGGCCAAGGCTGCGGCGTTTC
GTGACTTCCCCTCCCCTGAGCCCGGACCCTACCACGCCGACTACTTGACCTCCCCTACTGGCCTGTGGG
GACTTGCAGGTCACAGGCAGCGGCCACTGTCCCTACAGCACTGCCAGAAGGCGGTGGGCAAGGACAAC
TTTACCCTGATCCCCGAGGGTGTCAACGGGATAGAGGAGCGGATGACGGTCGTCTGGGACAAGGCGGTG
GCTACTGGCAAAATGGATGAGAACCAGTTTGTGCTGTCACCAGCACCAATGCAGCCAAGATCTTTAAC
CTGTACCCAAGGAAAGGCGGATTGCCGTGGGCTCGGATGCCGACGTGGTCACTGGGACCCCGACAAAG
TTGAAGACCATAACAGCCAAAAGTCACAAGTCGGCGGTGGAGTACAACATCTTCGAGGGTATGGAGTGC
CACGGCTCCCCTAGTGGTTCATCAGCCAGGGCAAGATCGTCTTTGAAGACGAAACATCAACGTCAAC
AAGGGCATGGGCCGCTTATTCCGCGGAAGGCGTTCGCGGAGCACCTGTACCAGCGCGTCAAAATCAGG
AATAAGGTTTTTGGATTGCAAGGGGTTCCAGGGGCATGTATGACGGTCTGTGTACGAGGTACCAGCT
ACACCCAAATATGCAACTCCCGCTCCTCAGCCAAATCTTCGCTTCTAAACACCAGCCCCACCCATC
AGAAACCTCCACCAGTCCAACCTCAGCTTATCAGGTGCCAGATAGATGACAACAATCCCAGGCGCACC
GGCCACCGCATCGTGGCGCCCCCTGGTGGCCGCTCCAACATCACCAGCCTCGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001288661

Insert Size: 1713 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288661.1](#)

RefSeq Size: 3278 bp

RefSeq ORF: 1713 bp

Locus ID: 1400

UniProt ID: [Q14194](#)

Cytogenetics: 4p16.2

MW: 62 kDa

Gene Summary:

This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) contains an alternate 5' exon, compared to variant 1. It differs in the 5' UTR and initiates translation at an alternate downstream AUG resulting in a protein (isoform 3) with a shorter N-terminus, compared to isoform 1.