

## Product datasheet for **SC336827**

### **BTBD15 (ZBTB44) (NM\_001301098) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BTBD15 (ZBTB44) (NM_001301098) Human Untagged Clone
Tag:	Tag Free
Symbol:	BTBD15
Synonyms:	BTBD15; HSPC063; ZNF851
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336827 representing NM\_001301098.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGTGTGAAAAATTTACTCATAGCTCCTCTCCACAGCCAGGAAATGCTTGAAAGCTAAATATG
CTGCGAAATGATGGACATTTTGTGATACACTATTCGTGTCCAGGACAAAATCTCCGGGCACATAAG
GTGGTACTAGCAGCTTGCACTGATTTCTTCGCACCAAACCTTGTAGGCCAAGCCGAGGATGAGAACAAG
AATGTGTTGGATCTGCATCATGTTACAGTGACTGGCTTTATACCTCTTTTAGAATATGCTTACACAGCC
ACTCTATCAATTAACACAGAAAATATTATTGATGTTCTAGCAGCAGCCAGCTATATGCAAAATGTTCACT
GTTGCCAGCACCTGCTCAGAGTTCATGAAATCAAGCATTTTATGGAATACACCAACAGCCAACCTGAA
AAGGGTCTAGATGCTGGACAAGAAAATAATTCTAACTGCAATTTTACTTCTCGAGATGGGAGCATTTCT
CCCGTGTCTCAGAGTGCAGTGTGGTAGAAGAACCATTCTGTCTGCCGAGAATCCCGGAGAAAGCGC
AAAAGCTACATTGTTATGTCTCCTGAAAGTCTGTAAAGTGTGGCACACAACAAGCTCACCCAGGTA
TTGAATTTCTCAGCTTCTACTCAGAAAATAGAAACCAACCAGTTGACTTCTCTTAGCTTTTCTTGG
ACTTTTCTTTTGGAAATTGATCGAAGGATTCAGCCTGAGAAAAGTTAAGCAAGCAGAAAATACCCGGACT
TTAGAATTACCTGGCCATCTGAGACCGGTAGAAGAATGGCTGATTATGTGACTTGTGAGAGCACAAAA
ACTACCTTGCCCTTAGGTACCGAAGAAGATGTCGGGTCAAAGTAGAAAGATTAAGTGTGAGGAGGTC
CATGAGGAAGTGTCCAGCCTGTGAGTGCATCTCAGAGTTCGCTGAGTGTACAGCAGACAGTTCCAGGA
AGTGAACAAGTCCAAGAGGACCTTCTGATTAGTCCACAGTCTTCTCTATAGGCTCAGTAGATGAAGGC
GTTTCTGAGGGCTTGCCACACTTCAAAGCACGTCTAGCACTAATGCTCCTCCGGATGATGATGATCGA
TTGGAAAATGTTCAAGTATCCCTACCAACTCTACATTGCTCCTCCACCAGCAGTACAGAGCGACCAAGT
CCAAATGGTCCCAGACAGACCTTTTCAAGTGTCCAACCTGCCGGGTGCGATTACCCGTATTCAGAACCTA
AAGCAGCACATGCTCATCCACTCAGGAATTAACCACTTTCAAGTGTGACCGCTGTGGGAAAAAGTTACC
AGGGCTTACTCGCTAAAGATGCATCGCCTAAAGCATGAAGGTAACGCTGTTTCCGGTGCCAGATATGT
AGTGCCACTTTCCTTTCGCGGAATATAAACACCACATGAGGGTTTCCCGGCACATTATCCGCAAG
CCTCGGATTTACGAGTGCAAAACATGTGGCGCCATGTTCCACCAACTCTGGAAATTTAATCGTGCACCTG
AGGAGTCTGAACCATGAAGCATCAGAGCTAGCAAACACTTCCAGAGCAGTGATTTCTAGTACCCGGAC
TACTTAAACCAGGAACAAGAAGAGACCCTGTTCATATGATCTTGGAGAACACGGTTTTGAAAGCAAC
TCCTCTGTTCAAATGCCTGTCAATTCACAGTACCCTCAAAGGCAAGGAACCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM\_001301098

**Insert Size:** 1713 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001301098.1](#)

**RefSeq Size:** 7203 bp

**RefSeq ORF:** 1713 bp

**Locus ID:** 29068

**UniProt ID:** [Q8NCP5](#)

**Cytogenetics:** 11q24.3

**MW:** 63.9 kDa

**Gene Summary:** May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.