

Product datasheet for **SC336836**

SLC44A3 (NM_001301079) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC44A3 (NM_001301079) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC44A3
Synonyms:	CTL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336836 representing NM_001301079.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCTACTCGGTGGTGGCTGGAGCCGCGGGAAGACTCCTCTTTGGCTATGACAGCTTTGGCAACATG
TGTGGCAAGAAGAACTCCCCCGTGAAGCTGGAAGTCAAAGGGTCTTCTCAGGGCAGGACATGACCCTAAAAAACAC
GTGTTCTTTATGAATTCCTGCAACTGGAAGTCAAAGGGTCTTCTCAGGGCAGGACATGACCCTAAAAAACAC
TTCAACTATACCCACAGTCCAAAAGCAGACTCACTGTGTCCCAGGCTACCAGTTCCTCCAAGCAAGTCA
TTTCCCTATTTAACCGATGTGTCCCTCAAACACCTGAGTGCTACTCCCTATTTGCATCTGTTTTGATA
AATGATGTTGACACCCTCCACCGAATTCTAAGTGAATCATGTCGGAAGAGATAACAATCCTTGGCCTG
TGTATCCTCGCATTAGCCTTGTCTTTGGCCATGATGTTTACCTCAGATTCATCACCACCCTTCTGGTT
CACATTTTCATTTCAATTGGTTATTTTGGGATTGTTGTTTGTCTGCGGTGTTTTATGGTGGCTGATTAT
GACTATACCAACGACCTCAGCATAGAATTGGACACAGAAAGGGAAAATATGAAGTGCCTGCTGGGTTT
GCTATCGTATCCACAGGCATCACGGCAGTCTGCTCGTCTTGATTTTTGTCTCAGAAAGAGAATAAAA
TTGACAGTTGAGCTTTTCCAAATCACAATAAAGCCATCAGCAGTCTCCCTTCTGCTGTTCCAGCCA
CTGTGGACATTTGCCATCCTCATTTTTCTTCTGGGTCTCTGGGTGGCTGTGCTGCTGAGCCTGGGAAC
GCAGGAGCTGCCAGGTTATGGAAGCGGCCAAGTGAATATAAGCCCTTTTCGGGCATTCGGTACATG
TGGTTCGTACCATTTAATTGGCCTCATCTGGACTAGTGAATTCATCCTTGCCTGCCAGCAAATGACTATA
GCTGGGGCAGTGGTTACTTGTATTTCAACAGAAGTAAAAATGATCCTCCTGATCATCCCATCCTTTCG
TCTCTCCATTCTCTTCTTACCATCAAGGAACCGTTGTGAAAGGGTCATTTTAATCTCTGTGGTG
AGGATCCGAGAATCATTGTGATGTACATGCAAAACGCACTGAAAGAACAGCAGCATGGTGCATTGTCC
AAGTACCTGTTCCGATGCTGCTACTGCTGTTTCTGGTGTCTTGACAAATACCTGCTCCATCTCAACCAG
AATGCATATACTACAACCTGCTATTAATGGGACAGATTTCTGTACATCAGCAAAAGATGCATTCAAAATC
TTGTCCAAGAAGTCAAGTCACTTTACATCTATTAAGTCTTTGGAGACTTCATAATTTTTCTAGGAAAG
GTGTTAGTGGTGTGTTTCACTGTTTTTGGAGGACTCATGGCTTTTAACTACAATCGGGCATTCCAGGTG
TGGGCAGTCCCTCTGTTATTGGTAGCTTTTTTGCCTACTTAGTAGCCATAGTTTTTTATCTGTGTTT
GAAACTGTGCTGGATGCACTTTTCTGTGTTTGTGTTGATCTGGAAACAAATGATGGATCGTCAGAA
AAGCCCTACTTTATGGATCAAGAATTTCTGAGTTTCGTAAGGAGCAACAAATTAACAATGCAAGG
GCACAGCAGGACAAGCACTCATTAAAGGAATGAGGAGGGAACAGAAGTCCAGGCCATTGTGAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001301079

Insert Size: 1722 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001301079.1](#)

RefSeq Size: 2236 bp

RefSeq ORF: 1722 bp

Locus ID: 126969

UniProt ID: [Q8N4M1](#)

Cytogenetics: 1p21.3

Protein Families: Transmembrane

MW: 64.5 kDa