

## Product datasheet for **SC336841**

### HCE (RNGTT) (NM\_001286426) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	HCE (RNGTT) (NM_001286426) Human Untagged Clone
Tag:	Tag Free
Symbol:	RNGTT
Synonyms:	CAP1A; hCAP; HCE; HCE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC336841 representing NM\_001286426.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCTCACACAAGATCCCGCCGGTGGCTGAACTGTCCCGCGCGGCCAGCCGGTGGCAGGAAGA
TTCTTACCTCTGAAGACAATGTTAGGACCAAGATATGATAGTCAAGTTGCTGAAGAAAATCGTTCCAT
CCAGCATGCTCTCAAATTACCTAAAGAGCCTAAAGGTTAAAATGGGCTTGTGGTGGACCTGACAAT
ACTTCAAGTTTCTATGACCGAAAATGACATAGAAAAAGAAGGAATCAAATATATAAAAATTCAAGTAAA
GGACATGGTGAAGTACCCTACCACTGAGAATACTGAGACCTTTATTCGTCTGTGTGAGCGGTTAATGAA
AGAAAATCCACCTGAACCTATAGGTGTTCACTGACTCATGGCTTCAATCGCACTGGTTTCTCATATGT
GCCTTTTTGGTGGAGAAAATGGATTGGAGTATCGAAGCAGCAGTTGCTACTTTTGCCCAAGCCAGACCA
CCAGGAATCTACAAGGTGATTATTTGAAGGAACTTTTTCGTGGTATGGTGACATAGAGGAAGCACCA
CCCCACCTCTATTGCCAGATTGGTGTGGTGGAGGATGATGAAGACGAAGATGAGGATGAGGATGAAAG
AAGGAATCAGAACCCGGGTCAAGTGCTTCTTTGGCAAAAGGAGAAAAGAACGGTTAAAACCTGGGCGCT
ATTTTCTTGAAGGTGTTACTGTTAAAGGTGTAACCTCAAGTAACAACACAACCAAAGTTAGGAGAGGTA
CAGCAGAAGTGTCAATTCTGTGGCTGGGAAGGGTCTGGATTCCCTGGAGCACAGCCTGTTTCCATG
GACAAGCAAAATATTAACCTTTAGACCTGAAGCCATACAAAAGTAAGCTGGAAAGCAGATGGTACTCGG
TACATGATGTTGATTGATGGCACAAATGAAGTTTTTATGATTGATAGAGACAATTCAGTATTTTCATGTT
TCAAATCTGGAATTTCCATTTTCGTAAGATCTTCGTATGCATTTATCAAATACTCTCTTGGATGGCGAG
ATGATTATTGACAGAGTAAATGGACAGGCTGTTCCCTAGATATTTGATATATGACATAATTAATTCAT
TCACAGCCCGTTGGAGATTGTGATTTTAAATGTTCTGCTGCAGTGTATAGAACGAGAAATATAAGTCCT
CGACACGAAAAAATGAAGACTGGGCTCATTGACAAAACACAGGAACCATTTAGCGTCAGAAATAAGCCG
TTTTTTGACATCTGTACTTCAAGAAAGAAATACAAACCTGGTCGATGTGATGATATTTTGAATGGAAG
CCTCCAGTCTGAATTTCTGTGGATTTTCGTCTAAAAATAACAAGAATGGGAGGAGAAGGGTTACTTCT
CAGAATGTTGGCCTCCTGTATGTTGGAGGTTATGAAAGACCTTTGCACAAATCAAGGTGACAAAAGAG
CTGAAACAGTATGACAACAAAATATAGAATGCAAATTTGAGAACAACAGCTGGGTCTTCATGAGACAG
AGAACAGACAAAAGTTTTCTAATGCCTACAACACTGCCATGGCTGTGTGTAACAGCATCTCAAACCT
GTCACCAAGGAGATGCTGTTGAGTTCATCGACAGATGTACTGCAGCTTCTCAAGGACAGAAGCGAAAA
CATCATCTGGACCCTGACACGGAGCTCATGCCACCACCCTCCAAAAGACCACGCCCTTAACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001286426

**Insert Size:** 1725 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286426.1](#)

**RefSeq Size:** 4772 bp

**RefSeq ORF:** 1725 bp

**Locus ID:** 8732

**UniProt ID:** [O60942](#)

**Cytogenetics:** 6q15

**Protein Families:** Druggable Genome, Phosphatase

**MW:** 66 kDa

**Gene Summary:**

Bifunctional mRNA-capping enzyme exhibiting RNA 5'-triphosphatase activity in the N-terminal part and mRNA guanylyltransferase activity in the C-terminal part. Catalyzes the first two steps of cap formation: by removing the gamma-phosphate from the 5'-triphosphate end of nascent mRNA to yield a diphosphate end, and by transferring the gmp moiety of GTP to the 5'-diphosphate terminus.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 3' coding region, compared to variant 1, resulting in an isoform (b) that is shorter than isoform a.