

Product datasheet for **SC336848**

CDH18 (NM_001291957) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDH18 (NM_001291957) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDH18
Synonyms:	CDH14; CDH14L; CDH24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336848 representing NM_001291957.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAAAACTACTAGCACATCTTGCATCTGTCCAGTCTAGTGTGTCTCTGTTTTGTGCAGAGGTGTTAT
GGAACGTCTCACCACAGCTCCATCAAGGTGATGAGAAACCAACCAACACATTGAAGGTGAAACCGAA
GTCCATCATCGTCCAAAAGGGGATGGGTATGGAATCAGTTCTTTGTTTTAGAAGAACATATGGGACCA
GATCCTCAGTATGTTGAAAGCTGCACTCCAATTCTGACAAAGGTGATGGATCTGTCAAGTACATCCTT
ACTGGAGAGGGTCTGGGACTATATTTATCATTGACGATACCACGGGTGATATCCACTCAACAAAAAGC
CTAGACAGAGAGCAGAAGACCCACTATGTGCTTCATGCTCAAGCTATTGATAGACGTACAAACAAACCT
CTTGAGCCTGAATCCGAGTTCATCATCAAAGTGAAGACATCAATGACAACGCTCCAAAATTCACAGAT
GGACCATACATTGTTACTGTGCCTGAAATGTCAGATATGGGTACCTCTGTTCTACAGGTGACAGCTACT
GATGCAGATGACCCTACCTATGAAACAGCGCTCGGGTGGTTACAGCATTCTCCAAGGACAACCCTAC
TTCTCCGTCGACCCTAAAACAGGAGTTATTAGAACGGCCTTACATAACATGGACAGAGAAGCCAGAGAA
CATTACTCCGTAGTCATTCAAGCCAAAGACATGGCTGGGCAAGTTGGAGGGCTTTTCAGGATCTACAACA
GTCAACATCACCTTAACCGATGTCAATGACAACCCACCACGCTTTCCTCAAAAACACTATCAGCTATAT
GTTCTGAGTCAGCTCAAGTTGGTTGAGCTGTTGGGAAAATCAAGGCAAAATGATGCTGACACTGGCTCA
AATGCTGACATGACCTACTCCATCATAAATGGTATGGCATGGGAATATTCTCAATCTCCACTGACAAA
GAGACCAGAGAAGGAATCCTTTCTTTAAAGAAGCCACTGAACTATGAGAAAAAGAAGTCATATACCTC
AACATAGAAGGAGCAAAACACATCTTGATTTTCGCTTTTCTCACTTGGGTCCTTTTAAAGATGCTACT
ATGCTGAAGATCATTGTTGGGGATGTAGATGAACCACCACTATTTCCATGCCTTCTACCTCATGGAA
GTCTACGAAAATGCCAAGATTGGGACCGTCGTTGGTACAGTTTTGGCACAAGATCCTGACAGTACTAAC
AGCTTAGTAAGATACTTCATCAACTACAATGTTGAAGACGACAGATTTTTCAACATTGATGCCAATACT
GGGACCATTAGGACTACAAAGTTCTCGACAGAGAAGAACTCCATGGTACAACATCACAGTCACTGCT
TCAGAAAATTGATAATCCTGATTTGCTGAGCCATGTCACAGTGGGTATTAGAGTTCTGGATGTCATGAC
AATCCACCCGAACCTGCCAGGGAATATGATATTATTGTATGTGAAAATCTAAGCCTGGCCAGGTTATT
CATACCATCAGTGCCACTGATAAAGATGATTTTGCCAATGGACCAAGGTTAACTTCTTTCTTGATGAA
CGCCTGCCTGTAATCCAACTTCACTCTGAAGGACAATGAAGCAGCAGCAGCACCCCTCACCATCAGGG
TTTGTGCATGCGAGAGAGATGGCGTGTGCGGACCTGCCATGCAGAAGCCTTCTGTCTCGGCTGGTT
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001291957

Insert Size: 1728 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291957.1](#)

RefSeq Size: 2974 bp

RefSeq ORF: 1728 bp

Locus ID: 1016

Cytogenetics: 5p14.3

Protein Families: Transmembrane

MW: 63.7 kDa

Gene Summary:

This gene encodes a type II classical cadherin from the cadherin superfamily of integral membrane proteins that mediate calcium-dependent cell-cell adhesion. Mature cadherin proteins are composed of a large N-terminal extracellular domain, a single membrane-spanning domain, and a small, highly conserved C-terminal cytoplasmic domain. Type II (atypical) cadherins are defined based on their lack of a HAV cell adhesion recognition sequence specific to type I cadherins. This particular cadherin is expressed specifically in the central nervous system and is putatively involved in synaptic adhesion, axon outgrowth and guidance. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2014]

Transcript Variant: This variant (4) differs in the 5' UTR and has an alternate splice site in the 3' coding region, which results in a frameshift, compared to variant 1. It encodes isoform 3, which has a shorter and distinct C-terminus, compared to isoform 1. Variant 4 and 8 encode the same protein (isoform 3). This isoform (3) may undergo proteolytic processing similar to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.