

Product datasheet for **SC336857**

TTC39A (NM_001297663) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39A (NM_001297663) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39A
Synonyms:	C1orf34; DEME-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336857 representing NM_001297663.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGACCTCGGCTGGCGGCCCGCCAGGAGCCCTGCCCGGGGACTCCTGAGAGCAGCTCCATGAGGCC
CTGGACCAGTGCATGACCGCCCTGGACCTCTTCCCTACCAACCAGTTCTCAGAAGCACTCAGCTACCTC
AAGCCCAGAACCAAGGAAAGCATGTACCACTCACTGACATATGCCACCATCCTGGAGATGCAGGCCATG
ATGACCTTTGACCCTCAGGACATCTGCTTGCCGGCAACATGATGAAGGAGGCACAGATGCTGTGTGAG
AGGCACCGGAGGAAGTCTTCTGTAACAGATTCTTTCAGCAGCCTGGTGAACCGCCACGCTGGGCCAA
TTCCTGAAGAGGAAATCCACGCTGAGGCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTGACC
TTCCTGCAGGACGAGAACATGGTGAAGTTCATCAAAGGCGGCATCAAAGTTCGAAACAGCTACCAGACC
TACAAGGAGCTGGACAGCCTTGTTCAGTCTCACAACTGCAAGGGTGAGAACCACCCGCACTTTGAA
GGAGGAGTGAAGCTTGGTGTAGGGGCTTCAACCTGACACTGTCCATGCTTCTACTAGGATCCTGAGG
CTGTTGGAGTTTGTGGGTTTTTCAGGAAACAAGGACTATGGGCTGCTGCAGCTGGAGGAGGAGCGTCA
GGGCACAGCTTCCGCTCTGTGCTGTGTGATGCTCCTGCTGTGCTACCCACCTTCTCACCTTGTG
CTCGGTACTGGGAACGTCAACATCGAGGAGGCCGAGAAGCTCTTGAAGCCCTACCTGAACCGGTACCCT
AAGGGTGCCATCTTCTGTTCTTTGAGGGGAGGATTGAAGTCATTAAGGCAACATTGATGCAGCCATC
CGGCGTTTCGAGGAGTGTGTGAGGCCAGCAGCACTGGAAGCAGTTCCACCACATGTGCTACTGGGAG
CTGATGTGGTGTTCACCTACAAGGCCAGTGAAGATGTCTACTTCTACGCCGACCTGCTCAGCAAG
GAGAACTGCTGGTCCAAGGCCACCTACATTTACATGAAGGCCCTACCTCAGCATGTTTGGGAAGGAG
GACCACAAGCCGTTCCGGGACGACGAAGTGAATTATTTTCAGCTGTGCCAGGCCGTAAGCTCAAGATT
GCTGGGAAATCTCTACCCACAGAGAAGTTGCCATCCGGAAGTCCCGGCGCTACTTCTCCTCAACCT
ATCTCGCTGCCAGTGCCTGCTCTGGAAATGATGTACATCTGGAACGGCTACGCCGTGATTGGGAAGCAG
CCGAAACTCACGGATGGGATACTTGAGATTATCACTAAGGCTGAAGAGATGCTGGAGAAAGGCCAGAG
AACGAGTACTCAGTGGATGACGAGTGTGGTGAATTGTTGAAAGGCCTGTGTCTGAAATACCTGGGC
CGTGTCCAGGAGGCCGAGGAGAATTTTAGGAGCATCTCTGCCAATGAAAAGAAGATTAATATGACCAC
TACTTGATCCCAAACGCCCTGCTGGAGCTGGCCCTGCTGCTTATGGAGCAAGACAGAAACGAAGAGGCC
ATCAAACCTTTGGAATCTGCCAAGCAAACTACAAGAATTACTCCATGGAGTCAAGGACACACTTTGCA
ATCCAGGCAGCCACACTCAAGCCAAGTCTTCCCTAGAGAACAGCAGCAGATCCATGGTCTCATCAGTG
TCCTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001297663

Insert Size: 1734 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297663.1](#)

RefSeq Size: 2738 bp

RefSeq ORF: 1734 bp

Locus ID: 22996

UniProt ID: [Q5SRH9](#)

Cytogenetics: 1p32.3

MW: 65.7 kDa