

Product datasheet for **SC336858**

CCDC138 (NM_001303105) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC138 (NM_001303105) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC138
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336858 representing NM_001303105.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGCCGAGGGTCGTCAAGCCACCGGGCAGGATTTAGTAGTGGAGAGTCTCAAAGCCGCTACGGA
CTCGGGGGCAGCTGCCCGACGAGTATGATTTTTCAAATTTTTATCAGTCTAAGTATAAGAGAAGAACT
CTAACCTCCCAGGTGATTTGGATATCTACTCTGGAGATAAAGTTGGTTCATCGTTAAAAATTCTGAT
GAAAGCAAGCATTGTAGAACACCATTGGGCAGCTTATTCAAGCACGTAATGTGAATTGCCTAGATGAT
GAACTGGATTCTTTCCATGATTTGAAGAAACAGGAAACAGAAGAAGGTTAATTGAAAATGATTATAGA
GTTAGTACCTCGAAAATAACCAAGCAGTCTTTAAAGAAATAGAAAAAGTTGCCTTGCCAACTAATACG
ACCTCATCGAGACCTCGGACTGAGTGTGTAGTGTGCAGGTGACTCTCCTTTGAAACCTGTCAGCTGT
CCAAAACTAAAGCATCAGACAAGCGGAGTTTACTTCCACATCAGATCAGTCAGATATATGACGAATTA
TTTCAGATACATCTGAAATTGCAGTGTGAAACTGCAGCACAACAGAAATTTGCTGAAAGAACTTCAAAG
CGAGAACGTTTTTACTTGAAAGAGAACAACACTGCTTTTCAGACATGAAAATGCCTTGAGTAAAATAAA
GGTGTTGAAGAAGAGGTTCTTACAAGATTTCAAATTTATAAAAGAGCAGCATGATGCAGAAGTTGAACAC
TTAACCGAAGTTCTTAAGGAAAAGAATAAAGAAACCAAGAGACTGAGGTCTCTTTTGATGCATTGAAA
GAATTGAATGATACCTTAAAAAACAGTTAAATGAAGCAAGTGAAGAAAACAGGAAGATAGACATTCAG
GCTAAAAGAGTTCAAGCTCGTTTACATAATTTACAGAGGAAGTACGAGTTTATGACAATACAGAGATTG
AAAGGAAGTTCCCATGCTGTTTCAAGTGAAGTGAAGTAAAACAAGAAAAGCACCAGTTTCAAAAACT
TACAAGGTACCACTTAATGGGCAAGTTTATGAACTTTTAACTGTCTTTCATGGACTGGATTTCCGATCAT
CATCTTAGCAAAGTGAACATGAAGAATCTGGAATGGATGGTAAAAACCACAACCTCAAATTTGCTTCC
CAGAGAAATGATATTCAGGAGAAGTGTGTAAGCTTTTGCCTCTAATGACAGAGCAGCTACAGTGGATG
CCATTTGTGAATATCAAACCTTACGAGCCTTTTGTAAAATTTATATATTGGTCCCTAAGGCAGCTAGAT
GCTGGAGCACAGCACTCGACTATGACATCAACATTGAGGAGATTGGGTGAAGACATTTTTAAAGGAGTG
GTAACAAAGGAATTCAGGATAATTCTCCACAGCATTCTGTGGAGAATAAACCAAGACAGCTGCTTTC
TTAAGAGCTCCAATTTGCCATTGAGATTTTTATCAACCTTAATTGTTCTAAAACAGTCACTCAAGCT
GATTACCTGGCTCAGGCATTTGATTCTCTTTGTTGGACTTGAAGACAGAAGAAGGAAAAACCTTGTTT
TTGGAGTATCAGGCTGTTCCAGTAATTAAGTCATCTAAGAATATCCAGTAAAGGACTCCTGTCTAAT
GTTATTGATAGTTTGTCCAGATGACGGTGAATCTAGAGTAATAAGAAGCTCTTTGAACTTTTACGA
TTCATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001303105

Insert Size: 1734 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001303105.1](#)

RefSeq Size: 2281 bp

RefSeq ORF: 1734 bp

Locus ID: 165055

UniProt ID: [Q96M89](#)

Cytogenetics: 2q13

MW: 66.4 kDa