

Product datasheet for **SC336861**

SLFN12 (NM_001289009) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLFN12 (NM_001289009) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLFN12
Synonyms:	SLFN3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336861 representing NM_001289009.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACATCAGTGTGTTGATTTGAAACGAATTATGCCGAGTTGGTCTAGATGTGGGAAGAGTCACTCTT
GGAGAGAACAGTAGGAAAAAATGAAGGATTGTAACCTGAGAAAAAGCAGAATGAAAGTGTCTCACGA
GCTATGTGTGCTCTGCTCAATTCTGGAGGGGAGTGATCAAGGCTGAAATTGAGAATGAAGACTATAGT
TATACAAAAGATGGAATAGGACTAGATTTGAAAAATTCTTTTAGTAACATTCTGTTATTTGTTCTGAG
TACTTAGACTTCATGCAGAATGGTAACTACTTTCTGATTTTGTGAAGTCATGGAGCTTGAACACCTCT
GGTCTGCCGATTACCACCTTGAGCTCCAATTTGTACAAAAGAGATATAACATCTGCAAAAGTCATGAAT
GCCACTGCTGCACTGGAGTTCCTCAAAGACATGAAAAAGACTAGAGGGAGATTGTATTTAAGACCAGAA
TTGCTGGCAAAGAGGCCCTGTGTGATATAACAAGAAGAAAATAACATGAAGGCCTGGCCGGGTTTTT
TTTGATAGAACAGAACTTGATCGGAAAGAAAAATTGACCTTTACTGAATCCACACATGTTGAAATTTAA
AACTTCTCGACAGAAAAGTTGTTACAACGAATTAAGAGATTCTCCCTCAATATGTTTCTGCATTTGCA
AACTACTGATGGAGGATATTTGTTCAATTGGTTTAAATGAAGATAAAAGAAATAATTGGCTTTAAAGCAGAG
ATGAGTGACCTCGATGACTTAGAAAAGAGAAATCGAAAAGTCCATTAGGAAGATGCCTGTGCATCACTTC
TGTATGGAGAAGAAGAAGATAAATTATTCATGCAAATTCCTGGAGTATATGATAAAGGAAGTCTTTGT
GGATATGTCTGTGCACTCAGAGTGGAGCGCTTCTGCTGTGCAGTGTGGTAAAGAGCCTGATTCCTGG
CATGTGAAAGATAACCGTGTGATGCAGTTGACCAGGAAGGAATGGATCCAGTTCATGGTGGAGGCTGAA
CCAAAAATTTCCAGTTCATATGAAGAGGTGATCTCTCAAATAAATACGTCATTACCTGCTCCCCACAGT
TGGCCTCTTTTGGAAATGGCAACGGCAGAGACATCACTGTCCAGGGCTATCAGGAAGGATAACGTATACT
CCAGAAAACCTTTGCAGAAAACCTGTTCTTACAACATGAAGGACTTAAGCAATTAATATGTGAAGAAATG
GACTCTGTGAGAAAAGGGCTCACTGATCTTCTCTAGGAGCTGGTCTGTGGATCTGGGCTTGAAGAGAAC
CACAAAGTCTCTGTGATGCTTCTGATTTCCAGGACAGTCCCTCCAGTCTATACACCTTCCACATG
GTACAGGATGAGGAGTTTAAAGGCTATTCTACACAACTGCCCTAACCTTAAAGCAGAAGCTGGCAAAA
ATTGGTGGTTACACTAAAAAGTGTGTGTCATGACAAAGATCTTCTACTTGAGCCCTGAAGGCATGACA
AGCTGCCAGTATGATTTAAGGTCGCAAGTAATTTACCCTGAATCCTACTATTTTACAAGAAGGAAATAC
TTGCTGAAAGCCCTTTTTAAAGCCTTAAAGAGACTCAAGTCTCTGAGAGACCAGTTTTCTTTGCAGAA
AATCTATACCAGATAATCGGTATAGATTGCTTTCAGAAGATGATAAAAAGATGTTTAAATCTTGTGCA
AGGCTCACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001289009

Insert Size: 1737 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289009.1](#)

RefSeq Size: 2465 bp

RefSeq ORF: 1737 bp

Locus ID: 55106

UniProt ID: [Q8IYM2](#)

Cytogenetics: 17q12

MW: 67 kDa