

Product datasheet for **SC336865**

LACTBL1 (NM_001289974) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LACTBL1 (NM_001289974) Human Untagged Clone
Tag:	Tag Free
Symbol:	LACTBL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336865 representing NM_001289974.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAAGACCCAGGATGGCTGGTGGCCAACTTGCTGAAAGTGAAGAAGAAGTGGCTCTCCCGGCCATC
TGCAGCTTCTTCTGCTGCTCAGTGGTCATGACTGGCTGCTTCTGTGGCAGTATCACCTCCCAAG
CTGAAGACCGGTTCCCTGGGACCAAGAGGACCTCTGCCCTGTGAGGATGTGCCCGGCACCCTGAG
CCCGTGCCCTGGCTCACCCACTCCCGTGTTGAAGGAGGCCCTGGAAAAGTGGACCAGATCCTGCGC
CAGGCAATGTCTGCCCCAGGCGTGGCTGCCATGTCTGCAGTTGTCATCCACAATGATACTGTGCTCTGG
ACAGGGAACTTTGGGAAGAAGAATGGCTCAGACCCGGCTTCTGGGGCCCCCAATGAGTACACCATGTAC
AGGATCTCCAGCATCTCCAAGATCTTTCCTGTCCTCATGTGTACCGTCTGTGGGAGGAGGGCATCGTG
GCCTCCCTAGATGACCCTCTGGAGCGGTATGCCAGCACCTTACCATTAAACAACCCGCTGGGCTGGCA
TCAGCCGAACAGCAGGGCTGATGGACGGGCTGGAGCAGGTGGGCCCGCCCCAAGGCCCTCACCTGTC
ACCCTTCAAGGATGGCCAGCCAGCTCTCAGGGCTGCCAGAAGGCTGCGCTCCACTTCACTGCTGTGG
AAGGGCAGCACCCAGGAGGCCCTGAACCTGCTCAAGGACGATGTGCTGGTGTGGACCCGGGAACACAGA
TGCCATTACAGCACGCTGGCCTTCTCGCTCCTGGCCACGTCTGGCAGCTCACACCCGCCAGGGCGAC
TACCAGCGCTGGGTCTCGGAGAAGCTGCTGGAGCCGCTGGGGATGGCAGACACGGGCTTTGACCTCACG
CCCAGCTGCGCGCGCCTGGCCGCGGGCTTCTACGGCAGCGGGCGGCCGGCCACTCTATGACCTG
GGCTGGTACCGGCCGTCGGGCCAGATGTACTCTACCGCCCGGACCTGGCCAAGCTGGCCGTGGCGCTC
CTGGGGCGCGGGCCCCGGCGGCTCCTGCGGCCGACGCGGCCAAGACGCTGCTGGCGCCGCTGCTGGCC
TGCCCGGGCGCCTACTTCGCCAACGAGACGGGCACGCCGTGGGAGTTCCACGCGCAGCGGGGTACCGC
GTGGTGGCAAGGACGGCGATCTGGACGGCTACGCCGCCACCTTCTCGCTGGTGGCCCGCTGCGTCTG
GGCTGGTGTGCTGCTGGCAGGGCCACGGCCGCGCCGGGCCGACCTGGTGGCGCGGGCCTACGATGAG
CTCCTGCCCGCCTGGAGCGCGCCTCCGGGAGGCCGAGCCGGCCGGCTCCGCCGCCACCGCGCAC
CCCTTCGCCGGCTACTTCACCTTCGCCAACCTGACCTTCTACGAGGTGCGCGCGGGCGCGGGCGAG
CTGCGCCTGCGCCAGTTCGGGCCGCGGTGGAGGCGTGGTGCCTCCAGCGTTCGCGACCCTGGCGCTG
CGCCACCTGCACGGCCGCTCTTCAGCTGCACGTGGCCACGAGTTCCTCCTGCGCACTGCCACTCGGC
GACCGTGGCTCTGCTCGAGGCCAGCACGGGCGAGTGGTCAACTTCTACCCCTGGATACACACGGG
CTGTCACCCGGCTTCGACGTGCCCGCCTCAACACGTACAGAGTGTGCGGCTGCGGGGCAAGCCGGTG
TTCAAGACCCAGTGA
AGCGGACCGACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_001289974

Insert Size: 1740 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289974.1](#)

RefSeq Size: 1740 bp

RefSeq ORF: 1740 bp

Locus ID: 646262

Cytogenetics: 1p36.12

MW: 63.5 kDa