

Product datasheet for **SC336870**

DGLUCY (NM_001286472) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DGLUCY (NM_001286472) Human Untagged Clone
Tag:	Tag Free
Symbol:	DGLUCY
Synonyms:	C14orf159
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336870 representing NM_001286472.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCTTCACACTCCACCTGAGGTCCCGCCTCCCTCGCCATAAAGAGTTTATTCTACAAAAGAAA
CCAAACATCAGAAATACATCCAGCATGGCTGGAGAGCTCCGACCAGCCAGCCTGGTGGCTCCTGCCAGG
TCCTTGCTCCAGCTTTTGAAGATTCTGCCAGGTCAACACTGGTCTCTACCCCTGCTGGGCCAGAGT
GAGCCAGAAAAGTGGATGCTGCCCCCTCAAGGTGCTATCTCAGAGACCAGGATGGCCATCCCCAGTTC
TGGAAATACGAGTTCCGTGCCTGCACCGGCAGCCTGGCTTCGCTGGAGCAGTACTCGGAGCAGTGAAG
GACATGGTGGCCTTCTTCTGGGCTGCAGCTTCTCCCTGGAGGAGGCCTTGGAGAAAGCGGGGCTCCCC
AGAAGAGACCCAGCAGGTACAGCCAGGCGGGTGCATACAAGACAACAGTGCCTTGTGTACCCATGCT
GGCTTCTGCTGCCCTCTGGTGGTCACGATGAGGCCATTCCCAAGGACAAGCTGGAAGGGCTGGTGGC
GCCTGCTGCTCCCTCGGAGGTGAGCAGGGGCAACCTGTTACATGGGCGACCCAGAAGTGTGGGAATC
AAAGAGCTTTCCAAACCTGCCTACGGGGATGCCATGGTGTGTCCCCAGGGGAGGTTCCAGTGTCTGG
CCTTCTCCGCTGACCAGTCTCGGAGCTGTGAGCAGCTGTGAGACCCCACTGGCTTTTGCCAGCATCCA
GGCTGCACAGTTATGACTGACCTGAAGGATGCAAAGGCTCCACCTGGTTGTCTCACCCAGAGAGAATT
CCAGAGGTCCATCACATTTCCAAAGATCCTCTGCACTACAGCATCGCGTCAGTCTCTGCTTCTCAGAAG
ATCAGAGAACTAGAGTCTATGATCGGCATAGACCCAGGGAACCGGGGATTGGGCACCTGCTCTGTAAA
GATGAGCTGCTGAAGGCCTCTCTCTCGCTGTCCCATGCCCGCTCAGTGCATCACCACTGGGTTCCCC
ACACATTTCAATCATGAGCCTCCAGAAGAGACAGATGGCCACCAGGAGCTGTTGCTCTGGTTGCCTTC
CTGCAGGCCTTGGAGAAGGAGTGCCTAATCGTTGACCAGAGAGCCTGGAAGTGCACCAGAAGATT
GTTGAAGATGCTGTTGAGCAAGGTGTTCTGAAGACGCAGATCCCGATATTAAGTACCAAGTGGATCA
GTGGAAGCTGCTCAGGCATTCTGTGCAAAAATGGGGACCCGCAGACACCTAGATTTGACCACCTGGTG
GCCATAGAGCGTGCCGGAAGAGCTGCTGATGGCAATTACTACAATGCAAGGAAGATGAACATCAAGCAC
TTGGTTGACCCCATGACGATCTTTTTCTTGCTGCGAAGAAGATTCTGGAATCTCATCAACTGGTGT
TCTAACTGGGAGGCTATGCCCTGGCCTGCGCACTCTACATCCTGTACTCATGTGCTGTCCACAGTCAG
TACCTGAGGAAAGCAGTCGGACCCTCCAGGGCACCTGGAGATCAGGCCTGGACTCAGGCCCTCCCGTCG
GTCATTAAGGAAGAAAAATGCTGGGCATCTTGGTGCAGCACAAAGTCCGGAGTGGCGTCTCGGGCATC
GTGGGCATGGAGGTGGATGGGCTGCCCTCCACAACCCACGCCGAGATGATCCAGAAGCTGGTGGAC
GTCACCACGGCACAGGTGTA
```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001286472

Insert Size: 1746 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286472.1](#)

RefSeq Size: 2832 bp

RefSeq ORF: 1746 bp

Locus ID: 80017

UniProt ID: [Q7Z3D6](#)

Cytogenetics: 14q32.11

MW: 62.8 kDa

Gene Summary:

D-glutamate cyclase that converts D-glutamate to 5-oxo-D-proline.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (8) differs in the 5' UTR, uses an alternate in-frame splice site, and lacks an alternate in-frame exon, compared to variant 1. The encoded isoform (d) is shorter than isoform a.