

Product datasheet for **SC336880**

KLHL32 (NM_001286250) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | KLHL32 (NM_001286250) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | KLHL32 |
| Synonyms: | BKLHD5; dJ21F7.1; KIAA1900; UG0030H05 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |




[View online »](#)

Fully Sequenced ORF: >SC336880 representing NM_001286250.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCGTCTGAACGCTGCCTCAGTATTCAAGAAATGCTGACAGGCCAGAGGCTCTGCCACTCCGAATCT
CACAATGACAGTGTCTGGCAGCGCTGAATCAGCAGAGGAGTGATGGCATCCTCTGCGACATCACCCCTG
ATTGCTGAGGAACAGAAATTCATGCTCACAAGGCAGTCTAGCAGCATGCAGTACTATTCCGGATT
TTGCTGGAGCCAGGTGTGATCCAGGATGTGCTAGCAGCGGGCAGTCACCTACAGCTGTTGGAGCTTCTC
AATTTATGCTCCCACTATCTCATCCAGGAATTAATAGCTTTAATTACTTGGATCTGTACAGACTTGT
GACCTCTTTAACCTCACTTTGTTGGAGAAGGCAGTGATCGATTTCTTAGTGAAACATCTCTGAACTC
CTGAAGAGCCGCCAGAAGAAGTTCTAACGCTTCCCTATTGCCTGCTTCAGGAGGTGCTGAAGAGCGAC
CGCCTGACCTCCCTGAGTGAAGAGCAGATCTGGCAGTAGCTGTGAGGTGGTTGGAACACAAGTCCAC
TACCAGTACATGGACGAGCTCTGCAATACATCCGCTTTGGCCTAATGGATGTGGATACTCTCCATACA
GTTGCCCTGTCCCACCCCTTGTCCAAGCAAGTGAGACTGCAACAGCCCTTGTCAACGAGGCCCTGGAA
TACCACCAGAGCATCTATGCACAGCCTGTCTGGCAGACTCGCAGGACCAACCACGATTCCAGTCAGAC
ACTCTGTATATCATTGGTGGGAAAAAGCGCGAGGTCTGCAAGGTCAAGGAACCTCGGTACTTCAATCCT
GTTGATCAGGAGAATGCTCTCATAGCTGCCATTGCCAAGTGGAGTGAGCTGGCTCCCATGCCGTGGGA
AGGAGCCACCATTTGTGTGGCAGTCATGGGGACTTCTGTTTGTGGCAGGAGGGGAAGTTGAGCATGCC
AGTGGCCGGACGTGTGCTGTGAGGACTGCCTGTGCTATGACCCCGCAGTAATTCCTGGGCAGAGATA
GCACCCATGAAAAACTGCCGGGAGCATTGTGCTGGGTGCCATGGAGGAATACCTCTATGCAGTTGGG
GGCAGAAATGAACTGCCCGAGTTCTGCCTACAGTTGAGCGATATTGCCCAAGAAGAACAATGGACT
TTTGTTCAGTCCTTTGACAGATCCCTTTCATGCCATGCTGGATATGTGGCTGATGGTCTTCTTTGGATA
TCAGGTGGAGTAACTAATACGGCACAATATCAGAACAGGCTAATGGTGTATGAACCTAACCAAGAATAAG
TGGATAAGCCGTAGCCCATGCTGCAGAGAAGGTCTACCATTCCATGGCTGCTGTACAAAAGGAAGCTT
TATGTTCTTGAGGCAATGACCTAGACTACAATAATGACCGGATCCTTGTGCGCCATATAGATTCTTAC
AACATAGACTGACAGTGGACACGTTGTAATTTCAACCTGCTGACTGGCCAAAATGAATCTGGAGTT
GCTGTCCATAATGGGAGAATATATTTAGTTGGTGGATTTCAATTTGGACAAATGAGCCTCTGGCTTGT
ATCCAGGTAAGTGTAAAGCAGAGAAGGCAAGAAGAAGTATTCTATGGGCCTACACTCCCTTTTGT
TCCAATGGAATAGCAGCATGCTTCTTCCAGCTCCATATTTACATGCCCTAACCTTCAAACCTCTTCAA
GTGCTCATCACAGGATTGGCACCATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001286250

Insert Size: 1755 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286250.1](#)

RefSeq Size: 3699 bp

RefSeq ORF: 1755 bp

Locus ID: 114792

UniProt ID: [Q96NJ5](#)

Cytogenetics: 6q16.1

MW: 66.6 kDa