

## Product datasheet for **SC336902**

### ZNF285 (NM\_001291489) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF285 (NM\_001291489) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF285  
**Synonyms:** ZNF285A  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC336902 representing NM\_001291489.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGATTAAGTTCAGGAAAGGGTGACATTCAGGATGTGGCTGTTGTCTTCACCAAGGAAGAGCTGGCA
CTATTGGATAAAGCCCAGATAAACCTGTACCAAGATGTGATGCTGGAAAACCTCAGGAACCTCATGTTA
GTGAGAGACGGGATTAACAACATTTTGAATCTTCAGGCAAAGGGGTTAAGTTACCTTCGCAAGAA
GTGCTTCATTGCTGCCAGATTTGAAACAAAGGATCCGGGATTTAACTGTGAGTCAGGATTATATCGTG
AACCTTCAAGAAGAGTGTCCCCACATTTAGAAGATGTTCCCTCAGTGAAGAGTGGCAGGCATTTCT
CTTCAGATTTCTGAAAATGAAAACCTATGTAGTAAATGCCATTATCAAAAATCAAGATATCACAGCATGG
CAAAGCCTGACACAGGTTCTTACCCAGAATCGTGGAGGAAAGCCAACATAATGACCGAGCCCCAGAAC
TCTCAGGGAAGATATAAGGGAATTTACATGGAAGAGAAATTGTACAGACGTGCTCAGCATGATGACAGC
CTCAGTTGGACCTCATGTGATCATCATGAGTCCCAAGAATGTAAGGAGAGAGGACCCTGGTAGACATCCC
AGCTGTGGGAAAAACTTGGGTATGAAATCAACGGTTGAAAAACGTAATGCGGCCCATGTATTACCACAG
CCTTTCCCATGTAATAACTGTGGGGTGGCCTTTGCAGATGATACAGATCCTCATGTCCATCACAGCACT
CACCTAGGAGAAAAATCTTATAAATGTGACCAGTATGAAAGAAGTCTTAGTCAGAGCCAAGATCTTATC
GTTCAATGTAAGTCACTCTGGCAAGACTCCCTATGAATCCACGAATGGCCTATGGGCTGCAAAACAG
AGCTCAGACCTTCCAGATATCAGAAAGTCTCCTCAGGAGACAAACCTACAAATGTAAGAATGTGGC
AAGGGCTTCAGGCGCAGCTCTTCCCTTCAACCATCATCGAGTCCACACAGGGGAGATGCCTACAAA
TGCGATGAATGTGGGAAAGGGTTTGATTTAGGTCACTTCTTTGTATTATCAGGGAGTACACACAGGG
AAAAAGCCCTATAAATGTGAAGAGTGTGGGAAAGGGCTTTGATCAGAGCTCCAACCTTCTTGTCCATCAG
AGAGTCCACACTGGAGAGAAGCCCTACAAATGCAGTGTGAGTGTGGCAAGTGTCTTGTAGTCAAGCTCCGTT
CTTCAAGTCCACTGGAGGTTTACACAGGGGAGAAACCATATAGGTGTGGTGTGAGTGTGAAAGGGCTTC
AGCCAATGTACACACCTTACATTCACCAGAGAGTCCACACAGGGGAGAAACCATACAAATGCAATGTG
TGTGAAAGGATTTGCGTATAGCTCTGTTCTTACACTCATCAGAGAGTTCACACTGGAGAAAAACCA
TATAAATGTGAAGTGTGTGAAAGTGTTCAGTTACAGTTCATATTTTCACTTACATCAAAGAGATCAC
ATCAGAGAGAAACCATATAAATGTGATGAGTGTGGTAAAGGCTTCAGCCGGAATTCAGATCTTAATGTT
CACCTCAGAGTCCACACAGGAGAGAGGCCCTATAAGTGTAAAGGCATGTGGTAAGGGCTTCAGTCGTAAT
TCATACCTCCTTGCCCATCAGAGAGTGCATATAGATGAGACACAGTACACACACTGTGAGCGTGGAAAG
GACCTTCTGACTCATCAAAGACTACATGAGCAGAGAGAAACATTATAA
  
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001291489

**Insert Size:** 1773 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001291489.1](#)

**RefSeq Size:** 2794 bp

**RefSeq ORF:** 1773 bp

**Locus ID:** 26974

**UniProt ID:** [Q96NJ3](#)

**Cytogenetics:** 19q13.31

**Protein Families:** Transcription Factors

**MW:** 68.1 kDa

**Gene Summary:**

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) uses an alternate splice site in the 5' UTR and lacks an alternate in-frame exon in the 5' coding region, compared to variant 1. The encoded isoform (b) is shorter than isoform 1. Variants 2 and 3 encode the same isoform (b).