

## Product datasheet for **SC336908**

### **TMEM132B (NM\_001286219) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	TMEM132B (NM_001286219) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM132B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC336908 representing NM\_001286219.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGAGCAAAGTGGACACGATTGTGAAC TTCACCCACCAGCACTTCACCTCCCAGTTCGAGGTCACT
GTCTGGGCACCCAGGCTCCCCGCAGATTGAGATCTCAGACACCGAGCTGAGCCAGATCAAGGGCTGG
AGGATCCCGGTTGCTGCCAACAGAAGGCCCTACCCGGGAAAGCGATGACGAGGACGATGAGGAGAAGAAG
GGACGAGGCTGCTCCCTGCAGTACCAGCACGCCACAGTGCCTGCTCCTCACCCAGTTTGTGGCCGAGTCA
CCTGACTTAGGGCAGCTGACCTACATGCTGGGCCCGACTGGCAGTTTGACATCACTGACCTTGTGACC
GAGTTCATGAAGGTGGAGGAGCCGAAAAATCGCTCAGTTACAGGACGGCAGGACCTGGCTGGTCGGGAG
CCGGGAATAACCCAGGTGCAGGTCTCTCGCCGTTGTCTGACTCCATCCTGGCTGAGAAGACGGTGATT
GTCCTGGATGACCGAGTACCATCGCGGAGCTGGGAGTGCAGCTCGTAGCTGGCATGTCTCTCCCTG
CAGCCACACCGAGCAGACAAAAGGGCCATCGTCTCCACAGCTGCTGCCCTGGATGTTCTTCAGTCCCCA
CAGCAGGAAGCAATAGTAAGTTCTTGGATTTTGTTCAGTATGGTTCGGTGACACCTTTAGACATTTAC
GATCCTAAGGATTATTCTGTTACTGTCTCATCATTGGATGAAATGGTGGTGTCTGTCCAGGCAAACCTT
GAGTCCAATGGCCAATTGTGGTTGCAGAGGGTGAAGGACAAGGGCCTTTGATTAAGTTAGAAAATGATG
ATAAGTGAACCTTGTGAGAAGACCAAGAGGAAGAGTGTCTTGCCGTGGGTAAAGGAAAATGTCAAGGTC
AAATTCGAACCAAGTAGTGATGAGCACCAAGGAGGCAGCAATGATATTGAGGGCATAAAATCGGGAATAT
AAAGACCACCTCAGTAATTCCATAGAGCGCGAAGGAAACCAGGAGAGAGCAGTCCAGGAATGGTTCCAC
CGTGGCACACCTGTTGGCCAAGAGGAAAGTACCAACAAAAGCACAACCCCCAGTCTCCCATGGAAGGG
AAGAATAAGTTACTCAAAGTGGTGGTCCAGATGCCTTTACAAGCTTCCCCTCAAGGGAAGTCACCG
GACCCCAATAATCCTAGTGACCTCACAGTGACCTCAAGGGGGCTAACGGACTTGGAGATTGGCATGTAT
GCCTTGCTCTGCGTCTTCTGTCTGGCCATTCTGGTCTTCTTGATCAACTGCGTGGCGTTTGCCTGGAAA
TACAGACACAAAAGGTTTGTGTGAGTGAGCAGGGCAACATCCCCATTCCCACGACTGGGTCTGGCTT
GGGAATGAAGTGGAACTTTGGAGAACCCTGTTGACATTACACTCCCATCAGAGGAGTGCACAACCATG
ATAGACAGGGCCTGCAGTTCGAGGAGAGGAACTTCTTCTGAATGGCAGTTCACAGAAAGACTTTTCAT
AGTCAACTACTCAGACCCTCTGACTATGTCTATGAGAAAGAAATTAATAATGAACCTATGAATCTTCG
GGCCAAAGAGGAAGAGAGTCAAGTTCACCTTCTACACCACCATCCTCCCAGAGGACGGCGGCCATAC
ACCAACTCCATCCTGTTTGACAGCGATGATAACATCAAGTGGGTCTGCCAAGATATGGGGCTGGGGAT
TCACAGGACTTTAGAGACTATATGAAAAGTCTGCAAGACCAGATGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001286219

**Insert Size:** 1773 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286219.1](#)

**RefSeq Size:** 9659 bp

**RefSeq ORF:** 1773 bp

**Locus ID:** 114795

**UniProt ID:** [Q14DG7](#)

**Cytogenetics:** 12q24.31-q24.32

**Protein Families:** Transmembrane

**MW:** 66.1 kDa