

Product datasheet for **SC336913**

ZNF570 (NM_001300993) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF570 (NM_001300993) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF570
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336913 representing NM_001300993.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGCATGTATTTGTTCTTTCTGGACTTCTCGATTGTCTGGGAGTGGGCTGAGGAGTGGCGTGTGGGTC
TCCGGAAGCTCGTCGAGGCCATCTGTGTGACTCCGGTCCGAGTGGAGTCATCTGAGGCCACTGCTATT
TCCCAAGAGAAGAGCCAGGAGGAAGAAAGAAATGGCTGTTGGGCTTCTTAAAGCCATGTACCAGGATTG
GTGACCTTCAGAGATGTGGCTGTAGACTTCTCCCAAGAGGAATGGGATTGTCTGGATTCTTCTCAAAGA
CATCTGTACAGTAATGTGATGCTAGAGAACTACAGGATCTTGGTATCACTGGGACTTTGCTTTTCCAAA
CCAAGTGTGATATTATTGTTGGAACAAGGAAAAGCACCCTGGATGGTGAAGAGAGAGCTGACAAAAGGC
TTGTGCTCAGGCTGGGAGCCTATATGTGAGACTGAAGAATTAACCCCAAAGCAGGATTTTATGAAGAA
CATCAATCCCAGAAGATAATAGAAACACTTACAAGCTATAACCTTGAATACTCCAGTTTGAGAGAAGAG
TGGAAATGTGAGGCTATTTTGAAGGCAACCAGGTAATCAGAAGGCGTGTTCAGGAAGAGATAATC
ACTCATGAAGAACCCTTTTGTGAGAGAGAACAAGAATATAAATCTTGGGGAAGTTTTCATCAGAAC
CCTGCTTTGTACACAAAAGATAATCCCAAGAGGAGAAAGTACATAAACATGACACACAAAAGAGA
AGCTTTAAAAAAATTTAATGGCTATTAAGCCCAAGAGTGTCTGTGCAGAGAAGAACTTTTGAAATGT
AATGACTGTGAAAAAGTCTTCAGCCAGAGTTCATCCCTTACTCTTCATCAAGAATTCTACTGGAGAG
AAACCCTATAATGTATAGAGTGTGAAAAGCCTTCAGCCAGAGATCAAACTTGTTCACATCAGAGG
ATTCATACTGGAGAAAAACCTATGAATGTAAGGAATGTAGGAAAGCCTTCAGTCAGAATGCACACCTA
GTTCAACATCTGCGAGTTCATACTGGAGAAAAACCTACGAATGTAAGGTATGTCGAAAAGCCTTCAGC
CAGTTTGCCTACCTTGCTCAACATCAGAGAGTTCACACGGGAGAGAAAACCTATGAATGTATCGAATGT
GGGAAAGCATTTAGCAACAGATCATCCATTGCTCAACACAGAGAGTTCATACAGGAGAGAAAACCTAT
GAATGTAATGTCTGTGGGAAAGCATTTAGCCTTCTGTCATACCTTACTGTACATCAGAGAATACATACT
GGAGAGAGACCCTATGAATGTAAGGAATGTGGGAAAGCCTTTAGCCAGAATTCACACCTTGCTCAACAT
CAGAGAATTCATACTGGAGAAAAACCTTATAAGTGTGAGGAATGTAGGAAAGCATTTCAGCCAGATTGCC
TACCTTGCTCAGCATCAAAGAGTTCATACTGGAGAGAAAACCTATGAATGTATTGAATGTGGGAAAGGCT
TTAGCAATGACTCGTCCCTTACTCAACATCAGCGAGTTCATACTGGAGAGAAAACCTTATGAATGTACT
GTTTGTGGAAAGGCTTTTAGTTACTGTGGATCCCTTGCCCAACATCAGAGAATTCATACTGGAGAGAGA
CCCTATGAATGTAAGGAATGCAAAAAACCTTCAGGCAGCATGCACACCTTGCTCATCACCAGAGAATT
CACATTGGGAGTCACTGTCACCACCAACCCAGTCAATCACCAGTCTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001300993

Insert Size: 1779 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001300993.2</u>
RefSeq Size:	5334 bp
RefSeq ORF:	1779 bp
Locus ID:	148268
Cytogenetics:	19q13.12
MW:	68.4 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.