

Product datasheet for **SC336915**

PRMT5 (NM_001282955) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRMT5 (NM_001282955) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRMT5
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336915 representing NM_001282955.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGCATGGCGGTGGGGGTGCTGGTGGGAGCCCGTGTCCAGCGGGAGGGACCTGAATTGCGTC
CCCGAAATAGCTGACACACTAGGGGCTGTGGCCAAGCAGGGGTTGATTTCTCTGCATGCCTGTCTTC
CATCCCGTTTTCAAGAGGGAGTTTATTAGGAACCTGCTAAGAATCGGCCGGTCCCGAGACACGATCA
GACCTACTGCTGTGAGGAAGGGCTTTCTGCTGCCCTTAATCAGGAAGATAACACCAACCTGGCCAGA
GTTTTGACCAACCACATCCACACTGGCCATCACTCTTCCATGTTCTGGATGCGGGTACCCTTGGTGGCA
CCAGAGGACCTGAGAGATGATATAATTGAGAATGCACCAACTACACACACAGAGGAGTACAGTGGGGAG
GAGAAAACGTGGATGTGGTGGCACAACCTCCGGACTTTGTGTACTATAGTAAGAGGATTGCAGTGGCT
CTTGAAATGGGGCTGACCTCCCATCTAATCATGTCTTGTGCTGCTGGTGGGGAGCCCATCAAAGCA
GCCATTCTCCCGACTAGCATTCTCTGACCAATAAGAAGGGATTTCTGTCTTTCTAAGATGCACCAG
AGGCTCATCTCCGGCTCCTCAAGTTGGAGGTGCAGTTTATCATCACAGGCACCAACCACACTCAGAG
AAGGAGTTCTGCTCTACCTCCAATACCTGGAATACTTAAGCCAGAACCCTCCTCCACCTAATGCCTAT
GAACTCTTTGCCAAGGGCTATGAAGACTATCTGCAGTCCCGCTTCAGCCACTGATGGACAACTTGAA
TCTCAGACATATGAAGTGTGAAAAGGACCCCATCAAATACTCTCAGTACCAGCAGGCCATCTATAAA
TGTCTGTAGACCGAGTACCAGAAGAGGAGAAGGATACCAATGTCCAGGACTGATGGTGTGGGAGCA
GGACGGGGACCCCTGGTGAACGCTTCCCTGCGGGCAGCCAAGCAGGCCGACCGGGCGGATAAAGCTGTAT
GCTGTGGAGAAAAACCAATGCCGTGGTACGCTAGAGAAGTGGCAGTTTGAAGAATGGGAAGCCAA
GTGACCGTAGTCTCATCAGACATGAGGGAATGGTGGCTCCAGAGAAAGCAGACATCATTGTCAGTGAG
CTTCTGGGCTCATTGCTGACAATGAATTGTCGCTGAGTGCCTGGATGGAGCCAGCATTCTAAAA
GATGATGGTGTGAGCATCCCCGGGGAGTACACTTCTTTCTGGCTCCCATCTTCTCTCAAGCTGTAC
AATGAGGTCCGAGCCTGTAGGGAGAAGGACCGTGACCCTGAGGCCAGTTTGAAGTGCCTTATGTGGTA
CGGCTGCACAACCTCCACCAGCTCTCTGCACCCAGCCCTGTTTACCTTCAGCCATCCCAACAGAGAT
CCTATGATTGACAACAACCGCTATTGCACCTTGGAAATTTCTGTGGAGGTGAACACAGTACTACATGGC
TTTGCCGGCTACTTTGAGACTGTGCTTTATCAGGACATCACTCTGAGTATCCGTCCAGAGACTCACTCT
CCTGGGATGTTCTCATGGTTTCCCATCTCTTCCCTATTAAGCAGCCATAACGGTACGTGAAGGCCAA
ACCATCTGTGTGCTTTCTGGCGATGCAGCAATCCAAGAAGGTGGTATGAGTGGGCTGTGACAGCA
CCAGTCTGTTCTGCTATTATAACCCACAGGCCGCTCATATAACCATTGGCCTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001282955

Insert Size: 1782 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282955.1](#)

RefSeq Size: 2409 bp

RefSeq ORF: 1782 bp

Locus ID: 10419

UniProt ID: [O14744](#)

Cytogenetics: 14q11.2

Protein Families: Stem cell - Pluripotency

MW: 67.7 kDa

Gene Summary:

This gene encodes an enzyme that belongs to the methyltransferase family. The encoded protein catalyzes the transfer of methyl groups to the amino acid arginine, in target proteins that include histones, transcriptional elongation factors and the tumor suppressor p53. This gene plays a role in several cellular processes, including transcriptional regulation, and the assembly of small nuclear ribonucleoproteins. A pseudogene of this gene has been defined on chromosome 4. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]

Transcript Variant: This variant (5) lacks an exon and uses an alternate splice site in the 5' coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (e) is shorter than isoform a.