

## Product datasheet for **SC336932**

### ZNF285 (NM\_001291488) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF285 (NM_001291488) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF285
Synonyms:	ZNF285A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336932 representing NM\_001291488.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATTAAGTTCAGATTCAAGGGATACATGTGCAGGAAAGGGTGACATTCAAGGATGTGGCTGTTGTC
TTCACCAAGGAAGAGCTGGCACTATTGGATAAAGCCAGATAAACCTGTACCAAGATGTGATGCTGGAA
AACTTCAGGAACCTCATGTTAGTGAGAGACGGGATTAACAACAACATTTGAATCTTCAGGCAAAGGGG
TTAAGTTACCTTCGCAAGAAGTCTTCATTGCTGGCAGATTTGAAACAAGGATCCGGGATTTAACT
GTGAGTCAGGATTATATCGTGAACCTTCAAGAAGAGTGTCCCCACATTTAGAAGATGTTCCCTCAGT
GAAGAGTGGCAGGCATTTCTTCAGATTTCTGAAAAAGAACTATGTAGTAAATGCCATTATCAAAA
AATCAAGATATCACAGCATGGCAAAGCCTGACACAGGTTCTTACCCAGAATCGTGGAGGAAAGCCAAC
ATAATGACCGAGCCCAGAACTCTCAGGGAAGATAAAGGGAATTTACATGGAAGAGAAATGTACAGA
CGTGCTCAGCATGATGACAGCCTCAGTTGGACCTCATGTGATCATCATGAGTCCCAAGAATGTAAGGA
GAGGACCTGGTAGACATCCCAGCTGTGGAAAAAATTGGGTATGAAATCAACGGTTGAAAAACGTAAT
GCGGCCATGTATTACCACAGCCTTCCCATGTAATAACTGTGGGGTGGCCTTTCAGATGATACAGAT
CCTCATGTCCATCACAGCACTACCTAGGAGAAAAATCTTATAAATGTGACCAGTATGAAAGAACCTT
AGTCAGAGCCAAGATCTTATCGTTCATTGTAAGAACTCACTCTGGCAAGACTCCCTATGAATTCACGAA
TGGCCTATGGGCTGCAAACAGAGCTCAGACCTTCCAGATATCAGAAAAGTCTCCTCAGGAGACAAACCC
TACAAATGTAAGAATGTGGCAAGGGCTTCAGGCGCAGCTCTCCCTTACAACCATCATCGAGTCCAC
ACAGGGGAGATGCCCTACAAATGCGATGAATGTGGAAAGGGTTTGGATTTAGGTCACTTCTTTGTATT
CATCAGGGAGTACACACAGGGAAAAAGCCCTATAAATGTGAAGAGTGTGGGAAGGGCTTTGATCAGAGC
TCCAACCTTCTTGTCATCAGAGAGTCCACACTGGAGAGAAGCCCTACAAATGCAGTGAAGTGGCAAG
TGCTTTAGTTCAAGCTCGTTCCTCAAGTCCACTGGAGGTTTACACAGGGGAGAAACCATATAGGTGT
GGTGAGTGTGAAAGGGCTTCAGCCAATGTACACACCTTACATTCACCAGAGAGTCCACACAGGGGAG
AAACCATACAAATGCAATGTGTGTGAAAGGATTTGCGTATAGCTCTGTTCTTCACTCATCAGAGA
GTTCACTGGAGAAAAACCATATAAATGTGAAGTGTGTGAAAGTGTTCAGTTACAGTTCATATTTT
CACTTACATCAAAGAGATCACATCAGAGAGAAACCATATAAATGTGATGAGTGTGGTAAAGGCTTCAGC
CGAATTCAGATCTTAATGTTACCTCAGAGTCCACACAGGAGAGAGGCCCTATAAGTGAAGGCATGT
GGTAAGGGCTTCAGTCGTAATTCATACCTCCTTGCCCATCAGAGAGTGCATATAGATGAGACACAGTAC
ACACACTGTGAGCGTGAAAGGACCTTCTGACTCATCAAAGACTACATGAGCAGAGAGAAACATTATAA
```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001291488

**Insert Size:** 1794 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001291488.1](#)

**RefSeq Size:** 2724 bp

**RefSeq ORF:** 1794 bp

**Locus ID:** 26974

**UniProt ID:** [Q96NJ3](#)

**Cytogenetics:** 19q13.31

**Protein Families:** Transcription Factors

**MW:** 68.9 kDa

**Gene Summary:**

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) encodes the longest isoform (a).