

## Product datasheet for **SC336939**

### PHF19 (NM\_001286840) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PHF19 (NM_001286840) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHF19
Synonyms:	MTF2L1; PCL3; TDRD19B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC336939 representing NM\_001286840.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGGTCTTGGTAATCCGTGGTCCCTATCCCAGCGCCCACTGTGTCAGGGGAAGCTGATGGAGAATCGA
GCTCTGGATCCAGGGACTCGGGACTCCTATGGTGCCACCAGCCACCTCCCAACAAGGGGGCCCTGGCG
AAGGTCAAGAACAACCTCAAGACTTGATGTCCAACTGACGGAGGGCCAGTATGTGCTGTGCCGTGG
ACAGATGGCCTGTACTACCTCGGGAAGATCAAGAGGGTCAGCAGCTCTAAGCAAAGCTGCCTCGTACT
TTCGAAGATAATTCAAAATACTGGGTCCTATGGAAGGACATACAGCATGCCGGTGTTCAGGAGAGGAG
CCCAAGTGCAACATCTGCCTAGGGAAGACATCAGGGCCGCTGAATGAGATCCTCATCTGCCGGAAGTGT
GGCCTGGGTTACCACCAGCAGTGCCACATCCCCATAGCGGGCAGTGCTGACCAGCCCTGCTCACACCT
TGGTCTGCCGACGCTGCATCTTCGCACTGGCTGTGCGGAAAGGGCGCGCTGAAGAAGGGCGCCATC
GCCAGGACGCTGCAGGCCGTGAAGATGGTGTCTCTACCAGCCCAGGAGCTCGAGTGGGACTCGCCC
CATCGCACCAACCAGCAGCAATGCTACTGCTACTGCGGGCGGCCGGAGAATGGTACCTGCGGATGCTG
CAATGTTACCGGTGCAGGCAAGTGTCCACGAGGCCTGCACCCAGTGCCTCAATGAGCCCATGATGTTT
GGAGACCGGTTTTACCTGTTCTTCTGCTCCGTGTGAACCAGGGCCAGAGTACATCGAGAGGCTGCC
CTGCGATGGGTGGATGTGGTTCACTGGCCCTCTATAATCTGGGGGTACAGAGCAAGAAGAAGTACTTT
GACTTTGAGGAGATTCTGGCCTTTGTCAACCACCCTGGGAGCTCCTGCAGCTTGGAAGCTCACCAGC
ACCCAGTGACAGATCGAGGACCACATCTCTCAACGCTCTGAACAGTTATAAAAGCCGGTTCCTCTGC
GGCAAGGAGATCAAGAAGAAGAAGTGCATCTCCGCCTGCGCATCCGCGTCCCACCAACCCGCCAGGG
AAGCTGCTGCCTGACAAAGGACTGCTGCCAAATGAGAACAGCGCCTCCTCTGAGCTGCGTAAGAGAGGA
AAGAGCAAGCCTGGTTTGTTCCTCACGAATCCAGCAGCAGAAAAGCGAGTTTATAGAAGAAAAGA
TCAAAGTTTTTGCTGGAAGATGCTATTCCAGTAGTGACTTACCTCAGCCTGGAGCACCAACCACCAC
CTGGCTAGCATATTTGACTTCACGCTGGATGAAATTCAAAGTTTAAAAAGTGCCAGCTCAGGCCAGACC
TTCTTCTCAGATGTCGACTCCACCGACGCTGCCAGCACCTCTGGCTCTGCCTCCACCAGCCTCTCTAT
GACTCCAGATGGACAGTGGGAGCCGAAAGAGGAAGCTGGCAGCCAAGGCATACATGCCCTGCGGGCA
AAGCGGTGGGAGCTGAGCTGGATGGACGCTGCCCTCGGACAGCAGTGCAGAGGGGGCTTCAGTCCCC
GAGCGGCCAGACGAAGGCATTGACAGCCACACATTTGAGAGCATCAGTGAAGATGACTCATCCCTGTCC
CACCTCAAGTATCTATACCAACTACTTTGGTGCAGCTGGGCGGTTGGCCTGTGGGAGAAAGTACCAG
GTGTTGGCTCGGAGGGTCACACCTGAGGGCAAGGTTACGTACCTGGTGGAGTGGGAAGGGACCACCCCT
TACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001286840

**Insert Size:** 1800 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001286840.1</a>
<b>RefSeq Size:</b>	4181 bp
<b>RefSeq ORF:</b>	1800 bp
<b>Locus ID:</b>	26147
<b>UniProt ID:</b>	<a href="#">Q5T6S3</a>
<b>Cytogenetics:</b>	9q33.2
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	67.6 kDa
<b>Gene Summary:</b>	<p>Polycomb group (PcG) that specifically binds histone H3 trimethylated at 'Lys-36' (H3K36me3) and recruits the PRC2 complex. Probably involved in the transition from an active state to a repressed state in embryonic stem cells: acts by binding to H3K36me3, a mark for transcriptional activation, and recruiting H3K36me3 histone demethylases RIOX1 or KDM2B, leading to demethylation of H3K36 and recruitment of the PRC2 complex that mediates H3K27me3 methylation, followed by de novo silencing. Recruits the PRC2 complex to CpG islands and contributes to embryonic stem cell self-renewal. Also binds dimethylated at 'Lys-36' (H3K36me2). Isoform 1 and isoform 2 inhibit transcription from an HSV-tk promoter. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (c).</p>