

Product datasheet for **SC336943**

RPS6KC1 (NM_001287220) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RPS6KC1 (NM_001287220) Human Untagged Clone
Tag:	Tag Free
Symbol:	RPS6KC1
Synonyms:	humS6PKh1; RPK118; RSKL1; S6K-delta-1; S6PKh1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC336943 representing NM_001287220.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTTAAAGCTCTGCCTTTGAAGAGTAGTCTTACTCCAAGTCTCAAGATGACAGCAACCAGGAAGAT
GATGGCCAAGATAGCTCTCCAAAGTGGCCAGATTCTGGTTCAAGTTCAGAAGAAGAAATGTACTACTAGT
TATTTAACATTATGCAATGAATATGGGCAAGAAAAGATTGAACCAGGGTCTTTGAATGAGGAGCCCTTC
ATGAAGACTGAAGGGAATGGTGTGATACAAAAGCTATTAAGAGCTTCCAGCACACCTTCTGCTGCTGAC
AGTGACAGCCCCAGCACAGCTGAGAGCTCACGAGCTGAAGTCTTCCCAACGATGACCCAGAAGCA
GTTAGTCTCCAAGAACATCAGATTCCTCAGTAGATCAAAAAATAGCCCATGGAATCTTTAGGATA
GACAGTAAGGATAGCGCAAGTGAACCTCTGGGACTTGACTTTGGAGAAAAATTGTATAGTCTAAATCA
GAACCTTTGAAACCATTCTTACTCTCCAGATGGAGACAGTCTTCTAGGAGTTTTAACTAGTGAA
AGCAAGGTAGAGTTTAAAGCTCAGGACACCATTAGCAGGGGCTCAGATGACTCAGTGCCAGTTATTTCA
TTTAAAGATGCTGCTTTTATGATGATGTCAGTGGTACTGATGAAGGAAGACCTGATCTTCTGTAAATTTA
CCTGGTGAATTGGAGTCAACAAGAGAAGCTGCAGCAATGGGACCTACTAAGTTTACACAACTAATATA
GGGATAATAGAAAATAAACTCTTGGAAAGCCCTGATGTTTTATGCCTCAGGCTTAGTACTGAACAATGC
CAAGCACATGAGGAGAAAGGCATAGAGGAACTGAGTGTCCCTCTGGGCCAAATCCTATAGTATAACA
GAGAAACACTATGCACAGGAGGATCCCAGGATGTTATTTGTAGCAGCTGTTGATCATAGTAGTTCAGGA
GATATGTCTTTGTTACCCAGCTCAGATCCTAAGTTTCAAGGACTTGGAGTGGTTGAGTCAGCAGTAACT
GCAAACAACACAGAAGAAAGCTTATTCGGTATTTGTAGTCCACTCAGGTGCTAATGAATATATTGCA
AGCACAGACACTTTAAAAACAGAAGAAGTATTGCTGTTTACAGATCAGACTGATGATTTGGCTAAAGAG
GAACCAACTTCTTTATTCAGAGAGACTCTGAGACTAAGGGTGAAGTGGTTTAGTGCTAGAAGGAGAC
AAGGAAATACATCAGATTTTGGAGACCTTGATAAAAAATTAGCACTAGCCTCCAGGTTTTACATCCCA
GAGGGCTGCATTCAAAGATGGGAGCTGAAATGGTGGTAGCCCTTGATGCTTTACATAGAGAGGGAATT
GTGTGCCGCGATTTGAACCCAAACAACATCTTATTGAATGATAGAGGACACATTCAGCTAACGATTTTT
AGCAGGTGGAGTGGGTTGAAGATTCCTGTGACAGCGATGCCATAGAGAGAATGTACTGTGCCCCAGAG
GTTGGAGCAATCACTGAAGAACTGAAGCCTGTGATTGGTGGAGTTTGGGTGCTGTCTCTTTGAACCT
CTCACTGGCAAGACTCTGGTTGAATGCCATCCAGCAGGAATAAATACTCACACTACTTTGAACATGCCA
GAATGTGTCTCTGAAGAGGCTCGCTCACTCATTCAACAGCTCTTGCAAGTCAATCCTCTGGAACGACTT
GGTGCTGGAGTTGCTGGTGTGAAGATATCAAATCTCATCCATTTTTACCCCTGTGATTGGGCAGAA
CTGATGAGATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001287220

Insert Size: 1806 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287220.1](#)

RefSeq Size: 4038 bp

RefSeq ORF: 1806 bp

Locus ID: 26750

Cytogenetics: 1q32.3

Protein Families: Druggable Genome, Protein Kinase

MW: 66.1 kDa

Gene Summary:

Sphingosine kinase catalyzes the formation of sphingosine 1 phosphate, a lipid cellular messenger. The protein encoded by this gene can bind to sphingosine kinase and to phosphatidylinositol 3-phosphate, suggesting a role in sphingosine 1 phosphate signaling. The encoded protein can also bind to peroxiredoxin-3 and may help transport it to mitochondria. [provided by RefSeq, Mar 2017]

Transcript Variant: This variant (5) differs in its 5' UTR and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (d) has a shorter N-terminus, compared to isoform a. Variants 5 and 22-30 all encode the same isoform (d).

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.