

Product datasheet for **SC336954**

ZNF169 (NM_003448) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF169 (NM_003448) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF169
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336954 representing NM_003448.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCACCAAGGACTCCTGACAACCAAGGAGGACATTGATGGCCTTCCGGGATGTGGCTGTGGCCTTC
ACCCAGAAGGAGTGGAAGCTATTGAGTTCTGCTCAGAGGACCTGTACAGGGAGGTGATGCTGGAGAAC
TACAGCCATCTGGTCTCCCTGGGAATTGCATTTTCCAAACCAAACTCATCGAACAGCTGGAGCAAGGC
GACGAACCTTGGAGAGAGGAGAACGAACATCTTCTGGACCTTTGTCCAGCAGAGCCTAGAACAGAATTC
CAGCCAAGTTTTCCACCTGGTGGCCTTTCTAGCTCGCAGCTCCTCAGACAATATGCGTAAGTGGC
CATCCACACAGATCTTCCCAAGCTCATCTGCAGGAGGTGACTTCCAAGTGAAGCTCCAAGATGCTCT
AGTGAAGGAGGAGAAAGTGAGAGACAGAAGGCCCGACAGCTCATTAGAAAGAGGCCAAGCAGAATT
TCTAGGACATTCTCAGCCACATCAAGGTGACCCAGTAGAATGGGTAGAAGGGAACAGAGAAGGAGGA
ACAGACCTTCGCTGGCCAAAGGATGAGTCTTGGGGGTGAGACACAATGTTGAAGGGAGCAGACACT
TCAGAATCTGGAGCAGTCATACGTGGAACATATAGACTGGGACTTAGCAAAAAGTCAAGCCTGTTACGC
CACCAGAAGCATCATGTGTGCCCTGAATGCGGGAGAGGCTTTGCCAGAGATCAGACCTTATCAAGCAC
CAGAGGACACACACCGGGGAGAAGCCATACCTGTGTCTGAGTGTGGCGCTCGGTTTAGCCAGAAGGCC
TCCCTCTCCATACACCAGAGGAAGCACTCGGGGAGAAGCCGTATGTGTGCAGGGAATGTGGGCGACAC
TTCAGGTATACATCTCTCTCACTAATCACAAGAGGATTCACTCCGGGAGAGGCCCTTTGTATGTCAG
GAGTGTGGGCGAGGCTTTCGCCAGAAGATAGCCCTCCTTCTACACCAGAGGACGCACTTGAGGAGAAG
CCCTTCGTGTGCTCTGAGTGTGGGAGAGGCTTTGCCAGAAGGCATCACTCCTCCAGCACCAGAGCTCA
CACACAGGGGAGAGGCCCTTCTGTGCCTTGTGTGGCGTAGCTTCAGGCAGCAGTCACTCCTCCTT
AGTCACCAAGTCAACACTCAGGAGAGAAGCCTTATGTCTGTGCTGAGTGTGGGCACAGCTTTCGCCAA
AAGGTCACTCTCATCAGGCACCAGAGGACACACAGGGGAGAAGCCTTACCTGTGCCCCAGTGTGGG
CGGGGTTTTAGCCAGAAGGTCAACCTCATTGGACACCAGAGGACACACAGGGGAGAAGCCCTACCTG
TGCCCTGATTGTGGGCGTGGCTTTGGTCAGAAGGTCAACCTCATCAGACACCAGAGGACACACAGGG
GAGAAGCCTTATCTGTGCCCAAGTGTGGGCGTGATTTGGCTTTAAGTCGCTCCTACCCGACACCAG
AGGACACACTCAGAGGAGGAGCTTTACGTAGACAGGTGTGTGGACAAGGACTTGCCAGAAGTCACAC
CTTATCTCTGACCAAGGACACACTCAGGAGAGAAGCCCTGCATTTGCGATGAATGTGGGCGCGGCTTT
GGCTTTAAGTCTGCCCTCATCCGACATCAGCGGACCCATTCTGGGAGAAGCCGTATGTCTGACGGGAG
TGTGGGCGTGGCTTTAGCCAGAAGTCTCACTTGCATAGACACAGGAGGACCAAGTCTGGTCATCAGCTC
CTACCCCAAGAGGTCTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_003448

Insert Size: 1815 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_003448.2</u>
RefSeq Size:	2409 bp
RefSeq ORF:	1815 bp
Locus ID:	169841
Cytogenetics:	9q22.32
Protein Families:	Transcription Factors
MW:	68.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.