

Product datasheet for **SC336984**

TTC39A (NM_001297665) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39A (NM_001297665) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39A
Synonyms:	C1orf34; DEME-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC336984 representing NM_001297665.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGCCAGAAGGGCCACAAAGACTCCCTTTACCTTTGTGGAGGGACTCCTGAGAGCAGCCTCCATGAG
GCCCTGGACCAGTGCATGACCGCCTGGACCTTCTCTCACCAACCAGTTCTCAGAAGCACTCAGCTAC
CTCAAGCCCAGAACCAAGGAAAGCATGTACCACTCACTGACATATGCCACCATCCTGGAGATGCAGGCC
ATGATGACCTTTGACCCTCAGGACATCCTGCTTCCGGCAACATGATGAAGGAGGCACAGATGCTGTGT
CAGAGGCACCGGAGGAAGTCTTCTGTAACAGATTCTTCAGCAGCCTGGTGAACCGCCCCACGCTGGGC
CAATTCCTGAAGAGGAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTG
ACCTTCTGCAGGGTTCTCACACGGAGGGGCAGTCAGGCCCAGAGCCTTGCATGATCCCTCTCACGCC
TGCAGTGCACCTGGGCCAGGCCGTGACATCTTTCTCTGCAGGACGAGAACATGGTGAGCTTC
ATCAAAGCGGCATCAAAGTTCGAAACAGCTACCAGACCTACAAGGAGCTGGACAGCCTTGTTCAGTCC
TACAATACTGCAAGGGTGAAGAACCCCGCACCTTGAAGGAGGAGTGAAGCTTGGTGTAGGGGCCCTTC
AACCTGCACTGTCCATGCTTCTACTAGGATCCTGAGGCTGTTGGAGTTTGTGGGGTTTTAGGAAAC
AAGGACTATGGGCTGCTGCAGCTGGAGGAGGGAGCGTCAGGGCACAGCTCCGCTCTGTGCTCTGTGTC
ATGCTCCTGCTGTGCTACCACACCTTCCACCTTCGTGCTCGGTACTGGGAACGTAACATCGAGGAG
GCCGAGAAGCTCTTGAAGCCCTACCTGAACCGGTACCCTAAGGGTGCCATCTTCTGTTCTTTCAGGG
AGGATTGAAGTCATTAAGGCAACATTGATGCAGCCATCCGGCGTTTTCGAGGAGTGTGTGAGGCCAG
CAGCACTGGAAGCAGTCCACCACATGTGCTACTGGGAGCTGATGTGGTGTCTCACCTACAAGGGCCAG
TGGAAGATGTCTACTTCTACGCCGACCTGCTCAGCAAGGAGAACTGCTGGTCCAAGGCCACCTACATT
TACATGAAGGCCGCTACCTCAGCATGTTTGGGAAGGAGGACCACAAGCCGTTTCGGGGACGACGAAGTG
GAATTATTTTCGAGCTGTGCCAGGCTGAAGCTCAAGATTGCTGGGAAATCTCTACCCACAGAGAAGTTT
GCCATCCGGAAGTCCCGGCGCTACTTCTCTCCAACCCTATCTCGCTGCCAGTGCCTGCTCTGAAATG
ATGTACATCTGGAACGGCTACGCCGTGATTGGGAAGCAGCCGAAACTCACGGATGGGATACTTGAGATT
ATCACTAAGGCTGAAGAGATGCTGGAGAAAGGCCAGAGAACGAGTACTCAGTGGATGACGAGTGTG
GTGAAATTGTTGAAAGGCTGTGCTGAAATACCTGGGCCGTGTCCAGGAGGCCGAGGAGAAATTTAGG
AGCATCTCTGCCAATGAAAAGAAGATTAATATGACCACTACTTGATCCCAAACGCCCTGCTGGAGCTG
GCCCTGCTGCTTATGGAGCAAGACAGAAACGAAGAGGCCATCAAATTTTGAATCTGCCAAGAAAAAC
TACAAGAATTACTCCATGGAGTCAAGGACACACTTTTGAATCCAGGCAGCCACACTCCAAGCCAAGTCT
TCCTAGAGAACAGCAGCAGATCCATGGTCTCATCAGTGTCTTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001297665

Insert Size: 1842 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297665.1](#)

RefSeq Size: 2866 bp

RefSeq ORF: 1842 bp

Locus ID: 22996

UniProt ID: [Q5SRH9](#)

Cytogenetics: 1p32.3

MW: 69.8 kDa