

Product datasheet for SC337001

DGLUCY (NM_001286471) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DGLUCY (NM_001286471) Human Untagged Clone
Tag:	Tag Free
Symbol:	DGLUCY
Synonyms:	C14orf159
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC337001 representing NM_001286471. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCCTCACACTCCACCTGAGGTCCCGCCTCCCTCTGCCATAAGGAGTTTGATTCTACAAAAGAAA
CCAAACATCAGAAATACATCCAGCATGGCTGGAGAGCTCCGACCAGCCAGCCTGGTGGTCTGCCAGG
TCCCTTGCTCCAGCTTTTGAAGATTCTGCCAGGTCAACTGGTCTCTACCCCTGCTGGGCCAGAGT
GAGCCAGAAAAGTGGATGCTGCCCCCTCAAGGTGCTATCTCAGAGACCAGGATGGGCCATCCCCAGTTC
TGGAAATACGAGTTCGGTGCCTGCACCCGACGCTGGCTTCGCTGGAGCAGTACTCGGAGCAGTGAAG
GACATGGTGGCCTTCTTCTGGGCTGCAGCTTCTCCCTGGAGGAGGCCTGGAGAAGCGGGGCTCCCC
AGAAGAGACCCAGCAGGTACAGCCAGACAACAGTGCCTTGTTGTTACCCATGCTGGCTTCTGCTGCCCT
CTGGTGGTCACGATGAGGCCATTCCCAAGGACAAGCTGGAAGGGCTGGTGCAGGCTGCTGCCCTC
GGAGGTGAGCAGGGGCAACCTGTTACATGGGCGACCCAGAAGTGTGGGAATCAAAGAGCTTTCCAAA
CCTGCCTACGGGGATGCCATGGTGTGTCCCCAGGGGAGGTTCCAGTGTCTGGCCTTCTCCGCTGACC
AGTCTCGGAGCTGTCAGCAGCTGTGAGACCCACTGGCTTTTCCAGCATCCCAGGCTGCACAGTTATG
ACTGACCTGAAGGATGCAAAGGCTCCACCTGGTTGTCTCACCCAGAGAGAATCCAGAGGTCCATCAC
ATTTCCAAGATCCTCTGCACTACAGCATCGCGTCAGTCTGCTTCTCAGAAGATCAGAGAAGTAGAG
TCTATGATCGGCATAGACCCAGGGAACCGGGGATTGGGCACCTGCTCTGTAAGATGAGCTGTAAG
GCCTCTCTCGCTGTCCATGCCGCTCAGTGCATCACCACTGGGTTCCCCACACATTTCAATCAT
GAGCCTCCAGAAGAGACAGATGGCCACCAGGAGCTGTTGCTCTGGTTGCCTTCTGCAGCCTTGAG
AAGGAGTGCATAATCGTTGACCAGAGAGCCTGGAAGTGCACCAGAAGATTGTTGAAGATGCTGTT
GAGCAAGGTGTTCTGAAGACGCAGATCCCGATTAATACTACCAAGGTGGATCAGTGAAGCTGCTCAG
GCATTCTGTGCAAAAATGGGGACCCGACACCTAGATTTGACCACCTGGTGGCCATAGAGCGTGCC
GGAAGAGCTGCTGATGGCAATTACTACAATGCAAGGAAGATGAACATCAAGCACTGGTTGACCCATT
GACGATCTTTTCTGCTGCGAAGAAGATTCTGGAATCTCATCACTGGAGTGGTGGATGGAGGCAAC
GAGCTTGGGATGGTAAAGTCAAGGAGGCTGTGAGGAGGCACATACGGCACGGGGATGTCATCGCCTGC
GACGTGGAGGCTGACTTTGCCGTATTGCTGGTGTCTAACTGGGAGGCTATGCCCTGGCCTGGCA
CTCTACATCCTGTACTCATGTGCTGTCCACAGTACCTGAGGAAAGCAGTCCGACCCCTCCAGGGCA
CCTGGAGATCAGGCTGGACTCAGGCCCTCCCGTCCGTCATTAAGGAAGAAAAATGCTGGGCATCTTG
GTGACGACAAAAGTCCGGAGTGGCGTCTCGGGCATCGTGGGATGGAGGTGGATGGGCTGCCCTCCAC
AACACCCACGCCGAGATGATCCAGAAGCTGGTGGACGTACCACGGCACAGGTGTA


```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001286471

Insert Size: 1851 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286471.1](#)

RefSeq Size: 3017 bp

RefSeq ORF: 1851 bp

Locus ID: 80017

UniProt ID: [Q7Z3D6](#)

Cytogenetics: 14q32.11

MW: 66.4 kDa

Gene Summary:

D-glutamate cyclase that converts D-glutamate to 5-oxo-D-proline.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (7) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, and 7 encode the same isoform (a).