

## Product datasheet for **SC337005**

### ZNF135 (NM\_001289402) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF135 (NM_001289402) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF135
Synonyms:	pHZ-17; pT3; ZNF61; ZNF78L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC337005 representing NM\_001289402.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGGACACCTTCAGGCTTCTGGTCTCTGTGGGACATTGGTTACCGAAGCCGAATGTCATCTCCCTG
CTGGAGCAAGAGGCAGAGCTGTGGCGGTGGAGTCTAGACTTCCCAAGGCGTGTACCCAGACTTGGAA
ACTAGACCCAAAGTCAAAGTGTAGTTCTAAAGCAAGGCATCTCTGAAGAAATATCCAACAGTGTATC
TTGGTAGAAAGATTCTGTGGGATGGTCTGTGGTACTGCAGGGGTGAGGACACTGAGGGCCACTGGGAA
TGGAGTTGTGAGAGTCTAGAGAGCCTGGCAGTGGCGGTGGCCTTACGCCTGTGAAGACGCCTGTTCTG
GAGCAGTGGCAGAGGAATGGGTTTGGGAAAACATAAGTCTGAACCTGATCTCCACATCAACCAATG
ACTCTGAAAGACAAAGCCCCACACATGGGGAACACGTGAAAAAGGGAGAAGCCAGACCTAAATGTT
TTACAGAAAACCTGTGTAAGAGAAAACCTACAATGTCAGGAATGCGGAAAGGCCTTTAGTCACAGC
TCAGCACTTATCGAACACCACCGGACGCACACAGGAGAGACCTTACGAATGTCACGAATGCTTAAAA
GGCTTCCGGAACAGCTCGGCACCTACCAAACACCAGAGAATCCACTGGGGAGAAAACCTATAAATGC
ACTCAGTGTGGGAGGACCTTCAACCAAATGCCCCACTGATCCAGCACAGAGAAGCTCACACAGGTGAG
AAGCCCTATGAATGCAGCGAATGTGGGAAATCCTTCAGTTTTAGGTCCTCCTTCAGCCAGCACGAGCGA
ACTCACACAGGCGAGAAGCCCTACGAGTGCAGTGAAGTGTGGGAAAGCCTTCCGGCAAAGCATCCACCTC
ACCCAGCATCTGCGAATCCACACTGGGGAGAAACCTATCAGTGTGGTGAAGTGTGGCAAGGCCTTCAGC
CACAGCTCATCCTTGACCAAACACCAGCGAATCCACACAGGGGAGAAGCCCTACGAGTGCATGAGTGT
GGAAAAGCCTTACCCAGATCACACCACTGATTCAGCACCAGAGGACCCACACAGGAGAAAAGCCCTAT
GAGTGTGGTGAAGTGTGGGAAAGCCTTCACTCAGAGCACACTCCTGACCGAGCATCGGAGGATTACACA
GGAGAGAAGCCCTATGGATGCAACGAGTGTGGGAAAACCTTCAAGCACAGCTCCTCACTCAGCCAGCAT
GAGCGGACACACAGGAGAGAAGCCCTATGAGTGCAGTCAAGTGTGGGAAAGCCTTCCGGCAGAGCACA
CACCTCACCAAACACCAGCGAATCCACACAGGGGAGAAGCCCTATGAATGCAATGACTGCGGCAAGGCA
TTCAGTACAGCTCGTCCCTCACAAACATCAGCGAATCCACACTGGGGAGAAGCCCTACGAATGCAAC
CAGTGTGGCAGAGCCTTACGCCAGCTTGTCCCCTCATTACAGATCAGAGGATCCACACAGGAGAGAAA
CCCTATGAATGTAACAGTGTGGCAGAGCCTTACGCCAGAGCTCCCTTCTCATCGAACACCAGAGGATT
CACACCAAGGAAAAGCCGATGGGTGCAATGAGTGTGGGAAATCCTTACGCCACAGCTCCTCGCTCAGC
CAGCACGAAAGGACGCACACTGGGAAAAGCCCTATGAGTGTACGATTGCGGAAAGTCTTTAGGCAG
AGCACCCACCTCACTCAGCACCGGAGGATCCACACAGGAGAGAAGCCATATGCATGCAGGGACTGTGGA
AAGGCCTTTACCCACAGCTCCTCCCTTACCAAGCACAGAGAAGCTCACACTGGATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001289402

**Insert Size:** 1851 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001289402.1](#)

**RefSeq Size:** 3740 bp

**RefSeq ORF:** 1851 bp

**Locus ID:** 7694

**UniProt ID:** [P52742](#)

**Cytogenetics:** 19q13.43

**Protein Families:** Transcription Factors

**MW:** 70.5 kDa

**Gene Summary:**

Plays a role in the regulation of cell morphology and cytoskeletal organization. May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (7) uses an alternate in-frame splice site in the 3' coding region, and includes an alternate exon that results in a distinct 5' UTR and uses of a downstream start codon, compared to variant 1. The encoded isoform (7) is shorter, compared to isoform 1.