

Product datasheet for **SC337006**

EIF4B (NM_001300821) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EIF4B (NM_001300821) Human Untagged Clone
Tag:	Tag Free
Symbol:	EIF4B
Synonyms:	EIF-4B; PRO1843
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337006 representing NM_001300821.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGGCCTCAGCAAAAAGAAGAATAAGAAGGGGAAGACTATCTCCCTAACAGACTTTCTGGCTGAG
GATGGGGTACTGGTGGAGGAAGCACCTATGTTTCCAAACCAGTCAGCTGGGCTGATGAAACGGATGAC
CTGGAAGGAGATGTTTCGACCACTTGGCACAGTAACGATGACGATGTGTATAGGGCGCCTCCAATTGAC
CGTTCCATCCTTCCCACTGCTCCACGGGCTGCTCGGGAACCCAATATCGACCGGAGCCGCTTCCCAAA
TCGCCACCCTACACTGCTTTTCTAGGAAACCTACCCTATGATGTTACAGAAGAGTCAATTAAGGAATTC
TTTCGAGGATTAATATCAGTGCAGTGCCTTACCACGTGAACCCAGCAATCCAGAGAGGTTGAAAGGT
TTTGTTATGCTGAATTTGAGGACCTGGATTCCCTGCTCAGTGCCTGAGTCTCAATGAAGAGTCTCTA
GGTAACAGGAGAATTCGAGTGGACGTTGCTGATCAAGCACAGGATAAAGACAGGGATGATCGTCTTTT
GGCCGTGATAGAAATCGGGATTCTGACAAAACAGATACAGACTGGAGGGCTCGTCTGCTACAGACAGC
TTTGATGACTACCCACCTAGAAGAGGTGATGATAGCTTTGGAGACAAGTATCGAGATCGTTATGATTCA
GACCGGTATCGGGATGGGTATCGGGATGGGTATCGGGATGGCCACGCCGGGATATGGATCGATATGGT
GGCCGGGATCGCTATGATGACCGAGGCAGCAGAGACTATGATAGAGGCTATGATCCCGGATAGGCAGT
GGCAGAAGAGCATTTGGCAGTGGGTATCGCAGGGATGATGACTACAGAGGAGCGGGGACCCTATGAA
GACCGATATGACAGACGGGATGATCGGTGCTGGAGCTCCAGAGATGATTACTCTCGGGATGATTATAGG
CGTGATGATAGAGGTCCCCCAAGACCCAACTGAATCTAAAGCCTCGGAGTACTCCTAAGGAAGAT
GATTCCTCTGCTAGTACCTCCAGTCCACTCGAGCTGCTTCTATCTTTGGAGGGCAAAGCCTGTTGAC
ACAGCTGCTAGAGAAAGAGAAGTAGAAGAACGGCTACAGAAGGAACAAGAGAAGTTGCAGCGTCAGCTG
GATGAGCCAAAACCTAGAACGACGGCTCGGGAGAGACACCAAGCTGGCGAAGTGAAGAACTCAGGAA
CGGGAACGGTCGAGGACAGGAAGTGAGTCATCACAACTGGGACCTCCACCACATCTAGCAGAAGTAAG
TCAGACCAGGATGCACGAAGGAGAGAGAGTGAGAAGTCTCTAGAAAATGAAACACTCAATAAGGAGGAA
GATTGCCACTCTCAACTTCTAAACCTCCCAAACCTGATCAGCCCTAAAGGTAATGCCAGCCCCTCCA
CCAAAGGAGAAATGCTTGGGTGAAGCGAAGTTCTAACCTCCTGCTCGATCTCAGAGCTCAGACACAGAG
CAGCAGTCCCCTACAAGTGGTGGGGAAAAGTAGCTCCAGCTCAACCATCTGAGGAAGGACCAGGAAGG
AAAGATGAAAATAAAGTAGATGGGATGAATGCCCAAAAGGCCAACTGGGAACCTAGCCGTGGTCCA
GGAGACGGAGGGAACAGAGACCACTGGAAGGAGTCAGATAGGAAAGATGGCAAAAAGGATCAAGACTCC
AGATCTGCACCTGAGCCAAAGAAACCTGAGGAAAATCCAGCTTCCAAGTTCAGTTCTGCAAGCAAGTAT
GCTGCTCTCTGTTGATGGTGAAGATGAAAATGAGGGAGAAGATTATGCCGAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001300821

Insert Size: 1851 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001300821.1](#)

RefSeq Size: 4076 bp

RefSeq ORF: 1851 bp

Locus ID: 1975

Cytogenetics: 12q13.13

Protein Pathways: mTOR signaling pathway

MW: 69.7 kDa

Gene Summary:

Required for the binding of mRNA to ribosomes. Functions in close association with EIF4-F and EIF4-A. Binds near the 5'-terminal cap of mRNA in presence of EIF-4F and ATP. Promotes the ATPase activity and the ATP-dependent RNA unwinding activity of both EIF4-A and EIF4-F. [UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.