

Product datasheet for **SC337015**

ZNF539 (ZNF254) (NM_001278677) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF539 (ZNF254) (NM_001278677) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF254
Synonyms:	BMZF-5; HD-ZNF1; ZNF91L; ZNF539
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337015 representing NM_001278677.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTAGAGAACTACAGAAACCTGGCCTTCCTGGGTATTGCTGTCTCTAAGCCAGACCTGATCACCTGT
CTGGAACAAGGGAAAGACCCCTGGAATATGAAGCGACATGAGATGGTGGATGAACCCCAAGTATGTGT
CCTCATTTTGTCAAGACCTTTGGCCAGAGCAGGGCATGGAAGATTCTTTTCAAAAAGCAATACTGAGA
AGATATGGAAAATATGGACATGAGAATTTACAGTTAAGAAAAGGCTGTAAGAGTGTGGATGAGTATAAG
GTGAACAAAGAAGGTTATAATGGACTTAACCAGTGTTCACAACCTGCCAGAGCAAAGTATTTCAATGT
GATAAATATTTGAAAGTCTTCTATAAATTTTAAATTCAAACAGACCTAAGATAAGACATACTGAAAAG
AAATCTTTCAAATGTAAAAACGTGTCAAATTATTTTGCATGCTTTCACATAAAACCCAACAAAAAGC
ATTTATCATAGAGAGAAGTCTACAAATGTAAAGAATGTGAAAAACCTTAATTGGTCTCAACCCTT
ACTAATCATAGGAAAATTTATACTGAAGAGAAACCTTACAAATGTGAAGAATATAACAAATCTCCTAAG
CAACTCTCAACCCTTACTACACATGAAATAATTCATGCTGGAGAGAACTCTACAAATGTGAAGAATGT
GGTGAAGCTTTTAAATCGATCCTCAAATCTTACTACACATAAGATAATTCATACTGGAGAGAAACCTTAC
AAGTGTGAAGAATGTGGCAAAGCATTATCTGGTCTCAACCCTTACTGAGCATAAGAAAATTCATACT
AGAAAAGAAACCTACAAGTGTGAAGAATGTGGCAAAGCATTATATGGTCTCAACACTAAGTAGACAT
AAGAGGATGCACACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCTTTAGCCAGTCTCA
ACCCTTACTACACATAAGATAATTCATACTGGAGAGAAACGCTACAAATGCTTAGAATGTGGCAAAGCT
TTAAGCAACTCTCAACTCTTACTACACATAAAATAATTCATGTTGGAGAGAACTCTACAAATGTGAA
GAATGCGGCAAAGGTTTTAATCGATCTTCAAATCTTACTACACATAAGATAATTCATACTGGAGAGAAA
CCTTACAAGTGTGAAGAATGTGGCAAAGCATTATCTGGTCTCAACCCTTACTAAACATAAGAGAATT
CATACTAGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCATTATATGGTCTCAACCCTAACT
AGACATAAGAGGATGCACACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAATCTTTTAGCCAA
TCCTCAACCCTTACTACACATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAGAATGTGGC
AAAGCCTTTAACTGGTCTCAACTCTTACTAAACATAAGATAATTCATACTGAAGAGAAACCTACAAA
TGTGAAAATGTGGCAAAGCCTTAAAGCAGTCTTCAATCCTTACTAACCATAAGAGAATTCATACTGGA
GAGAAACCTATAAATGTGAAGAATGTGGCAAATCTTTAACCGTCTTCAACTTTTACTAAACATAAG
GTAATTCATACTGGAGTAAAACCTACAAATGTGAAGAATGTGGCAAAGCATTTTCTGGTCTCAACC
CTAACTAAACATAAGAGAATTCATACTGGAGAGCAACCCTACAAATGGGAAAAATTTGGCAAAGCCTT
AATCGGTCTCGCACCTACCACAGATAAGATAACTATTGGAGAGAAATCTTACAAGTATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001278677

Insert Size: 1857 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278677.1](#)

RefSeq Size: 4334 bp

RefSeq ORF: 1857 bp

Locus ID: 9534

Cytogenetics: 19p12

Protein Families: Transcription Factors

MW: 72.4 kDa

Gene Summary:

Zinc finger proteins have been shown to interact with nucleic acids and to have diverse functions. The zinc finger domain is a conserved amino acid sequence motif containing 2 specifically positioned cysteines and 2 histidines that are involved in coordinating zinc. Kruppel-related proteins form 1 family of zinc finger proteins. See ZFP93 (MIM 604749) for additional information on zinc finger proteins.[supplied by OMIM, Jul 2002]

Transcript Variant: This variant (1) represents the longest transcript and encodes isoform a. Variants 1, 2, and 3 encode isoform a.