

Product datasheet for **SC337018**

SLC1A7 (NM_001287595) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC1A7 (NM_001287595) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC1A7
Synonyms:	AAAT; EAAT5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337018 representing NM_001287595.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGCCGCATGCCATCTTGGCACGGGGAGGGACGTGTGCAGGCGGAATGGACTCCTCATCCTGTCT
GTGCTGTCTGTGCATCGTGGGCTGCCTCCTCGGCTTCTTCTTGAGGACCCGGCGCCTCTACCACAGGAA
ATTAGTTACTTCCAGTTCCTGGAGAGCTCCTGATGAGGATGCTGAAGATGATGATCCTGCCACTGGTG
GTCTCCAGCTTGATGTCCGGACTTGCTCCCTGGATGCCAAGACCTCTAGCCGCCTGGGCGTCTCACC
GTGGCGTACTACCTGTGGACCACCTTCATGGTGTGCATCGTGGGCATCTTCATGGTCTCCATCATCCAC
CCAGGCAGCGCGGCCAGAAGGAGACCACGGAGCAGAGTGGGAAGCCATCATGAGCTCAGCCGATGCC
CTGTTGGACCTCATCCGGAACATGTTCCAGCCAACTAGTAGAAGCCACATTCAAACAGTACCGCACC
AAGACCACCCAGTTGTCAAGTCCCCAAGGTGGCACCAGAGGAGGCCCTCCTCGGCGGATCCTCATC
TACGGGGTCCAGGAGGAGAATGGCTCCCATGTGCAGAACTTCGCCCTGGACCTGACCCCGCCCGGAG
GTCGTTTACAAGTCAGAGCCGGGCACCAGCGATGGCATGAATGTGCTGGGCATCGTCTTCTCTGCC
ACCATGGGCATCATGCTGGGCGCATGGGTGACAGCGGGGCCCCCTGGTCAGCTTCTGCCAGTGCCTC
AATGAGTCGGTCAAGATCGTGGCGGTGGCTGTGTGGTATTTCCCTTCGGCATTGTGTTCTCATT
GCGGGTAAAGTCCCTGGAGATGGACGACCCAGGGCCGTGCGGAAGAAGCTGGGCTTCTACTCAGTCACC
GTGGTGTGCGGGCTGGTCTCCACGGGCTCTTTATCCTGCCCTGCTCTACTTCTTCATACCAAGAAG
AATCCCATCGTCTTCATCCGTGGCATCCTGCAGGCTCTGCTCATCGCGCTGGCCACCTCCTCCAGGAAG
AACCAAGTGAGAGCACTGCAGTGCCCGCGGTCTTACCTCCTGCCGAGCAGCGCTGCAGAGGCACAAA
CACTCCTGGCGGGTGGTGCCTGGGGAGGGGTGAGCTTCCCGACAGAGGCCCCAGCCATCCTGGAGGA
GGGTGCTGGGCAGACACCGGCCCTCTAACCGTAGCTCAGCCACACTGCCATCACCTTCAAGTGCCTG
CTGGAGAACAACCACATCGACCGGCGCATCGCTCGCTTCTGCTGCCCGTGGGTGCCACCATCAACATG
GACGGCACTGCGCTCTACGAGGCTGTGGCCGCCATCTTCATCGCCAGGTCAACAACACTACGAGCTGGAC
TTTGGCCAGATCATCACCATCAGTATCACAGCCACTGCAGCCAGCATTGGGGCAGCTGGCATCCCCCAG
GCCGGCTCGTCACCATGGTTCATCGTGTACCTCCGTGGGACTGCCACCGATGACATCACCCTCATC
ATTGCCGTTGACTGGGCTCTGGACCGTTTCCGCACCATGATTAACGTGCTGGGTGATGCGCTGGCAGCG
GGGATCATGGCCATATATGTCGGAAGGATTTGCCCGGGACACAGGCACCGAGAACTGCTGCCCTGC
GAGACCAAGCCAGTGAGCCTCCAGGAGATCGTGGCAGCCAGCAGAATGGCTGTGTGAAGAGTGTAGCC
GAGGCCTCCGAGCTCACCTGGGCCCCACCTGCCCCACCACGTCCCGTTCAAGTGGAGCAGGATGAG
GAGCTGCCCGCTGCGAGTCTGAACCACTGCACCATCCAGATCAGTGAGCTGGAGACCAATGCTGA
ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001287595

Insert Size: 1860 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287595.1](#)

RefSeq Size: 2873 bp

RefSeq ORF: 1860 bp

Locus ID: 6512

Cytogenetics: 1p32.3

Protein Families: Transmembrane

MW: 66.9 kDa

Gene Summary:

Transports L-glutamate; the L-glutamate uptake is sodium- and voltage-dependent and chloride-independent. Its associated chloride conductance may participate in visual processing.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.