

## Product datasheet for **SC337042**

### **ANTXRL (NM\_001278688) Human Untagged Clone**

#### **Product data:**

|                           |                                            |
|---------------------------|--------------------------------------------|
| Product Type:             | Expression Plasmids                        |
| Product Name:             | ANTXRL (NM_001278688) Human Untagged Clone |
| Tag:                      | Tag Free                                   |
| Symbol:                   | ANTXRL                                     |
| Mammalian Cell Selection: | Neomycin                                   |
| Vector:                   | pCMV6-Entry (PS100001)                     |
| E. coli Selection:        | Kanamycin (25 ug/mL)                       |



[View online »](#)

**Fully Sequenced ORF:** >SC337042 representing NM\_001278688.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGGAGCCATGAGTCCCTGGGGCCCTACTTCCTGGTCTTCTGCTGCTGCTGCTCTCTCCACCG
CTTTTATAGACAGGAAGCCTTCGGTACCATGGACCTGACTGGAGAATATTCACCGCCTGGCCCTGGGC
TCCAGGAGAGCCACCACCACCATGGCCAGGATGGAGGCAGCACTGGCGCCAGGGGCAAGCAGGTCAC
AGATGCCAGGGCTCATTTGACCTCTACTTCATCTTGACAAGTCTGGCAGCGTGAACAATAACTGGATT
GACCTTTATATGTGGTGGAGGAAACAGTGGCGAGGTTCCAAAGCCAAATATTCGGATGTGCTTCATC
ACCTACTCCACAGACGGCCAGACTGTCTTGCCACTCACCTCAGACAAGAATAGAATAAAAAACGGTCTT
GACCAACTTCAGAAAATTGTGCTGACGGTCACACATTTCATGCAGGCAGGATTTAGAAAGCAATTCAA
CAGATCGAAAGTTTCACTCCGGAACAAGGTTCCAGCATGATTATTGCTATGACTGATGGAGAACTG
GTGGCAGATGCATTTAGGACACTCTCAGAGAAGCTCAAAGGCTCGGAACTGGGGGCCAACGTTTAC
ACCCTGGGTGTGGCTGATTATAATCTGGATCAGATAACAGCAATTGCAGACAGCCCTGGCCACGTGTTT
GCAGTGGAGAATGGCTTCAAGGCCCTGAGAAGCACCATTTGATGCCCTCACGTCAAAGGTCTGTCTTGAT
GTGACATCGGTGGAGCCTTCTCTGAGTGTGTAGGAGAACCTACCATGTGGTTATTCATGGAAATGGC
TTTCAGAACTAAAGAAACGGGATGAAGTTATTTGAGATTTATCTTCAATGAAAGCACTATCATTGAT
GAAAAGCCAACAGTATCGACAATAATTCATGAATTGCCCTGGGCCAAAACCTAGAAAAACCTGGAGAG
GAGTACTCTATTGAAGTCAGCTTGAACAAAGGCAAAACATTCTTCAAGAGCAATGTCAGCATCACCAGC
ACCACATGTGGCATTTTCCGCAACTGGCTCTATTTGTGCCACTCCTGCTGCTTGTGCCACTGCTGCTG
TGTTGTGTCTGGCGCTGTGCCGCAAGCAGACTGTCAAGGAGCCACCACCTGTGCAGAAGCCAGAAAAG
GAGCCAGAGCAGGAAAAACCACCATCACCAACCACCGCCTCCGCCTCCACCACCTCCACTCCCGCCT
CCGCCCCCAGCTCCTGTAAACACCTGCCCACTGTGATTATTTGTTGCTGTGGATGCCAAGGAGTGGGC
GGGATGAGAAGGATAGAGGGCAATCTGGATACCTTTTGTGACCTCTCTCACGCAAGCTGCCACCAGGTG
CCATGGATGTGTGTGAGAGCAGGGACAGGGGAGGTACCTCAGCTTAGCCCTTGACAGTCCCAATAT
GCACAGGCTCCCTGCTGCCCAAGGATCTGCTTTCCACACAGCCAGGAGTGCCTTTCCCTACCACAGGCT
CCCTGCAGCCCAAGGATGTGCTGAGACACAGCCGGGAGTGCCTCGCCCTCAAACAGGCTCGCTGCAGC
CCAAACATCTGCCTGAGACACAGCCAACACAGCAGGGAGTGCCTTGCCCGCAAACAGGCTCCCTGCAGC
CCAAGGATCTGCCTGAGACACAGCCCGGAGTACTTTTCCAAAGCAGACTCTGTGAACCCAAAGAGC
TGCTTTCAACCCAGCCGGGAGTGCCTCCCCTCACCTGCTCCTCCAGGTGCCGCTCCCCCAGCTAGG
TGCTTGAGGCCTCCCTCCAGGATGCTGCCGCTGCTGTCCCACTGCTCAGGCACACGGCAGAACCCCT
TTGCTACTCCCCCCTCAGAGCCCAACTTCAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001278688

**Insert Size:** 1896 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001278688.2](#)

**RefSeq Size:** 2285 bp

**RefSeq ORF:** 1896 bp

**Locus ID:** 195977

**UniProt ID:** [A6NF34](#)

**Cytogenetics:** 10q11.22

**MW:** 70.6 kDa