

Product datasheet for **SC337043**

GOLGA8R (NM_001282484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA8R (NM_001282484) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA8R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337043 representing NM_001282484.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCAGAAGAACTCAACACAACAAATTGGCTGCAGCCAAGAAAAAGTTAAAAGAATATTGGCAGAAA
AACAGACCTAGAGTTCCAGCAGGAGTGAACAGGAACAGGAAACAAATGGCAGTATCCCTGAGACAGCC
ACTTCCGGTGGTTGCCAGCCACCTGGGGATTAGCAACAGGTTTCCACAGGGAAGGCCCTACATCATCT
GCTACCTGAAAGATCTGGAGAGCCCGTCCAAGAACGAGCAGTAGTCTGGATTCAACGTCCTGTAATA
ATCAGTCGACTGAAGAACACCATCAAACTCTTTGAAACAACAGAAGAAACAAGTGAACATCAGCTGGAA
GAAGAAAAGAAAGCAACACGAGAGACAGAAAGCCGAAAGGGTCTAGAGGTTCAAATCCAGACATTG
ATCATACAGAAAGAGGAATAAATACGGACCTGTACCACATGGAACGTTCTCTCAGATACTTTGAAGAA
GAGTCCAAGGACCTGGCTGTCCGCTGCAACATTCATTGCAGTGTAAAGGAGAGTTAGAGAGGGCTCTG
TCTGCTGTCATCGCCACAGAGAAGAAGGCAAAACAGTTGTCCAGCTGCAGCAAAGCACATACAGAG
TGGGAGTTAGAGCAGTCCCTACAGGACCAGGCACTGCTGAAAGCGCAGCTGACACAGTTGAAGGAGTCA
TTTCAACAACCTCAATTAGAAAGAGATGAGTGTGCTGAACATATAGAAGGAGAGAGGGCCCGGTGGCAT
CAGAGGATGAGTAAATGTGCGAGGAGATTGCACATTAAAGAAAAGAGAAGCAGGATATGCGTTGGGTA
GAGCAGCTGGAGTGGAGCTTGTCCAACTCAAAAACAGACGGCTGAACCTTGCCCCGGAGCCCCCA
GCAGTGCCTCTGAGGTGGAGCTGCAGCACCTGAGGAAGGAACAGAGAGAGTGGCAGGAGAGCTCCAG
TCCCAGGTCAAAAACATCAGCACATAAGTCTCCTGAACCGCGCAAGAAGAGAGGATTGCGGAACAG
GAAGAGAGGCTTCGGAAGCAGGAGGAGAGGCTTCAGGAGCAGCAGAGAAGCTTCGGCAGCTGGCCAAG
CCACAGAGCGTCTTCGAGGAGCTGAACAATGAGAACAAGAGCACACTGCAGTTGGAGCAGCAAGTAAAG
GAGCTACAGGAGAAGCTTGGCGAGGAGCACCTGGAAGTGGCCAGCCAGCAGAAACCAGCAGCTAACGGCC
CAGCTGAGCCTCATGGCTCTCCCTGGGGAAGGACACGAGGAGAAACATCTGGACAGTGAAGGGGAGGAG
GCACCTCAGCCCATGCCGAGTGTCCAGAGGACCCGAGAGCAGGAGGCCATGAGCAGCTTTATGGAC
CACCTGAAGGAGAAGGCAGACCTGAGTGAAGTGTGAAAAACAAGAAGTTCGTTTATTCAATACTGG
CAAGAGAGATGCCATCAGAAAATCCATCACCTTTATCAGAACAGGGGGCCGTGCCAAAGATGCAGCA
CTGGGAGGAGGACACCATCAGGCTGGAGCTCAGGGAGGAGATGAAGGTGAAGCTGCTGGAGCTGCAGCA
GATGGTATTGCGGCTTACAGCAACTACAACATGGGCACAGAAAATTCCTGGCCGCTGCCACAACCTCT
GCTGATGAGCCCGGTCCAGGAGCCCCAGCCCCCAGGAGCTTGGGGCTGCAGACAAGCATGGTGATCTT
CGTGAGGTGACCCTCACCTCTCTGCCAAGGAGAGGCCAGGGAGGATCCTCTCCTTGACAAGCCTACT
GCACAGCCGATCGTGCAGGACCACCAGGAGCACCAGGCTTGGGAGCAACTGCTGTGTGCCATTATTT
TGTTGGGCTTGGCTGCCAAGAAGAGGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001282484

Insert Size: 1896 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282484.1](#)

RefSeq Size: 5172 bp

RefSeq ORF: 1896 bp

Locus ID: 101059918

UniProt ID: [I6L899](#)

Cytogenetics: 15q13.2

MW: 71.5 kDa