

Product datasheet for **SC337048**

GOLGA8M (NM_001282468) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA8M (NM_001282468) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA8M
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337048 representing NM_001282468.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCAGAAGAACTCAACACAACAAATTGGCTGCAGCCAAGAAAAAGTTAAAAGAATATTGGCAGAAA
AACAGCCCTAGAGTTCCAGCAGGAGCGAACAGGAACAGGAAACAAATGGCAGTATCCCTCAGACAGCC
ACTTCTGGTGGTTGCCAGCCACCTGGGGATTACGCAACAGGTTTTTCACAGGGAAGGCCCTACATCATCT
GCTACCTGAAAGATCTGGAGAGCCCGTGCCAAGAACGAGCAGTAGTCTGGATTCAAGGTCCGTAGAA
ATCAGTCAACTGAAGAACACCATCAAACTCTTTGAAACAACAGAAGAAACAAGTGAACATCAGCTGGAA
GAAGAAAAAGCAACAACAAGAAACAGAAAGCCAAAGGGTCTAGAGGTTCAACTCCAGACATTG
AACATACAGAAAGAGGAATAAATACGGACCTGTACCACATGAAACGTTCTCTCAGATACTTTGAAGAA
AAGTCAAAGGATCTGGCTGTCCGCTGCAACATTCATTGCAGCGTAAAGGAGAGTTAGAGAGTGTCTC
TCTGATGTCATGGCCACACAGAAGAAGGCAAAACAGTTGTCCAGCCCCAGTAAAGCAGGTACGGAG
TGGAAGTTAGAGCAGTCCATGCGGGAGGAGGCACTACTGAAAGTGCAGCTGACACAGTTGAAGGAGTCA
TTTCAACAAGTCCAATTAGAAAGAGATGAGTATTCTGAACATCTAAAAGGAGAGAGGGCCCGGTGGCAG
CAGAGGATGAGAAAAATGTCGCAGGAGATTTGCACATTAAAGAAAGAGAGAGCAGCAAGATATGCGTCGG
GTAGAGAAGCTGGAGAGGAGCTTGTCCAACTCAAAAACAGATGGCTGAACCTTGCCCCCGGAGCCC
CCAGCAGTGCCCTCTGAGGTGGAGCTGCAGCACCTGAGGAAGGAAGTAGAGAGAGTGGCAGGAGAGCTC
CAGGCCAGGTCAAAAACAATCAGCGCATAAAGTCTCTGAACAGCGACAAGAAGAGAGGATTGCGGAG
CAGGAAGAGAGGCTTCGGAAGCAGGAGGAGAGGATTGAGGAGCAGCACAAGAGCCTTCAGCAGCTGGCC
AAGCCACAGAGCGTCTTCGAGGAGCCGAACAATGAGAACAAAGAGCAGTGCAGTTGGAGCAGCAAGTA
AAGGAGCTACAGGAGAAGCTTGGCGAGGAGCACCTGGAAGCTGCCAGCCAGCAGAACCAGCAGCTAACG
GCCCAGCTGAGCCTCATGGCTCTCCCTGGGGAAGGACACGGAGGAGAACAATCTGGACAGTGAAGGGGAG
GAGGCACCTCAGCCCATGCCGAGTGTCCAGAGGACCCGAGAGCAGGAGGAGGATGAGCAGCTTTATG
GACCACCTGGAGGAGAAGGCAGACCTGAGTGAAGTGTGAAGAAACAAGAACTTCGCTTCATTCAATAC
TGGCAAGAGAGATGCCATCAGAAAATCCATCACCTTTTATCAGAACAGGGGGCCGTGCCAAGATGCG
GCACTGGGAGGAGACACCATCAGGCTGGAGCTCAGGGAGGAGATGAAGGTGAAGCTGCTGGAGCTGCA
GCAGATGGTATTGCGGCTTACAGCAACTACAACAATGGGCACAGAAAATTCCTGGCCGCTGCCACAAC
TCTGCTGATGAGCCCGTCCAGGAGCCCCAGCTCCCCAGGAGCTTGGGGTGCAGACAAGCATGGTGAT
CTTTGTGAGGTGAGCCTCACCTCTCTGCCAAGGAGAGGCCAGGAGGATCCTCTCTTGACAAGCCT
ACTGCACAGCCGATCGTGCAGGACCACAGGAGCACCCAGGCTTGGGCAGCAACTGCTGTGTGCCATTC
TTTTGCTGGGCTTGGCTGCCAAGAAGAAGGAGATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001282468

Insert Size: 1899 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282468.1](#)

RefSeq Size: 5270 bp

RefSeq ORF: 1899 bp

Locus ID: 653720

UniProt ID: [H3BSY2](#)

Cytogenetics: 15q13.1

MW: 71.5 kDa