

Product datasheet for **SC337056**

p66 alpha (GATAD2A) (NM_001300946) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	p66 alpha (GATAD2A) (NM_001300946) Human Untagged Clone
Tag:	Tag Free
Symbol:	GATAD2A
Synonyms:	p66alpha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337056 representing NM_001300946.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGACCGAAGAAGCATGCCGAACACGGAGTCAGAAACGAGCGCTTGAACGGGACCCAACAGAGGACGAT
GTGGAGAGCAAGAAAATAAAAATGGAGAGAGGATTGTTGGCTTCCAGATTTAAACTGACGGAGACATG
AGGGTGACACCTGAGCCGGGAGCAGGTCCAACCAAGGATTGCTGAGGGCAACAGAGGCCACGGCCATG
GCCATGGGCAGAGGCGAAGGGCTGGTGGCGATGGGCCCGTGGACATGCCACCTCACACAGTGACATG
AAGTCCGAGAGGAGACCCCTCACCTGACGTGATTGTGCTCTCCGACAACGAGCAGCCCTCGAGCCCG
AGAGTGAATGGGCTGACCACGGTGGCTTGAAGGAGACTAGCACCAGGGCCCTCATGAAAAGCAGTCCT
GAAGAACGAGAAAGGATGATCAAGCAGCTGAAGGAAGAATTGAGGTTAGAAGAAGCAAACCTCGTGTG
TTGAAAAAGTTGCGGCAGAGTCAAATACAAAAGGAAGCCACCGCCAGAAGCCCACAGTTCTGTGGG
AGCACCCTGACCACCCCTCCCGCTTGTTCGGGGCACTCAGAACATTCTGCTGGCAAGCCATCACTC
CAGACCTTTCAGCTCGGATGCCCGGCAGTGTACATACCCCGCCCTGGTCCGAGGTGGGCAGCAGGCG
TCCTCGAAGCTGGGGCCACAGGCCAGCTCACAGTCGTATGCCCCCACTCGTCAGGGGGGCTCAGCAA
ATCCACAGCATTAGGCAACATCCAGCAGAGGGCCACCGCCCTCCTCCTGGCCCCCGGGCGTCGGTG
CCCAGTGTGAGATTAGGGACAGAGGATCATCCAGCAGGGCCTCATCCGCGTCGCCAATGTTCCCAAC
ACCAGCCTGCTCGTCAACATCCCACAGCCACCCAGCATCACTGAAGGGGACAACAGCCACCTCCGCT
CAGGCCAACTCCACCCCACTAGTGTGGCTCTGTGGTCACTCTGCCGAGTCTCCAGCAAGCCGACAG
GCGGCCCGCAAGCTGGCGCTGCGCAAACAGCTGGAGAAGACGCTACTCGAGATCCCCCACCAAGCCC
CCAGCCCAGAGATGAATTCCTGCCAGCGCCGCAACAACAGTTCATCTACCTGGTCGGCCTGGAG
GAGGTGGTGAGAACCTACTGGAGACACAAGCAGGCAGGATGTGGCCGCCACTGTGCTGTCCCGGGAG
CCCTACATGTGTGACAGTGCAAGACGGACTTCACGTGCCGCTGGCGGGAGGAGAAGAGCGCGCCATC
ATGTGTGAGAAGTGCATGACAACCAACCAGAAGAAGGCGCTCAAGGTGGAGCACACCAGCCGGCTGAAG
GCCGCTTTGTGAAGGCGCTGCAGCAGGAACAGGAGATTGAGCAGCGGCTCCTGCAGCAGGGCACGGCC
CCTGCACAGGCCAAGGCCGAGCCACCGCTGCCCCACACCCCGTGTGAAGCAGGTCATAAAACCCCGG
CGTAAGTTGGCGTTCGCTCAGGAGAGGGCCGCGACTGGAGTAACGGGGCTGTGCTACAGGCCTCCAGC
CAGCTGTCCCGGGTTCGGCCACGACGCCCCGAGGTGTCTGCACAGTTCAGTCCGTACCCAACTG
CAGAACTCAGCCTCGGCCACAGCCCTGGTCAGCAGGACCGGCAGACATTCTGAGAGAACCGTGAGCGCC
GGCAAGGGCAGCGCCACCTCCAACCTGGAAGAAGACGCCCCCTCAGCACAGGCGGGACCCCTTGCCTTGTG
AGCCCAAGCCTGGCGGTGCAAGAGCTCCTCGGCCGTGGACCGCCAGCAGAGTACCTCCTGGACATG
ATCCACCCCGCTCCATCCCCAGTCAGCCACGTGGAATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001300946

Insert Size: 1905 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001300946.1](#)

RefSeq Size: 5386 bp

RefSeq ORF: 1905 bp

Locus ID: 54815

UniProt ID: [Q86YP4](#)

Cytogenetics: 19p13.11

Protein Families: Transcription Factors

MW: 68.1 kDa

Gene Summary:

Transcriptional repressor. Enhances MBD2-mediated repression. Efficient repression requires the presence of GATAD2B.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.