

Product datasheet for **SC337059**

ZNF133 (NM_001283005) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF133 (NM_001283005) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF133
Synonyms:	pHZ-13; pHZ-66; ZNF150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337059 representing NM_001283005.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGCTGCCGGCCTTCAGGGCGGCCTGGTGGCGCCGAGCATCCTGGGAGGTGTGGTCCGGGCGCT
GTGCCGAGGATTCGGGGTAGTGTAGTCTTGGCGCCCGCTGGAGGAGCTCCCAGCTGGTGGCTGGAGC
GACCCCTTGTCTTTGGTGGCGCTGATCCAGAGCCAGAGCTCTACCTCGATCCTTTCTGCCCTCCGGGT
TTCTCCAGTCAGAAATCCCATGCAGCATGTGCTGTGTAATCATCCCCCTGGATCTTCACATGCTTG
TGTGCAGAAGGTAACATCCAGCCTGGGGATCCGGGCCAGGGGACCAGGAGAAGCAGCAACAAGCCTCT
GAGGGGAGACCCTGGAGTGATCAAGCAGAAGGTCTGAGGGAGAAGGTGCCATGCCTTTGTTTGAAGA
ACCAAGAAAAGGACTCTGGGAGCGTTCTCCAGGCCACCCAGAGGCAGCCAGTCAGCTCTCGAACGGC
CTCAGAGGGGTGGAGTTAGAAGCCAGCCAGCTCAGACAGGGAACCTGAGGAAACAGACAAATTGTTG
AAGAGGATAGAAGTCTTAGGATTTGGAACAGTCAACTGTGGAGAGTGTGGACTGAGCTTCAGCAAGATG
ACAAACCTGCTCAGTCACCAGCGGATACACTCAGGGGAGAAGCCCTACGTGTGTGGGGTATGTGAGAAG
GGCTTCAGCCTAAAGAAGAGCCTCGCCAGACACCAAGAAGGCACACTCGGGGGAGAAGCCAATTGTGTGC
AGGGAGTGTGGACGAGGCTTTAACCGGAAGTCAACGCTAATCATACACGAACGGACACACTCCGGTGAG
AAACCTTACATGTGCAGTGAGTGTGGGCGAGGCTTCCAGCCAGAAGTCAAACCTCATCATAACACAGAGG
ACACACTCAGGGGAAAAGCCTTACGTGTGCCGGGAATGTGGCAAAGGCTTCCAGCCAGAAGTCAGCTGTG
GTGAGACACCAGAGGACACACTTGGAGGAGAAGACCATCGTGTGCAGTGACTGTGGCCTGGGCTTCCAGC
GACAGGTCAAACCTCATCTCCACCAGAGGACGCACTCTGGGGAAGCCCTACGCCTGCAAGGAGTGT
GGGCGATGCTTCAGGCAGAGGACCACCTTGTCAACCACCAGAGGACACACTCAAAGGAGAAGCCCTAT
GTGTGCGGGGTGTGTGGCACAGCTTCAGCCAGAATTCAACCCTCATCTCACAGGCGGACACACACT
GGGAGAAGCCGTATGTTTGTGGGGTGTGTGGCGAGGCTTTAGTCTCAAGTCACACCTCAACAGACAC
CAGAACATACACTCAGGAGAGAAGCCATTGTGTGCAAGGACTGTGGCCGGGCTTCCAGCCAGCAATCC
AACCTCATCAGACACCAGAGGACGCACTCAGGCGAGAAGCCATGGTGTGTGGGAGTGCGGGCGAGGC
TTCAGCCAGAAGTCAAACCTTGTGACACCAGAGGACGCACTCAGGGGAGAGGCCGTATGTGTGCCGA
GAGTGCGGGCGAGGCTTTAGCCACCAGGCCGTCTCATCAGGCACAAGCGGAAGCACTCGAGGGAGAAG
CCCTACATGTGCAGGCAGTGTGGACTGGCTTTGGCAATAAGTCAGCTCTAATTACACACAAGCGGGCT
CACTCGGAAGAGAAGCCTTGTGTGTGCAGAGAGTGTGGCCAAGGCTTTCTCCAAAAGTCACACCTCAC
TTACATCAAATGACACATACGGGGGAGAAGCCATATGTGTGCAAGACGTGTGGGCGGGCTTCCAGCCTC
AAGTCTCACCTCAGCAGACACAGGAAGACCACGTCTGTCCACCACAGACTGCCAGTGCAGCCGACCTT
GAGCCGTGTGCAGGGCAACCTTCGGATTCTTATACTCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001283005

Insert Size: 1908 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001283005.1](#)

RefSeq Size: 2264 bp

RefSeq ORF: 1908 bp

Locus ID: 7692

Cytogenetics: 20p11.23

Protein Families: Transcription Factors

MW: 70.2 kDa

Gene Summary:

May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (13) lacks five internal exons and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (f) has a distinct N-terminus and is shorter than isoform a.