

Product datasheet for **SC337069**

POMGNT1 (NM_001290129) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POMGNT1 (NM_001290129) Human Untagged Clone
Tag:	Tag Free
Symbol:	POMGNT1
Synonyms:	gnT-I.2; GNTI.2; GnT I.2; LGMD2O; LGMDR15; MEB; MGAT1.2; RP76
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337069 representing NM_001290129.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGGGACCAAGGAAATCAAGGTCTCAAATAGAAGAATACAGGAGCACTGGACAGGGCCGTGCTT
TTCCCTGCTGGTGACTGTCATTGTCAATATCAAGTTGATCCTGGACTCGCGAGCCATCAGTGAAGCC
AATGAAGACCCAGAGCCAGAGCAAGACTATGATGAGGCCCTAGGCCGCTGGAGCCCCACGGCGCAGA
GGCAGTGGTCCCCGGCGGGTCTGGACGTAGAGGTGATTCAAGTCGCAGCAAAGTATATGTGGCAGTG
GATGGCACCACGGTGCTGGAGGATGAGGCCCGGGAGCAGGGCCGGGGCATCCATGTCATTGCTCAAC
CAGGCCACGGCCACGTGATGGCAAACGTGTGTTTGACACGTAACCTCATGAGGATGAGGCCATG
GTGCTATTCTCAACATGGTAGCGCCCGGGCAGTGTCTCATCTGCACTGTCAAGGATGAGGGCTCCTTC
CACCTCAAGGACACAGCAAGGCTCTGCTGAGGAGCCTGGGCAGCAGGCTGGCCCTGCCCTGGGCTGG
AGGGACACATGGGCTTCTGTTGGACGAAAAGGAGTCTGTCTTGGGGAGAAACATTCTAAATCACCT
GCCCTCTCTTCTGGGGGACCCAGTCTGCTGAAGACAGATGTGCCATTGAGCTCAGCAGAAGAGGCA
GAGTGCCACTGGGACACAGAGCTGAACCGTGCAGCCGGCGCTTCTGCAGAAAGTTGAGGGCTAT
GGAAGTGTATGCAGCTGAAGGACCCACACCCATCGAGTTCAGCCCTGACCCACTCCCAGACAACAAG
GTCCTCAATGTGCTGTGGCTGTCATTGACAGGAAACCGACCAATTACCTGTACAGGATGCTGCGCTCT
CTGCTTTCAGCCAGGGGGTGTCTCCTCAGATGATAACAGTTTTTATTGACGGCTACTATGAGGAACCC
ATGGATGTGGTGGCACTGTTGGTCTGAGGGGCATCCAGCATACTCCCATCAGCATCAAGAATGCCCGC
GTGTCTCAGCACTACAAGGCCAGCCTCACTGCCACTTTCAACCTGTTCCGGAGGCAAGTTTGTGTG
GTTCTGGAAGAGGACCTGGACATTGCTGTGGATTTTTTTCAGTTTCTGAGCCAATCCATCCACCTACTG
GAGGAGGATGACAGCCTGTACTGCATCTGCTGGAATGACCAGGGGTATGAACACACGGCTGAGGAC
CCAGCACTACTGTACCGTGTGGAGACCATGCCTGGGCTGGGCTGGGTGCTCAGGAGGTCCTGTACAAG
GAGGAGCTTGAGCCCAAGTGGCCTACACCGGAAAAGCTCTGGGATTGGGACATGTGGATGCGGATGCTT
GAACAACGCCGGGGCCGAGAGTGATCATCCCTGACGTTTCCCGATCCTACCACTTTGGCATCGTCGGC
CTCAACATGAATGGCTACTTTACAGAGGCTACTTCAAGAAGCACAAGTTCAACACGGTTCCAGGTGTC
CAGCTCAGGAATGTGGACAGTCTGAAGAAAGAAGCTTATGAAGTGAAGTTACAGGCTGCTCAGTGAG
GCTGAGGTTCTGGACCACAGCAAGAACCCTTGTGAAGACTCTTCTGCCAGACACAGAGGGCCACACC
TACGTGGCCTTTATTGAAATGGAGAAAGATGATGACTTACCACCTGGACCCAGCTTCCAAGTGCCTC
CATATCTGGGACCTGGATGTGCGTGGCAACCATCGGGGCTGTGGAGATTGTTTCGGAAGAAGAACCAC
TTCTGTTGGTGGGGGTCCCGGCTTCCCTACTCAGTGAAGAAGCCACCCTCAGTACCCCAATTTTC
CTGGAGCCACCCCAAGGAGGAGGGAGCCCCAGGAGCCCCAGAACAGACATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001290129

Insert Size: 1917 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290129.1](#)

RefSeq Size: 2400 bp

RefSeq ORF: 1917 bp

Locus ID: 55624

Cytogenetics: 1p34.1

Protein Pathways: O-Mannosyl glycan biosynthesis

MW: 72.4 kDa

Gene Summary:

This gene encodes a type II transmembrane protein that resides in the Golgi apparatus. It participates in O-mannosyl glycosylation and is specific for alpha linked terminal mannose. Mutations in this gene may be associated with muscle-eye-brain disease and several congenital muscular dystrophies. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Feb 2014]

Transcript Variant: This variant (3) differs in the 5' and 3' UTRs, lacks a portion of the 5' coding region, and initiates translation at a downstream start codon compared to variant 1. The resulting protein (isoform 3) is shorter and has a distinct N-terminus compared to isoform 1.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.