

Product datasheet for **SC337072**

USHBP1 (NM_001297703) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USHBP1 (NM_001297703) Human Untagged Clone
Tag:	Tag Free
Symbol:	USHBP1
Synonyms:	AIEBP; MCC2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337072 representing NM_001297703.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTCCACCCGACTGACAAGAAGATGGATGGGGGCTCTGGCAGGGAAGTGGCCTCAGCCCCGAAGTG
CCCCACAACCTGCAGTGGAGGCCACCAAGCCCCAGAAGCAGCCCTACAGTACAAGGAGACTGTGCC
CCTGGGAATGGGGCCCCGATGTGTTTCAGACCTCCAGCACACTCTGAGCTCCCTGGAGGCAGCGGCT
GCAGCCTGGCGCCACCAGCCCCCAGCCATTCTGGGCCGATGGAGTTCGAAGGCACAAGCGAAGGGGA
GCAGGTTCCCTTGGGAAGCAGGAAGGGGCAGGAGCTGCCAGCGAGAGGCAGCTCGCCTGGCCGAGAGG
AATGCCTGGCTCCGGCTGGCCCTGAGTAGCCGAGAGGATGAGCTGGTCCGACGCAGGCCCTCCCTGGAG
GCCATCCGAGCTGAGAAGGAGACGCTGCAGAAAGAGGTTCAAGAGCTGCAGGATTCCTCTTGAGGCTG
GAGCCCTGCCACATCTTCCACAACCAAGCAGGTGGCTCAGGCAGTGGTTCTAGCAGCTCTGAGGCA
GACAGGGAACCTGGGAGACTCAGGACTCCTTCTCCCTGGCTCACCCCTGCTTCGGCGCCTCCGCAGC
CATTCCAGCACCCAGATCCTTGGGTCTCTTCCCAACCAGCCCTCAGTCTGAGATGCACATCATGGAA
GCCAGATGGAGCAACTCCGGGGGAGCATTGAGAAGCTCAAATGCTTTAATCGTCTGCTATCAGCTGTG
CTACAGGGATACAAGGGCCGCTGTGAAGGCCTCAGCATGCAGCTAGGCCAGCGGGAGGCTGAGGCCACA
GCATTGCATCTGGCCTTGCAGTACAGTGAACACTGTGAAGAGGCATACAGGGTTCTGCTTGTCTGCGG
GAGGCCGACTCAGGAGCAGGAGACGAAGCCCCCATGAGTGACCTGCAAGCAGCTGAAAAGGAAGCATGG
AGGCTGCTGGCACAAGAGGAGGCTGCCATGGATGCAGGAGCACAGCAGAAATCCACAGCCAAGCCCTGAA
GGCAGCAGTGTGGATAAGCCCACCCACAGGAAGTGGCTTCCAGCTCCGAAGCTATGTCCAGCGTCTC
CAGGAGCGCGCTTCTCTAATGAAGATTCTCTCAGAGCCTGGCCCCACCTTGGCACCCATGCCACTGTG
CCCCGTGCAGAAGCCATGGTGCAGGCCATTCTGGGGACCCAGGCTGGCCAGCTCTTCCCGACTGGAG
AAGACACAAATTCAGCAGGACCTGGTGGCCGCAAGGGAGGCCCTGGCGGACCTGATGCTTCGGCTGCAG
CTGGTGGCGCTGAGAAGCGGGGCTAGAGCTGCGGGAGGCTGCCCTCCGAGCCCTGGGTCCAGCTCAC
GTGCTCCTGCTGGAGCAGCTGCGGTGGGAACGGGCAGAGCTCCAGGCTGGTGGGGCAAACAGCAGTGGC
GGACATAGCAGCGGAGGTGGCAGCAGCGGGACAGGAAGAGTGGTATCAGGGCCTTCTGCTGTCCCT
GGTGGCACCAGTGGCATTGATGGTGGCCAGGTGGGCAGAGCCTGGGACCCAGAGAAGTTAGCCCAGGAA
CTGGCAGCATCGCTCACCAGGACACTGGACCTGCAGGAGCAGCTGCAGTCTCTGCGCAGGAGCTGGAA
CAGGTGGCTCAGAAGGGCGAGCCAGACGCTCTCAGAGTCCGAGCTGAACAGGGATTTATGCAAGGCC
CACAGCGCCTGGTCTGGCCTTCCGAGGAGCCACCGAAGCAGGAAGCAACGCCGGAAGTTGGAG
CAGCAGATGGCACTCATGGAGGCTCAGCAGGCCGAGGAGGTGGCGGTGCTCGAAGCCACTGCAAGGGCC
CTGGGGAAGCCAGGCCTCCCCTCCCGCCTCCCGAGCTTGGGGACACCTTTCTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001297703

Insert Size: 1920 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297703.1](#)

RefSeq Size: 2354 bp

RefSeq ORF: 1920 bp

Locus ID: 83878

UniProt ID: [Q8N6Y0](#)

Cytogenetics: 19p13.11

MW: 69.6 kDa