

Product datasheet for **SC337073**

PIB5PA (INPP5J) (NM_001284288) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIB5PA (INPP5J) (NM_001284288) Human Untagged Clone
Tag:	Tag Free
Symbol:	INPP5J
Synonyms:	INPP5; PIB5PA; PIPP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337073 representing NM_001284288.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCCCTCCCAAGGCTTGGCACACAGAGTACAGGGCCCTGGCAGGTGCCTGAGCCCCAACCTTCAGGCC
CAAGAAGCCCCAGCCCCAGTACCACCTCCTCTTACATCCACCCTGTATCCTCCCCTTGGTCAGCT
CAGCCTACCTGGAAGAGCGACCCCGCTTCCGGATCACTGTGGTACATGGAACGTGGGCACTGCCATG
CCCCAGACGATGTACATCCCTCCTCCACCTGGGCGGTGGTGACGACAGCGACGGCGCAGACATGATC
GCCATAGGGTTGCAGGAAGTGAAGTCCATGCTCAACAAGCGACTCAAGGACGCCCTTTCACGGACCAG
TGGAGTGAGCTGTTTATGGATGCGCTAGGGCCCTTCAACTTCGTGCTGGTGGTTCGGTGGAGTGCAG
GGTGTATCCTGCTGCTGTTTCGCAAGTACTACCACCTGCCCTTCTGCGAGACGTGCAGACCGACTGC
ACGCGCACTGGCCTGGGCGGCTACTGGGTAAAGGGTGGCGTGAGCGTGCCTGGCGCCCTTCGGG
CACATGCTCTGCTTCTGAACTGCCACTTGCTGCGCATATGGACAAGGGCGAGCAGCGCAAGACAAC
TCCAGACCATCCTCAGCTCCAGCAGTCCAAGGGCCGGCGCACAGGGCATCCTGGATCATGACCTC
GTGTTCTGGTTCGGGGACCTGAACTTCCGATTGAGAGCTATGACCTGCACTTTGTCAAGTTTGCCATC
GACAGTGACCAGCTCCATCAGCTCTGGGAGAAGGACCAGCTCAACATGGCCAAGAACACCTGGCCATT
CTGAAGGGCTTTTCAAGGGGGCCCTCAACTTCGCTCCCACCTTCAAGTTTGATGTGGGTACCAACAAA
TACGATACCAGTGCCAAGAAACGGAAGCCAGCTTGACAGACCGTATCCTATGGAAGGTCAAGGCTCCA
GGTGGGGTCCCAGCCCCCAGGACGGAAGAGCCACCGACTCCAGGTGACGCAGCACAGCTACCGCAGC
CACATGGAATACACAGTCAGCGACCACAAGCCTGTGGCTGCCAGTTCCTCCTGCAGTTTGCCTTCAGG
GACGACATGCCACTGGTGGCTGGAGGTGGCAGATGAGTGGGTGCGGCCGAGCAGGCGGTGGTGGAG
TACCGCATGGAACAGTGTTCGCCCGCAGCTCCTGGGACTGGATCGGCTTATACCGGGTGGGTTTCCGC
CATTGCAAGGACTATGTGGCTTATGTCTGGGCCAAACATGAAGATGTGGATGGAATACCTACCAGGTA
ACATTCAGTGAGGAATCACTGCCAAGGGCCATGGAGACTTCATCCTGGGCTACTATAGTACAACCAC
AGCATCCTCATCGGCATCACTGAACCTTCCAGATCTCGCTGCCTTCTCGGAGTTGGCCAGCAGCAGC
ACAGACAGCTCAGGCACCAGCTCAGAGGGAGAGGATGACAGCACACTGGAGCTCCTTGCACCCAAGTCC
CGCAGCCCCAGTCTGGCAAGTCCAAGCGACACCCGAGCCGAGCCGGGACTGGCCAGGTTCCCTGGG
CTTGCCCTACGGCCCTCATCCCCTGAACGCCGTGGTGGCAGCCGTAGCCCTCACCCAGAGCCCGCCG
CTGTCCCAGTGGCTCCTGACAGGAGCAGTAATGGCAGCAGCCGGGCGAGTGTGAAGAGGGGCCCTCT
GGTTGCCTGGCCCTGGGCTTCCACCAGCTGTGCTCGAAGCCTGGGCTGTGCCCCGCTTGGC
CTAGAGACTGTAGACCCTGGTGGTGGTGGCTCCTGGGGACCTGATCGGGAGGCCCTGGCGCCCAACAGC
CTGTCTCCTAGTCCCCAGGGCCATCGGGGGCTGGAGGAAGGGGGCCTGGGGCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001284288

Insert Size: 1920 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001284288.1](#)

RefSeq Size: 2481 bp

RefSeq ORF: 1920 bp

Locus ID: 27124

UniProt ID: [Q15735](#)

Cytogenetics: 22q12.2

Protein Families: Druggable Genome

Protein Pathways: Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system

MW: 70.5 kDa

Gene Summary: Inositol 5-phosphatase, which converts inositol 1,4,5-trisphosphate to inositol 1,4-bisphosphate. Also converts phosphatidylinositol 4,5-bisphosphate to phosphatidylinositol 4-phosphate and inositol 1,3,4,5-tetrakisphosphate to inositol 1,3,4-trisphosphate in vitro. May be involved in modulation of the function of inositol and phosphatidylinositol polyphosphate-binding proteins that are present at membranes ruffles (By similarity).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (5) differs in its 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (e) is shorter at the N-terminus, compared to isoform a. Both variants 5 and 6 encode isoform e.