

Product datasheet for **SC337076**

ZNF343 (NM_001282497) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF343 (NM_001282497) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF343
Synonyms:	dj734P14.5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337076 representing NM_001282497.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGTTGCCTTATCCTTCAGCACTGGGAGATCAATACTGGGAAGAGATTTTGCTTCAAAGAATGGG
GAAAATGTAGAGACTATGAAGAAATTGACCCAAAATCATAAAGCGAAAGGCTTGCCTTCTAATGATACT
GACTGCCCCAGAAAAGGAGGGAAAGGCCAAAATAGTGGTACCAGTTACATTCAGGGATGTGACTGTG
ATCTTCACAGAAGCAGAATGGAAGAGACTGAGTCCAGAGCAGAGGAATCTATACAAGAAGTGATGCTG
GAGAATTACAGGAATCTTCTCTCATTGGATGGAGTTTCGCTCTTGTGCCCAGGCTAGAGTGCAATCGT
GCAATCTCGGCTCACCACAACCTTTGCCTCCCGGGTTCAAGCAATTCTCCTGCCACAGCCTCTCGAGTA
GTTGGGATTACAGCAGAACCAAAGCCAGAAATCTACACTTGTTCCTCCTGCCTTCTGCCTTCTCCTGT
CAGCAGTTCCTCAGTCAACATGTACTTCAGATCTTCTGGGCTTATGTGCAGAAAATCACTTCCATCCA
GGGAATTCAGCCAGGGCATTGGAAACAGCAGGGGCAGCAGTATCCCATGTAAGCTGTTGGTTTGAA
AATGCAGAAGGTCAGGAGAGAGGAGGTGGCTCAAACCCTGGTCTGCAAGGACAGAGGAGAGAGAAAACC
TCAAGGGCATTCCCAGCCCACTCAAAGACAGTCAGCAAGTCTAGAAAAGGCAACATGGTGGTAGAA
ACAGAGCCAGCTCAGCCCAAAGACCAAACCTGTGCAGCTAGACAAAGGCTTGAAGGAATTAGAAAACC
TTGAGATTTGGAGCAATCAACTGTAGAGAGTATGAACCGGACCATAACCTGGAATCAAACCTTTATTACA
AACCCGAGGACCTCTTAGGGAAGAAGCCCTACATTTGCAGTGATTGTGGGCGAAGCTTTAAGATAGA
TCAACCCTCATCAGACACCATCGTATACACTCGATGGAGAAGCCTTATGTGTGCAGTGAGTGCAGGCGA
GGTTTTAGCCAGAAGTCCAACCTCAGCAGACACCAGAGAACACATTCAGAAGAGAAGCCTTATTTGTGC
AGGCCTATGTTTGCAGCGAGTGTGGGCGAGGCTTTAGCGAGAAGTCATCCTTTCATCAGACACCAGAGG
ACACACTCCGGTGAGAAAACCTATGTGTGCCTGGAGTGTGGACGAAGCTTTTGTGATAAGTCAACCTC
AGAAAACACCAGAGGATACACTCAGGGGAGAAGCCTTATGTTTGCAGGGAGTGTGGGCGAGGCTTTAGC
CAGAACTCAGATCTCATCAAACACCAGAGGACACACTTGGATGAGAAGCCTTATGTTTGCAGGGAGTGT
GGGCGAGGCTTTTGTGACAAGTCAACCCTCATCATACAGAGCGGACGCACTCTGGAGAGAAGCCTTAT
GTGTGTGGTGTGAGTGTGGCCGAGGCTTTAGTCGAAATCACTCCTCCTTGTCCACCAGAGGACACACTCA
GGGAGAAGCATTATGTCTGCAGGGAGTGTAGGCGAGGTTTTAGCCAGAAGTCAAATCTCATCAGACAC
CAGAGGACGCACTCAAATGAGAAGCCTTATATTTGCAGGGAATGTGGGCGAGGCTTTTGTGACAAGTCA
ACCCTCATTGTACATGAGAGGACACACTCAGGAGAGAAGCCTTACGTGTGCAGTGAGTGTGGCCGAGGC
TTTAGCCGGAATCACTCCTCCTTGTCCACCAGAGGACACACTCAGGGGAGAAGCATTATGTTTGTAGG
GAGTGTGGGCGAGGCTTTAGTCATAAGTCAAATCTCATCAGACACCAGAGGACACACTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001282497

Insert Size: 1923 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282497.1](#)

RefSeq Size: 3767 bp

RefSeq ORF: 1923 bp

Locus ID: 79175

UniProt ID: [Q6P1L6](#)

Cytogenetics: 20p13

Protein Families: Transcription Factors

MW: 73.4 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) contains an alternate in-frame exon in the central coding region, compared to variant 1. The encoded isoform (2) is longer, compared to isoform 1. Both variants 4 and 8 encode the same isoform (2).