

## Product datasheet for **SC337083**

### **PDE4 (PDE4B) (NM\_001297440) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PDE4 (PDE4B) (NM_001297440) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE4B
Synonyms:	DPDE4; PDEIVB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337083 representing NM\_001297440.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCTTTGATGTGGAAAATGGCCCTTCCCAGGTCCGAGTCCACTGGATCCCCAGGCCAGCTCTTCC
GCTGGGCTGGTACTTCACGCCACCTTTCCTGGGCACAGCCAGCGCAGAGAGTCATTTCTCTACAGATCA
GACAGCGACTATGACTTGTCAACAAAGGCGATGTCGAGAAACTCTTCTTCCAAGCGAGCAACACGGC
GATGACTTGATTGTAACCTTTTGGCCAGTCTTGGCAGCTTGCAGTGTGAGAAACAACCTTCACT
ATACTGACAAACCTTCATGGTACATCTAACAAGAGGTCCCAGCTGCTAGTCAGCCTCCTGTCTCCAGA
GTCAACCACAAGAAGAATCTTATCAAAAATTAGCAATGGAAACGCTGGAGGAATTAGACTGGTGTTTA
GACCAGCTAGAGACCATAACAGACCTACCGGTCTGTGAGTGGCTTCTAACAAGTTCAAAGAATG
CTGAACCCGGGAGCTGACACACCTCTCAGAGATGAGCCGATCAGGGAACCGGTGTCTGAATACATTTCA
AATACTTTCTTAGACAAGCAGAATGATGTGGAGATCCCATCTCCTACCCAGAAAGACAGGGAGAAAAAG
AAAAAGCAGCAGCTCATGACCCAGATAAGTGGAGTGAAGAAATTAATGCATAGTTCAAGCCTAAACAAT
ACAAGCATCTCACGCTTTGGAGTCAACACTGAAAATGAAGATCACCTGGCCAAGGAGCTGGAAGACCTG
AACAAATGGGGTCTTAACATCTTAAATGTGGCTGGATATTCTACAATAGACCCCTAACATGCATCATG
TATGCTATATTCCAGGAAAGAGACCTCCTAAAGACATTCAGAATCTCATCTGACACATTTATAACCTAC
ATGATGACTTTAGAAGACCATTACCATTCTGACGTGGCATATCACAACAGCCTGCACGCTGCTGATGTA
GCCAGTCGACCCATGTTCTCCTTTCTACACCAGCATTAGACGCTGTCTTACAGATTTGGAGATCCTG
GCTGCCATTTTTCAGCTGCCATCCATGACGTTGATCATCCTGGAGTCTCCAATCAGTTTCTCATCAAC
ACAAATTCAGAACTTGCTTTGATGTATAATGATGAATCTGTGTTGGAAAATCATCACCTTGCTGTGGGT
TTCAAACTGCTGCAAGAAGAACAACACTGTGACATCTTCATGAATCTCACAAGAAGCAGCGTCAGACACTC
AGGAAGATGGTTATTGACATGGTGTAGCAACTGATATGTCTAAACATATGAGCCTGCTGGCAGACCTG
AAGACAATGGTAGAAACGAAGAAAGTTACAAGTTCAGGCGTTCTTCTCCTAGACAATATACCGATCGC
ATTCAGTCTTTCGCAACATGGTACACTGTGCAGACCTGAGCAACCCACCAAGTCTTGGAAATGTAT
CGGCAATGGACAGACCGCATCATGGAGGAATTTTCCAGCAGGGAGACAAAGAGCGGGAGAGGGGAATG
GAAATTAGCCCAATGTGTGATAAACACACAGCTTCTGTGGAAAAATCCCAGGTTGGTTTCATCGACTAC
ATTGTCCATCCATTGTGGGAGACATGGGCAGATTTGGTACAGCCTGATGCTCAGGACATTCTCGATACC
TTAGAAGATAACAGGAACTGGTATCAGAGCATGATACCTCAAAGTCCCTCACCACCTGGACGAGCAG
AACAGGGACTGCCAGGTCTGATGGAGAAGTTTCAGTTTGAAGTACTCTCGATGAGGAAGATTCTGAA
GGACCTGAGAAGGAGGGAGAGGGACACAGCTATTTTCAGCAGCACAAAGACGCTTTGTGTATTGATCCA
GAAAACAGAGATTCCTGGGAGAGACTGACATAGACATTGCAACAGAAGCAAGTCCCCCGTGGATACA
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001297440

**Insert Size:** 1935 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001297440.1](#)

**RefSeq Size:** 4070 bp

**RefSeq ORF:** 1935 bp

**Locus ID:** 5142

**UniProt ID:** [Q07343](#)

**Cytogenetics:** 1p31.3

**Protein Families:** Druggable Genome

**Protein Pathways:** Progesterone-mediated oocyte maturation, Purine metabolism

**MW:** 73.2 kDa

**Gene Summary:** This gene is a member of the type IV, cyclic AMP (cAMP)-specific, cyclic nucleotide phosphodiesterase (PDE) family. The encoded protein regulates the cellular concentrations of cyclic nucleotides and thereby play a role in signal transduction. Altered activity of this protein has been associated with schizophrenia and bipolar affective disorder. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014]

Transcript Variant: This variant (e) differs in the 5' UTR and lacks an alternate exon in the 5' coding region and initiates translation at an alternate downstream start codon, compared to variant 1. The encoded isoform (4) has a shorter and distinct N-terminus when compared to isoform 1.