

## Product datasheet for **SC337085**

### SUV420h1 (KMT5B) (NM\_001300908) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SUV420h1 (KMT5B) (NM_001300908) Human Untagged Clone
Tag:	Tag Free
Symbol:	KMT5B
Synonyms:	CGI-85; CGI85; MRD51; SUV420H1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337085 representing NM\_001300908.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTACTTAGACATGGAGAAAACGACTTCAGTGTATGTACTCCACAAGAAAACTGTGCTCACTC
TGGTGGGTCTGCTGCGTTTATAAACCATGATTGCAGACCTAATTGTAAGTTTGTGCACTGGTCGA
GATACAGCATGTGTGAAGGCTCTAAGAGACATTGAACCTGGAGAAGAAATTTCTGTTATTATGGAGAT
GGTTCCTTTGGAGAAAATAATGAGTTCTGCGAGTGTACACTTGCAGAAAGACGGGCACTGGTGCTTTT
AAATCCAGAGTGGGACTGCCTGCGCCTGCTCCTGTTATCAATAGCAAATATGGACTCAGAGAAACAGAT
AAACGTTTAAATAGGCTTAAAAAGTTAGGTGACAGCAGCAAAAATTCAGACAGTCAATCTGTCAGCTCT
AACACTGATGCAGATACCACTCAGGAAAAAACAATGCAACTTCTAACCGAAAATCTTCAGTTGGCGTA
AAAAAGAATAGCAAGAGCAGAACGTTAACGAGGCAATCTATGTCAAGAATTCAGCTTCTCCAACCTCT
ACCTCATCTAAGCTAACTCATATAAATAATTCCAGGGTACCAAAGAACTGAAGAAGCCTGCAAAGCCT
TTACTTTCAAAGATAAAATTTGAGAAATCATTGCAAGCGGCTGGAGCAAAGAATGCTTCAAGAAAACCTC
GAAATGGGAAACTTAGTACTGAAAGAGCCTAAAGTAGTCTGTATAAAAAATTTGCCATTAAAAAGAT
AAGGAGCCAGAGGGACCAGCCAAAGCCGAGTTGCCAGCGGGTGCTTGACTAGACACGCGCGGAGAGAA
CACAGACAGAATCCTGTGAGAGGTGCTCATTGCGAGGGGAGAGCTCGCCCTGCACCTACATAACTCGG
CGGTCACTGAGGACAAGAACAACTCTGAAGGAGGCCTCTGACATCAAGCTTGAACCAAATACGTTGAAT
GGCTATAAAAGCAGTGTGACGGAACCTTGCCCCGACAGTGGTGAACAGCTGCAGCCAGCTCCTGTGCTG
CAGGAGGAAGAACTGGCTCATGAGACTGCACAAAAGGGGAGGCAAAAGTGCATAAGAGTGACACAGGC
ATGTCCAAAAGAAGTCAAGACAAGGAAAACCTGTGAAACAGTTTGCAAAAATAGAGGAATCTACTCCA
GTGCACGATTCTCTGGAAAAGACGACGCGGTACCAGATTTGATGGGTCCCATTCTGACCAGGGTGAG
CACAGTGGCACTGTGGGCGTGCCTGTGAGCTACACAGACTGTGCTCCTTACCCGTCGGTTGTTCAAGTT
GTGACATCAGATAGCTTCAAACAAAAGACAGCTTTAGAACTGCAAAAAGTAAAAAGAAGAGGCGAATC
ACAAGGTATGATGCACAGTTAATCCTAGAAAATAACTCTGGGATCCCAAATGACTCTTCGTAGGCGT
CATGATAGCAGCAGCAAAAACAAATGACCAAGAGAATGATGGAATGAACTTCCAAAATAAGCATCAAG
TTAAGCAAAGACCATGACAACGATAACAATCTCTATGTAGCAAAGCTTAATAATGGATTTAACTCAGGA
TCAGGCAGTAGTTCTACAAAATAAAAATCCAGCTAAAACGAGATGAGGAAAATAGGGGGTCTTATACA
GAGGGGCTTCATGAAAATGGGGTGTGCTGCAGTGATCCTTTCTCTCTGGAGTCTCGAATGGAGGTG
GATGACTATAGTCAGTATGAGGAAGAAAGTACAGATGATTCTCTCTCTGAGGGCGATGAAGAGGAG
GATGACTATGATGATGACTTTGAAGACGATTTTATTCCTCTCTCCAGCTAAGCGCTTGAGGTTAATA
GTTGAAAAGACTCTATAGATATTGACATTTCTTCAAGGAGAAGAGAAGATCAGTCTTTAAGGCTTAAT
GCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001300908

**Insert Size:** 1938 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001300908.1](#)

**RefSeq Size:** 5802 bp

**RefSeq ORF:** 1938 bp

**Locus ID:** 51111

**Cytogenetics:** 11q13.2

**Protein Families:** Druggable Genome

**Protein Pathways:** Lysine degradation

**MW:** 71.6 kDa

**Gene Summary:**

This gene encodes a protein that contains a SET domain. SET domains appear to be protein-protein interaction domains that mediate interactions with a family of proteins that display similarity with dual-specificity phosphatases (dsPTPases). The function of this gene has not been determined. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]

Transcript Variant: This variant (4) uses an alternate first exon and an alternate splice junction at the 5' end of an exon, and uses a downstream start codon compared to variant 1. The resulting isoform (4) is shorter at the N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.