

Product datasheet for **SC337088**

ZNF567 (NM_001300979) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF567 (NM_001300979) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF567
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337088 representing NM_001300979.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCTCAGGGATCAGTGTCTTTCAATGATGTGACTGTGGACTTCACTCAGGAGGAGTGGCAGCACCTG
GATCATGCTCAGAAGACTCTATATATGGATGTGATGTTGGAAAACCTATTGCCACCTCATCTCTGTGGGG
TGTACATGACCAAACCTGATGTGATCCTCAAGTTGGAACGAGGAGAAGAGCCATGGACATCATTGCA
GGTCATACCTGCTTGAAGAAAACCTGGAAGCTGAAGACTTTTTAGTGAATTCAAGGAACCAAGAG
AAGTATTCTAGATCAGTTGTAAGCATCAACCACAAAAAAGCTGGTGAAGGAGAAGAGTAAATATATGAA
AAGACATTTACTCTAGGCAAAAACCTGTGAATCAAAAAATCTACCTCTGAATATGATACTCATGGA
AGGATTTTGA AAAATGTTTCAGAATTAATCATCAGTAATCTAAATCCTGCAAGAAAGAGACTTAGTGAG
TATAATGGATATGGGAAATCACTCCTGAGTACTAAACAAGAGACTACTCATCCTGAAGTCAAATCCCAT
AATCAAAGTGCCAGAGCTTTCAGTCATAATGAAGTTCTTATGCAGTATCAGAAAACGAAACTCCAGCA
CAGTCATTTGGATATAATGACTGTGAGAAATCATTCTTCAAAGGGGAGGCCTGATTACACATAGTAGA
CCTTACAAAGGAGAAAACCCATCTGTATATAATAAAAAAGAAGAGCAACCAATATTGAAAAAACAT
ACATGCAATGAATGTGGGAAATCTTTCTGCAGGAAATCAGTATTGATTCTGCATCAGGGAATTCACTCA
GAAGAAAAACCCATCAATGTCATCAATGTGGAATGCATTTAGAAGGAAATCATATCTCATTGATCAT
CAGAGAACTCACACAGGAGAGAAACCCTTTGTGCAATGAATGTGGTAAGTCTTCCGCCTCAAGACA
GCCCTCACTGATCATCAGAGAACACACAGGGGAGAAATCGTATGAATGTCTGCAATGTAGGAATGCC
TTCAGATTGAAGTCACACCTCATTGTCATCAGAGAACTCACACGGGAGAGAAACCATATGAGTGAAT
GACTGTGGGAAGTCCTTCCGCCAGAAGACAACACTCTCTACATCAGAGAAATCCATACAGGTGAGAAA
CCCTACATTTGTAAAGAATGTGGGAAGTCCTTTCACCAGAAGGCAAATCTTACTGTACATCAGAGAACT
CATACAGGGGAAAAGCCCTATATTTGTAATGAATGTGGGAAATCCTTCTCCAGAGACAACCCCTTGCT
CTTCATGAGAAAACCTATAATGAGGAGAAACCTATATTTGATGTAATGTGGAAGTCCTTCCGCCAG
AAGACAACCCCTGTAGCACATCAGAGAACATACAGGGGAGAAATCTTATGAATGTCTCACTGTGGG
AAGGCCTTTAGAATGAAGTCATACCTCATTGATCATCACCGAACTCACACAGGAGAGAAACCATATGAA
TGTAATGAATGTGGTAAATCATTAGTCAAAAGACAAATCTCAATCTACATCAGAGAAATTCATACAGGG
GAGAAACCCATGTTTGTAAATGAATGTGGGAAGTCCTTTCGCCAGAAAGCAACCCCTCACTGTACATCAG
AAAATACATACCGGCCAGAAATCCTATGAATGTCCTCAGTGTGGGAAAGCCTTTAGCAGGAAGTCATAT
CTCATTATCATCAAAGAACTCATACGGGAGAGAAACCATATAAATGTAGTGAATGTGGAAGTGCTTC
CGCCAGAAAGACAAATCTTATTGTACATCAGAGAACTCACACAGGTGAGAAACCCATGTTTGTAAATGAG
TGTGGTAAAGTCCTTTCAGTTATAAGAGAAACCTCATTGTCCATCAAAGAACTCACAAGGGAGAAAACATT
GAAATGCAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001300979

Insert Size: 1944 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001300979.1</u>
RefSeq Size:	2871 bp
RefSeq ORF:	1944 bp
Locus ID:	163081
UniProt ID:	<u>Q8N184</u>
Cytogenetics:	19q13.12
Protein Families:	Transcription Factors
MW:	75.2 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]