

## Product datasheet for **SC337090**

### ZNF202 (NM\_001301779) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF202 (NM_001301779) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF202
Synonyms:	ZKSCAN10; ZSCAN42
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

**Fully Sequenced ORF:** >SC337090 representing NM\_001301779.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTACAGCCGTGGAACCCAGAGGACCAGGATCTTTGGGAAGAAGAGGGAATTCTGATGGTGAACTG
GAAGATGATTTACCTGTCCGCCAGAGTCTGTCTTACAGAGGGATGACCCGGTCTGGAACCTCCCAC
CAGAACTCCGACGCTTCCGCTACCAGGAGGCAGCAAGCCCTAGAGAAGCTCTCATCAGACTCCGAGAA
CTTTGTCACCAGTGGCTGAGACCAGAGAGGCGGACAAAGGAGCAGATCCTAGAGCTGCTTGTGCTGGAA
CAATTTCTTACCGTCTACCTGGAGAACTACAGAGCTGGGTGCGGGGCCAACGGCCAGAAAAGTGGCGAG
GAGGCAGTGACGCTGGTGGAGGGTTTGCAGAAACAACCCAGGAGACCAAGGCGGTGGGTGACTGTCCAT
GTTACAGGCCAGGAAGTCTGTGAGAGGAGACGGTGCATTTAGGAGTGGAGCCTGAGTCACCTAATGAG
CTGCAGGATCCTGTGCAAAGCTCGACCCCGAGCAGTCTCCTGAGGAAACCACACAGAGCCAGATCTG
GGGGCACCCGAGAGCAGCGTCCACACCAGGAAGAGGAGCTCCAGACCCTGCAGGAGAGCGAGGTCCCA
GTGCCCGAGGACCAGACCTTCTGCAGAGAGGAGCTCTGGAGACTCAGAGATGGTTGCTTCTTACT
GCTCTGTCACAGGACTGGTAACGTTCAAGGATGTGGCCGTATGCTTTCCAGGACCAGTGGAGTGAT
CTGGACCAACACAGAAAGATTCTATGGAGAATATGTCTTGGAAGAAGACTGTGGAATTGTTGTCTCT
CTGTCATTTCAAATCCCCAGACCTGATGAGATCTCCAGGTTAGAGAGGAAGAGCCTTGGTCCCAGAT
ATCCAAGAGCCTCAGGAGACTCAAGAGCCAGAAAATCCTGAGTTTTACCTACACAGGAGATAGGAGTAAA
GATGAGGAAGAGTGTCTGGAGCAGGAAGATCTGAGTTTGGAGGATATACACAGGCCTGTTTTGGGAGAA
CCAGAAATTCACCAGACTCCAGATTGGGAAATAGTCTTTGAGGACAATCCAGGTAGACTTAATGAAAGA
AGATTTGGTACTAATATTTCTCAAGTGAATAGTTTTGTGAACCTTCGGGAAACTACACCCGTCACCCC
CTGTTAGGGAGGCATCATGACTGTTCTGTGTGTGGAAAGAGCTTCACTTGTAACCTCCACCTTGTAGA
CACCTGAGGACTCACACAGGAGAGAAACCCTATAAATGTATGGAATGTGGAAAAAGTTACACACGAAGC
TCACATCTTGCCAGGCACAAAAGGTTCAACAAGATGAACGCGCCTTACAAATATCCCCTAAACCGGAAG
AATTTGGAAGAGACCTCCCTGTGACACAGGCTGAGAGAACTCCATCAGTGGAGAAAACCTATAGATGT
GATGATTGCGGAAAGCACTTCCGCTGGACTTCCAGACTTGTGAGACATCAGAGGACACATACTGGAGAA
AAACCCTTCTTTGACTATTTGTGGCAAAAGCTTCCAGCCAGAAAATCTGTGTTAACAACACACCAAAGA
ATCCACCTGGGAGGCAAACCTACTTGTGTGGAGAGTGTGGTGGAGACTTCAAGTGAACACAGGCGGTAC
CTGGCGCACCGGAAGACGCACGCTGCTGAGGAACTCTACCTCTGCAGCGAGTGCGGGCGCTGCTTACC
CACAGCGCAGCGTTCGCCAAGCACTTGGAGGACACGCCTCAGTGAAGCCCTGCCGATGCAACGAATGT
GGGAAGAGCTTCAAGTGCAGGGACCACCTCGTCAGGCATCAGAGAACACACACTGGGAGAAACCATTC
ACGTGCCCTACCTGTGAAAAAGCTTCCAGAGGATATCACTTAATTAGGCATCAGAGGACCCACTCA
GAAAAGACCTCCTAG
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001301779

**Insert Size:** 1947 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001301779.1](#)

**RefSeq Size:** 4393 bp

**RefSeq ORF:** 1947 bp

**Locus ID:** 7753

**UniProt ID:** [O95125](#)

**Cytogenetics:** 11q24.1

**Protein Families:** Transcription Factors

**MW:** 74.7 kDa

**Gene Summary:**

Transcriptional repressor that binds to elements found predominantly in genes that participate in lipid metabolism. Among its targets are structural components of lipoprotein particles (apolipoproteins AIV, CIII, and E), enzymes involved in lipid processing (lipoprotein lipase, lecithin cholesteryl ester transferase), transporters involved in lipid homeostasis (ABCA1, ABCG1), and several genes involved in processes related to energy metabolism and vascular disease.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same isoform (1).