

## Product datasheet for **SC337099**

### **GGA3 (NM\_001291641) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	GGA3 (NM_001291641) Human Untagged Clone
Tag:	Tag Free
Symbol:	GGA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC337099 representing NM\_001291641.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAAGAAGTGTGGGAGGAGATTTATAACGAAGTGGGAAGTCCGCTTTTTGAATGAGTTAATCAAA
GTCGTCTCTCAAAGTACCTGGGGACAGGGTGTCTGAGAAAGTGAAGACCAAGTTATTGAGCTGCTG
TACAGCTGGACCATGGCCCTGCCAGAAGAAGCAAAGATCAAAGACGCCTACCACATGCTGAAGAGACAG
GGCATAGTGCAGTCTGACCCACCAATTCCTGTGGATAGGACGCTGATCCCCTCTCCACCCTCGTCCC
AAAAACCTGTTTTGATGATGAGGAGAAGTCCAAGCTTTAGCCAAGCTGCTGAAAAGCAAAAACCCA
GATGACCTGCAGGAGGCCAACAGCTCATCAAGTCCATGGTGAAGGAAGACGAGGCACGGATCCAGAAG
GTGACCAAGCGTCTGCACACGTTAGAGGAAGTTAAACAACAGTGAAGTCTGCTCAGTGAAGTCTGCTT
CATTACAGCCAGGAGGACTCTTCGGACGGGACAGAGAGCTGATGAAGGAGCTGTTTGTATCAGTGTGAG
AACAGAGGCGGACTTTATTTAACTCGCCAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGA
CTGCAAGCCAGTGACAACCTCTCCGGGTCACTCAACTTTACAAAACAATTATTGAAGGGCAGGTATC
AATGGCGAGGTGGCTACCTTAACCTGCCTGACTCGGAAGGAAACAGTCAAGTGAAGTGAAGTGAAGTGA
CTCATCGACCTTGCGGAGCTGGACACGACCAACAGTTTGTCTCCGTGTTGGCCCCAGCCTACTCCA
CCCTCCTCAGGCATCCCAATCCTCCCTCCACCACCCAGGCCTCAGGACCTCCACGGAGCCGCTCCTCT
AGCCAGGCCGAGGCCACCTGGGGCCAGCAGCACAAAGCAACGCCTCTCCTGGCTGGACGAGGAGCTA
CTCTGCTTGGGCTCGCCGACCCAGCCCTAATGTTCTCCCAAAGAGTCAAGTGAAGTGAAGTGAAGTGA
CACCTGCTCCAGAGGGAACAGTCCGACCTGGACTTCTTACGCCCCAGGCCGGGACCGCTGCCTGTGGC
GCCTCCGATGCTCCTCTGCTCCAGCCCTCAGCCCCCTCCTCAAGCAGTCCCAAGCTCCACTGCCGCT
CCCTTCCAGCTCCTGTGGTCCCAGCCAGTGTTCCTGCCCCAGTGGGGCTCCTCCTTTGTTTTCTACT
GGAGTGGCCCCAGCCTTGCCCCAAAAGTTGAGCCCGCAGTCCCTGGGCACCATGGCTTGCGTTGGGC
AACAGCGCGCTGCACCACCTGGATGCCCTCGATCAGTCTTAGAAGAGGCCAAAGTGACCTCGGGCTTG
GTGAAACCCACTACCTCCCTCTCATCCCCACCACCACCCAGCCAGGCCCTCCTGCCCTTCTCCACG
GGGCCCGCAGCCGCTCTTCCAGCCACTGAGTTTCCAGTCCCAGGGCAGCCCCCGAAGGGGCTGAG
CTCTCCCTGGCCAGCATCCACGTGCCCTGGAATCGATCAAGCCTAGCAGTGCCTTCTGTGACAGCC
TACGATAAAAACGGCTTCCGCATCCTCTTCCACTTTGCCAAGGAGTGTCCCCAGGACGACCTGACGTG
CTGGTGGTGGTGTCCATGCTGAACACGGCTCCCTTACCTGTCAAGAGCATCGTGTGACAGGCTGCA
GTGCCAAGTCAATGAAAGTGAAGTTGCAGCCACCTTCTGGGACAGAAGTCTCTCCATTTAGCCCCATC
CAGCCACCTGCAGCCATACCCAGGTATGTTGCTGGCCAATCCACTGAAGGAGAAGGTGCGGCTTCGG
TATAAGCTGACCTTCGCCCTGGGGGAGCAGCTGAGCACAGAGGTGGGCGAGGTGGACCAGTTCCTCCT
GTGGAACAGTGGGGGAACCTATGA
ACCGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

## Plasmid Map:



ACCN: NM\_001291641

Insert Size: 1956 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001291641.1](#)

RefSeq Size: 4118 bp

RefSeq ORF: 1956 bp

Locus ID: 23163

UniProt ID: [Q9NZ52](#)

Cytogenetics: 17q25.1

Protein Families: Druggable Genome

Protein Pathways: Lysosome

**MW:** 70.1 kDa

**Gene Summary:**

This gene encodes a member of the Golgi-localized, gamma adaptin ear-containing, ARF-binding (GGA) family. This family includes ubiquitous coat proteins that regulate the trafficking of proteins between the trans-Golgi network and the lysosome. These proteins share an amino-terminal VHS domain which mediates sorting of the mannose 6-phosphate receptors at the trans-Golgi network. They also contain a carboxy-terminal region with homology to the ear domain of gamma-adaptins. Multiple alternatively spliced transcript variants have been identified in this gene. [provided by RefSeq, Feb 2010]

Transcript Variant: This variant (3) contains an additional exon in the 5' region, and it thus differs in its 5' UTR and initiates translation at a downstream in-frame start codon, compared to variant long. The encoded isoform (1) is shorter at the N-terminus, compared to isoform long. Both variants 1 and 3 encode isoform 1.