

Product datasheet for **SC337100**

SLC26A6 (NM_001281733) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC26A6 (NM_001281733) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC26A6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337100 representing NM_001281733.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGCTGGCGGATGCGTCGGGACCGAGGGACACACAGGCACTGCTGTCTGCAACACAAGCAATGGAC
CTGCGGAGGCGAGACTACCACATGGAACGGCCGCTGCTGAACCAGGAGCATTGGAGGAGCTGGGGCGC
TGGGGCTCAGCACCTAGGACCCACCAAGTGGCGGACCTGGTTGCAGTGCTCCCGTGTCTGGGCTATGCC
CTTCTGTCTCAACACCTCCCGGTTTTGGTCTGGTTACCCCGGTATCCTGTGCGTGACTGGCTCCTGGGT
GACCTGTTATCCGGCCTGAGTGTGGCCATCATGCAGCTTCCGACAGGGCTTGGCCTACGCCCTCCTGGCT
GGATTGCCCCCGTGTGGCCTCTATAGCTCCTTCTACCCTGTCTTCTACTTCTCCTGTTGGCACT
TCCCGGCACATCTCCGTGGTGTGGAGGCTGCTGGAAGCTGCCCCAGAGCAAGTTGGCACCCTGGTC
ACTGCAGCTGTGGCTGGGGTGGTCTCGTGGTGGTGAAGCTGTTGAATGACAAGCTGCAGCAGCAGCTG
CCCATGCCGATACCCGGGAGCTGCTCAGCTCATCGGGCCACAGGCATCTCCTATGGCATGGGTCTA
AAGCACAGATTTGAGGTAGATGTCGTGGGCAACATCCCTGCAGGGCTGGTCCCCCAGTGGCCCCAAC
ACCCAGCTGTTCTCAAAGCTCGTGGGACAGCCCTTACCATCGCTGTGGTTGGGTTTGGCATTGGCATC
TCACTGGGAAGATCTTCCGCCCTGAGGCACGGCTACCGGGTGGACAGCAACCAGGAGCTGGTGGCCCTG
GGCCTCAGTAACCTTATCGGAGGCATCTTCCAGTGTCTCCCGTGTGCTTATGTCTCGGAGCCTG
GTACAGGAGAGCACCGGGGCAACTCGCAGGTTGCTGGAGCCATCTTCCCTTTTCATCCTCCTCATC
ATTGTCAAACCTGGGGAACCTTCCATGACCTGCCAAGGCGGTCTGGCAGCCATCATATTGTGAAC
CTGAAGGGCATGCTGAGGAGCTCAGCGACATGCGCTCCCTCTGGAAGGCCAATCGGGCGGATCTGCTT
ATCTGGCTGGTGACCTTACGGCCACCATCTTGTGAACCTGGACCTGGCTTGGTGGTGGTGGTGGTGGT
TTCTCCCTGTGCTCGTGGTGGTCCGGACACAGATGCCCACTACTCTGTCTGGGGCAGGTGCCAGAC
ACGGATATTTACAGAGATGTGGCAGAGTACTCAGAGGCCAAGGAAGTCCGGGGGTGAAGTCTTCCGC
TCTCGGCCACCGTGTACTTTGCCAATGCTGAGTTCTACAGTGATGCGCTGAAGCAGAGGTGTGGTGTG
GATGTCGACTTCTCATCTCCAGAAGAAGAACTGCTCAAGAAGCAGGAGCAGCTGAAGCTGAAGCAA
CTGCAGAAAGAGGAGAAGCTTCCGAAACAGGCTGCCTCCCCAAGGGCGCCTCAGTTTCCATTAATGTC
AACACCAGCCTTGAAGACATGAGGAGCAACAACGTTGAGGACTGCAAGATGATGGTGGAGCTCAGGAGAT
AAGATGGAAGATGCAACAGCCAATGGTCAAGAAGACTCCAAGGCCCCAGATGGGTCCACACTGAAGGCC
CTGGCCCTGCCTCAGCCAGACTTCCACAGCCTCATCCTGGACCTGGGTGCCCTCTCCTTTGTGGACT
GTGTGCCTCAAGAGCCTGAAGAATATTTCCATGACTTCCGGGAGATTGAGGTGGAGGTGTACATGGCG
GCCTGCCACAGCCCTGTGGTCAAGCAGCTTGGAGCTGGGCACTTCTTCGATGCATCCATACCAAGAAG
CATCTCTTTGCCCTGTCCATGATGCTGTACCTTTGCCCTCCAACACCCGAGGCCTGTCCCGACAGC
CCTGTTTCGGTACCAGACTTGA
```

Restriction Sites: Sgfl-Mlul

```
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Plasmid Map:


ACCN: NM_001281733

Insert Size: 1956 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001281733.1](#)

RefSeq Size: 2338 bp

RefSeq ORF: 1956 bp

Locus ID: 65010

UniProt ID: [Q9BXS9](#)

Cytogenetics: 3p21.31

Protein Families: Transmembrane

MW: 71.5 kDa

Gene Summary:

This gene belongs to the solute carrier 26 family, whose members encode anion transporter proteins. This particular family member encodes a protein involved in transporting chloride, oxalate, sulfate and bicarbonate. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Aug 2013]

Transcript Variant: This variant (6) has multiple difference in the coding region compared to variant 1. The resulting protein (isoform 6) is shorter but has the same N- and C-termini compared to isoform 1.