

Product datasheet for **SC337107**

ANUBL1 (ZFAND4) (NM_001282905) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANUBL1 (ZFAND4) (NM_001282905) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANUBL1
Synonyms:	ANUBL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337107 representing NM_001282905.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAACCTGAAAATGATTATTGCTTGAATGATTACAACATTTTCAAGGGTGTACCTTGAAGCTAGTT
TTGGCTATGCGTGGCGGACCTATAAATACTAGAAGAGTTCTACAGACGATCCACTTAGGAAGATGGCA
GAGTACTTGGATTCCAGTCGAGTTGAGGTCTGGGAGAAGACCTCCTGCAGCAAACAAGTTACATTTTTG
GTATACCAAGAAGGAGATCAATTGAATTTCTTCTGCAGTAGATAGGGGAGATGGCACTTTAACACCG
TTATCTGACTCTTCAAAGAAAATAGATTTTCATTTGCATGTCCTACGGAGAAAAGGAGAACATCGTATG
AGTGGTGGTTCTATGTATAATTCAGATACAGATGAGGATGAAGAACTGAGCCTTCTTCTCTGGGCAA
CAGATAATTGAAAATCAATAACTATGAATAAGATGAAGCTGCTGAAGGCTAAGATGAAGAACATGAAT
CTCAGCAAAAAGCCTAAGAAAGCTGTCAAGATAAAACCTACCCACCTGTAGCTCCTCGACCTTCTAGT
GGTTCCTGCACCATCTCGCCACCGATTGTTAAGGGTCTCCCAACATTGGTCAATCTTGTTACCT
GCTTTTGGGAATGCATATCCGCCGAAAATCTCCAGGAATGGAATTTCAAGTTTGGCAACTCAACTCTCT
GCTGAGAGATACATATCTTCTATCACTGGTGAATTTCTAAGGAAGATAATAGCTGGGAGAATAACACA
CTGTCTCACTTCACTAGCAACGTCAAACTACCTCCTCAGATACCCCATTTGGAGTTGGGAAATGACCAG
GAGCTTGCTGACTCTGTTCTCCATCTTGGATCATCCCTGCCTAGGCAAACAAAACATTTTTTAGGAAAC
TTGCCATCTAGTAATGGGAACATTGCTTACCTCAGAAGAAATGTGTAACCGAACAACTACTTCTACCT
AAAGTGGGCTCACTGGCTTCATTTGAGAAGGAAATGCTGATGAACAGAGCAGCGGCTTAGAAGGTGGC
TGCAAAGTGAATCTGGAGTTGCTACTACTAATGCTGACAAAGGGTTGAAAGCTCCAGAGCAGCATCTC
AAGCATGTTGCAGGAGTCTGAATGGGAATCAGTAGAGACTTCAGTTCTTAACACCGGAATTAAGT
CCTCACAAGAATAGACTCTTGTCACTCTTCGCTGTTCTGCACCAATGTCGCTACATAATCTCTGGTG
AAACCAGAGAGACAGTCCAAATGTTTTGAGTTTGGGAAGCTACAGCCTTCTTCTTCTCAGTCACTGGAT
GTTCAAAACATAACTGATTCTTCTTCTTAGGACTACTGCTTCAAGGTGTTAAAGTTGATTCATT
GGGAAAAGATCTGATGTTATTTCCAAAGTTGAGGCTCGGGATATCACAGAAATGACTAACAAGGCTTCC
AAAGAGCTGTTGGTTGTGTAATAATATCAGTTTTTCTGCTCACTGGCCGGGAGCACAAGCAGAAAT
AGATTACAGAGCACACGTGGGCGAGGAGGCTTCAAGACTCTGGAAGTGGGCTGTCTACAAACCTCCAG
CATTTTCAGGAAGAAAACCTTAGGAAAAGTTCTCCCAAGTTAGAACATACAGGAGTTTTTTGTCTACC
CATGGTGTGGAATGAATGGAAAATATGCAGCAGCAGGAAAAGCGTAGGAGAATGTACTACTCATCAC
CTCCACCTGTGAAAGCCCTCTTCCAGCAAAAGAAAACAACAAATCATTGTTTTCTTTGTGGAAG
AAAACAGGACTGGCTAGTAGCTACGAATGCAGATGTGGAACAACCTTCTGTGCATCTCATGTTATGCA
GAAACTCATGGCTGTACCTATGATTACAAGAGTGCAGGGAGGAGATACTTGCACGAGGCAAACTCTGTG
GTTAATGCACCAAAGCTTCCAAAATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001282905

Insert Size: 1962 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282905.1](#)

RefSeq Size: 3821 bp

RefSeq ORF: 1962 bp

Locus ID: 93550

Cytogenetics: 10q11.22

Protein Families: Druggable Genome

MW: 71.5 kDa