

Product datasheet for **SC337108**

ANUBL1 (ZFAND4) (NM_001282906) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANUBL1 (ZFAND4) (NM_001282906) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANUBL1
Synonyms:	ANUBL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337108 representing NM_001282906.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAACCTGAAAATGATTATTGCTTGAATGATTACAACATTTTCAAGGGTGTACCTTGAAGCTAGTT
TTGGCTATGCGTGGCGGACCTATAAATACTAGAAGAGTTCTACAGACGATCCACTTAGGAAGATGGCA
GAGTACTTGGATTCCAGTCGAGTTGAGGTCTGGGAGAAGACCTCCTGCAGCAAACAAGTTACATTTTTG
GTATACCAAGAAGGAGATCAATTGAATTTCTTTCTGCAGTAGATAGGGGAGATGGCACTTTAACACCG
TTATCTGACTCTTCAAAGAAAATAGATTTTCATTTGCATGTCCTACGGAGAAAAGGAGAACATCGTATG
AGTGGTGGTTCTATGTATAATTCAGATACAGATGAGGATGAAGAACTGAGCCTTCTTCTTCTGGGCAA
CAGATAATTGAAAATCAATAACTATGAATAAGATGAAGCTGCTGAAGGCTAAGATGAAGAACATGAAT
CTCAGCAAAAAGCCTAAGAAAGCTGTCAAGATAAAACCTACCCACCTGTAGCTCCTCGACCTTCTAGT
GGTTCCACTGCACCATCTGCCACCGATTGTTAAGGGTCTCCCCAACATTGGTCAATCTTGTTACCT
GCTTTTGGGAATGCATATCCGCCGAAAATCTCCAGGAATGGAATTTCAAGTTTGGCAACTCAACTCTCT
GCTGAGAGATACATATCTTCTATCACTGGTGAATTTCTAAGGAAGATAATAGCTGGGAGAATAACACA
CTGTCTCACTTCACTAGCAACGTCAAACTACCTCCTCAGATACCCCATTTGGAGTTGGGAAATGACCAG
GAGCTTGCTGACTCTGTTCTCCATCTTGGATCATCCCTGCCTAGGCAAACAAAACATTTTTTAGGAAAC
TTGCCATCTAGTAATGGGAACATTGCTTACCTTCAGAAGAAATGTGTAACCGAACAACTACTTCTACCT
AAAGTGGGCTCACTGGCTTCATTTGAGAAGGAAATGCTGATGAACAGAGCAGCGGCTTAGAAGGTGGC
TGCAAAGTGAATCTGGAGTTGCTACTACTAATGCTGACAAAGGGTTGAAAGCTCCAGAGCAGCATCTC
AAGCATGTTGCAGGAGTCTGAATGGGAATCAGTAGAGACTTCAGTTCTTAACACCGGAATTAAGT
CCTCACAAAGAATAGACTCTTGTCACTCTTCGCTGTTCTGCACCAATGTCGCTACATAATCTCTGGTG
AAACCAGAGAGACAGTCCAAATGTTTTGAGTTTGGGAAGCTACAGCCTTCTTCTTCTCAGTCACTGGAT
GTTCAAAACATAACTGATTCTTCTTCTTAGGACTACTGCTTCAAGGTGTTAAAGTTGATTCATT
GGGAAAAGATCTGATGTTATTTCCAAAGTTGAGGCTCGGGATATCACAGAAATGACTAACAAGGCTTCC
AAAGAGCTGTTGGTTGTGTAATAATATCAGTTTTTCTGCTCACTGGCCGGGAGCACAAGCAGAAAT
AGATTACAGAGCACACGTGGGCGAGGAGGCTTCAGAATCTGGAAGTGGGCTGTCTACAAACCTCCAG
CATTTTCAGGAAGAAAACCTTAGGAAAAGTTCTCCCCAGTTAGAACATACAGGAGTTTTTTTGTCTACC
CATGGTGTGGAATGAATGGAAAATATGCAGCAGCAGGAAAAGCGTAGGAGAATGTACTACTCATCAC
CTCCACCTGTGAAAGCCCTCTTCAGACAAAGAAGAAAACAACAAATCATTGTTTTCTTTGTGAAAG
AAAACAGGACTGGCTAGTAGCTACGAATGCAGATGTGGAACAACCTTCTGTGCATCTCATGTTATGCA
GAAACTCATGGCTGTACCTATGATTACAAGAGTGCAAGGAGGAGATACTTGCACGAGGCAAACTCTGTG
GTTAATGCACCAAAGCTTCCAAAATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001282906

Insert Size: 1962 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282906.1](#)

RefSeq Size: 3473 bp

RefSeq ORF: 1962 bp

Locus ID: 93550

UniProt ID: [Q86XD8](#)

Cytogenetics: 10q11.22

Protein Families: Druggable Genome

MW: 71.5 kDa