

Product datasheet for **SC337110**

TTC14 (NM_001288582) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC14 (NM_001288582) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC14
Synonyms:	DRDL5813; PRO19630
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337110 representing NM_001288582.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGACCGGGACCTTTTTCGGCAGTCGCTAAATTGCCACGGGTCGTCTTTGCTCTCTCTACTTCGGAGC
GAACAGCAGGACAATCCACACTTCGGTAGCCTCCTGGGTCGGCCGCCGAGCCAGCCCGGGGCCCGCCG
CCCCAGCACCCGTTGCAGGGCAGAAAAGAGAAGAGAGTTGACAACATCGAGATACAGAAATTCATCTCC
AAAAAAGCGGATCTGCTTTTTCGACTTTCTGGAATCAGATGCACCTGCAACTTCTGAAATTAATGAA
GACAGTGAAGATCATTATGCAATCATGCCACCTTTAGAGCAATTCATGGAGATACCTAGTATGGATCGG
AGAGAGCTGTTTTCCGAGATATTGAGCGTGGTGATATAGTGATTGGAAGAATTAGTTCTATTGGGAA
TTCGGTTTTTTCATGGTGTGATCTGTTTAGGAAGTGGTATCATGAGAGATATAGCCCACTTAGAAATC
ACAGCTCTTTGTCCTTAAGAGATGTGCCTTCTCACAGTAACCATGGGGATCCTTATCATATTACCA
ACTGGTGACATCATTCGAGCTGGAATCAAGGATATTGACAGATACCATGAAAAGCTAGCAGTATCTCTG
TATAGCTCTTCTCTCCACCACCTATCTGGTATTAATAGGTGAATTAGCTCTGAAGAGCTTCTCT
TTATACTACAGGAGAAGTGTGAGCTAAATAGCAATTCTTGGAGTCCTATGAAAATGTCATGCAGAGT
TCCTTGGGATTTGTTAATCCAGGAGTAGTTGAATTCCTCTAGAAAACTAGGAATAGATGAATCTAAT
CCACCATCTTTAATGAGAGGCCTACAAAGCAAAAATTTCTCTGAAGATGATTTTGCTTCTGCATTGAGA
AAAAACAATCCGCATCTTGGGCTTTAAATGTGTGAAGATCGGAGTTGACTATTTTAAAGTTGGACGC
CATGTGGATGCTATGAATGAATACAATAAAGCTTTGGAATAGACAAACAAACGTGGAAGCTTTGGTA
GCTCGTGGAGCATTATATGCGACAAAAGGAAGTTTGAACAAAGCAATAGAAGATTTTGAGCTTGCAATTA
GAAAACGTCCAACCTCACAGAAATGCAAGAAAATACCTCTGCCAGACACTTGTAGAGAGAGGAGGACAG
TTAGAAGAAGAAGAAAAGTTTTTAAATGCTGAAAGTTACTATAAGAAAGCTTTGGCTTTGGATGAGACT
TTTAAAGATGCAGAGGATGCTTTGCAGAACTTCATAAATATATGCAGAAATCTTTGGAATTAAGAGAA
AAACAAGCTGAAAAGGAAGAAAAGCAGAAAACAAAGAAAATAGAAACAAGTGCAGAAAAGTTGCGTAAG
CTCTTAAAGAAGAGAAGAGGCTAAAGAAGAAAAGAAAATCAACTTCTTCTCAAGTGTCTTCTCT
GCTGATGAATCAGTGTCTTCATCATCATCTCTTCTCTCTGTCACAAAAGGCATAAGAAACATAAG
AGGAACCGTTTCAGAGTCTTCTCGCAGTTCCAGAAGGCATTCATCTAGGGCATCCTCAAATCAGATAGAT
CAGAATAGGAAAGATGAGTGCTACCCAGTTCCAGCTAATACTTCAGCATCTTTTCTTAACCATAACAA
GAAGTGGAGAACTACTGGGAAGCAGGATAGGTTACAGTATGAAAAGACACAGATAAAAGAGAAGAT
AGATGCCCTCTCTTTCATCTTCACTTGAAATACCGGATGATTTTGGAGTGTACTCTATTTATTTAAA
AAGTTAACTATAAAACAGCCTCAGGCAGGTCCTTCAGGAGATATTCCAGAAGAGGGCATTGTATCATA
GATGACAGCTCCATTCATGTTACTGACCCTGAAGACCTTCAAGTGGGACAAGATATGGAGGTGGAAGAC
AGTGGTATTGATGATCCTGACCACGGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001288582

Insert Size: 1962 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288582.1](#)

RefSeq Size: 2881 bp

RefSeq ORF: 1962 bp

Locus ID: 151613

UniProt ID: [Q96N46](#)

Cytogenetics: 3q26.33

MW: 73.8 kDa