

## Product datasheet for **SC337114**

### ZNF133 (NM\_001282999) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF133 (NM_001282999) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF133
Synonyms:	pHZ-13; pHZ-66; ZNF150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337114 representing NM\_001282999.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCATT CAGGGATGTGGCTGTGGATTTACCCAGGATGAGTGGAGGCTGCTGAGCCCTGCTCAAAGG
ACTCTGTACAGAGAGGTGATGCTGGAGAACTACAGCAACCTGGTCTCACTGGGAATTTCAATTTCTAAA
CCAGAACTCATCACCCAGCTGGAGCAAGGAAAGAGACCTGGAGAGAGGAAAAAATGTTACCCGGCA
ACCTGTCCAGCAGATCCAGAGCCAGAGCTCTACCTCGATCCTTTCTGCCCTCCGGGTTTCTCCAGTCAG
AAATTCCCCATGCAGCATGTGCTGTGTAATCATCCCCCTGGATCTTACATGCTTGTGTGCAGAAGGT
AACATCCAGCCTGGGGATCCGGGCCAGGGGACCAGGAGAAGCAGCAACAAGCCTCTGAGGGGAGACCC
TGGAGTGATCAAGCAGAAGTCTGAGGGAGAAGGTGCCATGCCTTTGTTTGAAGAACCAAGAAAAGG
ACTCTGGGAGCGTTCTCCAGGCCACCCAGAGGCAGCCAGTCAGCTCTCGGAACGGCCTCAGAGGGGTG
GAGTTAGAAGCCAGCCCAGCTCAGACAGGGAACCTGAGGAAACAGACAAATTGTTGAAGAGGATAGAA
GTCTTAGGATTTGGAACAGTCAACTGTGGAGAGTGTGGACTGAGCTTCAGCAAGATGACAAACCTGCTC
AGTCACCAGCGGATACACTCAGGGGAGAAGCCCTACGTGTGTGGGTATGTGAGAAGGGGTTACAGCCTA
AAGAAGAGCCTCGCCAGACACCAGAAGGCACACTCGGGGAGAAGCCAATTGTGTGCAGGGAGTGTGGA
CGAGGCTTTAACCGAAGTCAACGCTAATCATACACGAACGGACACACTCCGGTGAGAAACCTTACATG
TGCAGTGAGTGTGGCGAGGCTTACGCCAGAAGTCAAACCTCATCATACACCAGAGGACACACTCAGGG
GAAAAGCCTTACGTGTGCCGGGAATGTGGCAAAGGCTTACGCCAGAAGTCACTGTGCGTGCAGACACCAG
AGGACACACTTGGAGGAGAAGACCATCGTGTGCAGTACTGTGGCCTGGGCTTACGCGACAGGTCAAAC
CTCATCTCCACCAGAGGACGCACTCTGGGGAGAAGCCCTACGCCTGCAAGGAGTGTGGCGATGCTTC
AGGCAGAGGACCACCCTTGTCAACCACCAGAGGACACACTCAAAGGAGAAGCCCTATGTGTGCGGGGTG
TGTGGGCACAGCTTCAGCCAGAATTCAACCCTCATCTCTCACAGGCGGACACACTGGGGAGAAGCCG
TATGTTTGTGGGTGTGTGGGCGAGGCTTTAGTCTCAAGTCACACCTCAACAGACACCAGAACATACAC
TCAGGAGAGAAGCCATTGTGTGCAAGGACTGTGGCCGGGCTTACGCCAGCAATCCAACCTCATCAGA
CACCAGAGGACGCACTCAGGCGAGAAGCCATGGTGTGTGGGAGTGCAGGCGAGGCTTACGCCAGAAG
TCAAACCTTGTGCACACCAGAGGACGCACTCAGGGGAGAGGCGTATGTGTGCCGAGAGTGCAGGCGA
GGCTTTAGCCACCAGGCCGTCTCATCAGGCACAAGCGGAAGCACTCGAGGGAGAAGCCCTACATGTGC
AGGCAGTGTGACTGGGCTTTGGCAATAAGTCAGCTTAATTACACACAAGCGGGCTCACTCGGAAGAG
AAGCCTTGTGTGTGCAGAGAGTGTGCCAAGGCTTCTCCAAAAGTCACACCTCACCTTACATCAAATG
ACACATACGGGGAGAAGCCATATGTGTGCAAGACGTGTGGGCGGGCTTACGCCTCAAGTCTCACCTC
AGCAGACACAGGAAGACCAGTCTGTCCACCACAGACTGCCAGTGCAGCCGACCCTGAGCCGTGTGCA
GGGCAACCTTCGGATTCCTTATACTCTCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

## Plasmid Map:



ACCN: NM\_001282999

Insert Size: 1965 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001282999.1](#)

RefSeq Size: 2503 bp

RefSeq ORF: 1965 bp

Locus ID: 7692

UniProt ID: [P52736](#)

Cytogenetics: 20p11.23

Protein Families: Transcription Factors

MW: 73.4 kDa

**Gene Summary:**

May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (7) differs in the 5' UTR and uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. The encoded isoform (c) is 1 aa longer than isoform a. Variants 5, 6, 7, 8 and 9 all encode isoform c.