

Product datasheet for **SC337123**

EXOC7 (NM_001282313) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EXOC7 (NM_001282313) Human Untagged Clone
Tag:	Tag Free
Symbol:	EXOC7
Synonyms:	2-5-3p; BLOM4; EX070; EXO70; Exo70p; EXOC1; NEDSEBA; YJL085W
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337123 representing NM_001282313.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGTCTATCTTATCATCCTTTGAGAGCCGCCTTATGAAGCTGGAGAAGCTCCATCATCCCTGTGCAC
AAGCAGACGGAGAATCTGCAGCGGTGCAGGAGAATGTTGAGAAGACGCTGTCTGCCTGGACCATGTC
ATCAGCTACTACCATGTGGCCAGTGACACTGAGAAGATCATCAGAGAGGGCCCCACAGGTAGGCTGGAA
GAGTACCTGGGAAGCATGGCCAAGATTCAGAAGGCAGTGGAGTATTTCCAGGACAACAGCCAGACAGC
CCGGAAGCTCAACAAAGTAAACTGCTCTTTGAGCGCGGAAGGAGGCCCTGGAGTCCGAATTTGCGAGC
CTGATGACGCGGCACAGTAAGGTCGTCTCGCCGTGCTCATCTTGGATCTGATCAGTGGTACGATGAT
CTGGAGGCCAGGAGGACGTGACCTGGAGCACCTGCCGAGAGCGTGTCCAGGATGTCATTGCGATC
TCCCGTGGCTGGTGAATATGGCCGAACCAAGATTTTCATGAACGTCTACTACCAGATACGCTCCAGC
CAGCTGGACCGCTCCATCAAAGGACTGAAGGAGCATTTCATAAGAGCAGTTCTTCTCTGGGGTCCC
TACTCCCTGCTATCCCCAACAAGAGGAAAGACACACCTACCAAGAAGCCAGTCAAGCGGCCAGGGACG
ATCCGTAAAGGCTCAGAACCTTCTGAAACAGTATCCCAGCATGGTCTAGATGGGAAAAAGGGGGGCTCT
AACCTCATTCTCTGGAAGGGAGAGATGACATGCTGGACGTGGAGACCGATGCCTACATCCACTGCGTC
AGTGCCTTCGTAAGCTGGCGCAGAGCGAGTACCAGCTGCTGGCCGACATCATCCCCGAGCACCACCAG
AAGAAGACCTTCGACTCCCTGATACAGGATGCCCTGGATGGGCTGATGCTGAAGGGGAGAACATCGTG
TCTGCTGCCCGGAAGGCCATTGTGCGACACGACTTCTCCACGGTGTCTCACCGTCTCCCCATCTGCGA
CACCTCAAGCAGACCAAGCCTGAGTTTGACCAGGTGCTCCAGGGCACGGTGCAGCACAAGAACAAG
CTGCCTGGCCTCATCACATCCATGGAGACCATCGGTGCCAAAGCGCTGGAGGACTTCGCAGACAACATC
AAGAATGACCCGGACAAGGAGTACAACATGCCGAAGGACGGCACCGTACACGAGCTCACCAGCAATGCC
ATCCTCTTCTGACGAGCTTTTGGACTTCCAGGAGACGGCAGGCCATGCTGGCCTCCCAAGTTCTT
GGGACACATAACAATATTCCTTTAGACCCCGAGAGACCAGCTCTTCGGCCACCAGCTACAGCTCTGAG
TTCAGCAAGCGGCTGCTAAGCACCTATATCTGTAAGTGTGGGCAACCTGCAGTTGAAGTGTGAGC
AAGTCCAAGGTGTACGAGGACCCAGCTCTGAGCGCCATCTTCTGCACAACAACATAAATTACATCCTC
AAGTCCCTGGAGAAGTCTGAAGTATCCAGCTGGTGGCAGTGACACAGAAGACTGCTGAGCGCTCCTAC
CGGAGCACATTGAGCAGCAGATCCAGACCTACCAGCGCAGCTGGTTAAAGGTGACTGATTACATCGCA
GAGAAGAATCTACCTGTGTTCCAGCCGGGAGTCAAGCTCCGGGACAAGGAGCGGCAGATTATCAAGGAG
CGTTTTAAGGGCTTCAATGATGGCTCGAAGAAGTGTGCAAAATCCAGAAGGCCTGGGCTATTCCAGAC
ACAGAGCAGAGGGACAGGATTCGCCAGGCCAGAAAGACATTGTCAAGGAGACCTACGGGGCCTTTCTA
CAGAAGTTTGGCAGCGTGCCCTTACCAAGAACCCTGGAGAAGTACATCAAGTACGGGGTGGAGCAGGTG
GGCGACATGATCGATCGCCTTTTCGACACCTCTGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001282313

Insert Size: 1971 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282313.1](#)

RefSeq Size: 4834 bp

RefSeq ORF: 1971 bp

Locus ID: 23265

UniProt ID: [Q9UPT5](#)

Cytogenetics: 17q25.1

Protein Families: Druggable Genome

Protein Pathways: Insulin signaling pathway

MW: 74.7 kDa

Gene Summary: The protein encoded by this gene is a component of the exocyst complex. The exocyst complex plays a critical role in vesicular trafficking and the secretory pathway by targeting post-Golgi vesicles to the plasma membrane. The encoded protein is required for assembly of the exocyst complex and docking of the complex to the plasma membrane. The encoded protein may also play a role in pre-mRNA splicing through interactions with pre-mRNA-processing factor 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Nov 2011]

Transcript Variant: This variant (3) differs in the 5' UTR and has multiple coding region differences, compared to variant 1. These differences cause translation initiation at a downstream AUG, and the resulting isoform (3) has a shorter N-terminus and is shorter compared to isoform 1.