

Product datasheet for **SC337128**

SLCO6A1 (NM_001289004) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLCO6A1 (NM_001289004) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLCO6A1
Synonyms:	CT48; GST; OATP-I; OATP6A1; OATPY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337128 representing NM_001289004.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTTCGTAGGCGTCGCCCGCAGCTCTGGGAGCCAGGATGAAGTCTCAAGGGGAGTAGAGCCGCTGGAG
GCCGCGCGGGCCAGCCTGCTAAGGACAGGAGGCCAAGGGAACCCGAAGTCCTCGAAGCCCGGGAAA
AAACACCGGTATCTGAGACTACTTCCAGAGGCCTTGATAAGGTTTCGGCGGTTTCCGAAAAAGGAAAAA
GCCAAGTCCTCAGTTTCCAAGAAGCCGGGAGAAGTGGATGACAGTTTGGAGCAGCCCTGTGGTTTGGGC
TGCTTAGTCAGCACCTGCTGTGAGTGTTGCAATAACATTCGCTGCTTCATGATTTTCTACTGCATCCTG
CTCATATGTCAAGGTGTGGTGTTTGGTCTTATAGATGTCAGCATTGGCGATTTTCAGAAGGAATATCAA
CTGAAAACCATTTGAGAAGTTGGCATTGAAAAGAGTTACGATATTTTCATCTGGCCTGGTAGCAATATTT
ATAGCATTCTATGGAGACAGAAAAAAGTAATATGGTTTGTAGCTTCCTCCTTTTAAATAGGACTTGGG
TCACTTTTATGTGCTTTTCCATCCATTAATGAAGAAAAATAACAAAGTAAGGTAGGAATTGAAGGTATT
GCAGAATGTACATCAATGATTGGATATGCTCTGGGTTATGTGCTAGGAGCACCAGTAAAGTCCCT
GAGAATACTACTTCTGCAACAAACACTACAGTCAATAATGGTAGTCCAGAATGGCTATGGACTTGGTGG
ATTAATTTTCTTTTGGCGCTGTCGTTGCATGGTGACATTAATACCATTGTCATGCTTTCCAAACAAT
ATGCCAGGTTCAACACGGATAAAAGCTAGGAAACGTAACAGCTTCATTTTTTTCAGCAGCAGACTTAAA
GATCTGAAACTTGGAACTAATATCAAGGATTTATGTGCTGCTCTTTGGATTCTGATGAAGAATCCAGTG
CTCATATGCCTAGCTCTGTCAAAAGCTACAGAATTTTAGTTATTATTGGAGCTTCTGAATTTTTCCT
ATATATTTAGAAAATCAGTTTATATTAACACCCACTGTGGCAACTACACTTGCAGGACTTGTTTTAATT
CCAGGAGGTGCACTTGGCCAGCTTCTGGGAGGTGTCATTGTTCCACATTAGAAATGTCTGTAAAGCC
CTTAGAGATTATAATGGTTACATCTGTGATATCACTTATACTGCTTGTGTTTATTATTTTGTACGC
TGTAATCCAGTGCAATTTGCTGGGATCAATGAAGATTATGATGGAACAGGGAAGTTGGGAAACCTCACG
GCTCCTTGCAATGAAAAATGTAGATGCTCATCTTCAATTTATTCTTCTATATGTGAAGAGATGATATT
GAATATTTTCTCCCTGCTTTCAGGGTGACATATTCTAAAGCACAAAACCAAAAAAGATGTACTAC
AATTGTTCTTGCAATTAAGAAGGATTAATAACTGCAGATGCAGAAGGTGATTTTATTGATGCCAGACCC
GGGAAATGTGATGCAAGTGCTATAAGTTACCTTTGTTTCATTGCTTTTATCTTTTCTACACTTATATTT
TCTGGTTTTTCTGGTGACCAATCGTCTTGCCATGACGCGGGTGTACCTGACAACTGCGTTCTCTG
GCCTTGGGTGTAAGCTATGTGATTTTGAAGATTTTGGGACTATTCCTGGACCATCAATCTTTAAATG
TCAGGAGAACTTCTTGATTTTACGGGATGTTAATAAATGTGGACACACAGGACGTTGTTGGATATAT
AACAGACAAAAATGGCTTTCTTATTGGTAGGAATATGTTTTCTTGCAAACTATGCACTATCATCTTC
ACTACTATTGCATTTTTCATATACAAACGTCGTCTAAATGAGAACACTGACTTCCAGATGTAACGTG
AAGAATCCAAAAGTTAAGAAAAAGAAGAACTGACTTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001289004

Insert Size: 1974 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289004.1</u>
RefSeq Size:	2503 bp
RefSeq ORF:	1974 bp
Locus ID:	133482
UniProt ID:	<u>Q86UG4</u>
Cytogenetics:	5q21.1
Protein Families:	Transmembrane
MW:	72.6 kDa