

Product datasheet for **SC337131**

WDR93 (NM_001284395) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR93 (NM_001284395) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR93
Synonyms:	C1d-87; CFAP297; FAP297
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC337131 representing NM_001284395.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTCATTCCCAAGAGGAAGTCAGACCCAGAAAATAAAGCACCCATTGGTACACGAAAGGGACCATTG
GAGGTGCCACCCCAACAGAGAAGGACTGGCCTAAAGACGATGAACAGGATCATGTCCTCGTGGATCCA
GATGAGGAGCTGGATTCCCTTGCCTCAGCCTTATCGAATGATCAACAAGCTGGTGAACCTTCTGTTTGAC
CAGTCTTGGGAAATTATTGAAGAGAGAAACGCACTGAGGGAAGCTGAGAGCAGCCAGATCCAGCCACC
GTCTACCTCCACTTGGAGAAATCCAGCTCAACAAAATGCCAAATTGTATGGCTGTTTCCCAAGACTAT
GTGTTTATTGGAGGAGCCAAAGGATTCTCAATTTATAATCTGTACAGTGCTAAACAAATATATGCGTGG
GAGAAGCTTAAGTTGATGTCACTTCCATCTGGCCACAGATTTAGGGAATGAAATACTATTGCTCCT
GTGGATGAAATGGGTATCATTAGACTCTTTATTTTTATAAGGAAGGACTTTACCTAGTCAAAGCCATC
AATGAAGTGGATGATACCAGCAAGCAAACTACCTGTATAAAGATGGAGATCTCTCAAGGAGGGGACTTT
GCAGCCTTCTCTACAAGGAGCCGAGATATTTGGCTGGATGTGTATAAATTGCCAAGGAGACTTGG
CTCAAGAAACTAGAGCACCCCAACTCACTTCAAATCCAAAGAAAAAAGTCAGACAGCCGCAACTGAAC
TCCCTTGGTCCCATTCTGCAGATCCTTTAGAAATGGATGCAAAATGTTAGTTTTAAAGGAGACATTTAAA
TTGAGTCTTCCAGTTTACATAATGAAGATCAAACCACTAAGCCTGTTACAGGCACCACATTCAAAGC
CCCCTGGAAGTCTTTGCAAAGATCAAGGACTGCTACGACTGGGCTCCGGGCAGAATCATTTATCAAG
GACAGTCAGTGGGAGCAGCAAGCTGAGATCTTCAACGCTTCTACAAGAAGTACCTAGATAGGGAGTGG
GAGGAAGAGCCACTCAGTACGGCCACCTTCTATTTCTTCTTCTAGCTGCCTATTTGCAATGCCACCG
GAAGTCAAGGGCCCTCAGGAATGGCCTGTGCTTGGTATACACTGGACCAGAAGTCACAATTTCTTC
CTGTATTCATAAACCGAACTCTAAAGGATAAAGCTGACCCCGAGGGTGTGTGGCCCTGTGCTGCGCCC
ATTGCAGTCTCTCAGCTCAGCTGCTCCTCCTACCTGGTGTGCTGCGCTGCGAGGATGGTGTGCTCACG
CTGTGGGACCTGGCCAAAGGACAGAATATGTATCCTGAAGGTCAAGTGAATCCCAAATGAAATGTGTG
GTGCTGTGCACAGACGCCTCCCTCCATCTGGTGGAGGCTAGCGGGACCCAAGGACCCACCATCAGTGTG
CTTGTTGAGAGGCCTGTAAAGCACCTGGATAAAACCATCTGTGCCGTGGCCCCAGTCCCAGCCTTACCT
GGCATGGTGTCTATTTTTCCAAGAATGGCTCTGTGTGCCTTATGGATGTGGCAAGCGTGAAATCATC
TGTGCTTTGCCCTCCGGGAGCCTTCTCTGGAGGTCCCTGGAAGCCAGTGTGCTGTGTCTCCA
GACCATCCATGTTTCTGCTCCGAGGAGACTATTCACATGAAACTGCGTCCACCGACGATGCTGGAATC
CAATATTCGTGTTTCTATTTAATTTTGGGCCTGCCACTCCTGAAAAATATCTCAAAAAATTGTACC
ATTCTCAAAGGGACTTGGATAACATGGCCTTCCCCCAAGCACTGCCACTGGAGAAGAGATGTGAGCGT
TTCTCCAGAAGAGCTATCGGAAGCTGGAGAAGACCCAGAGAAGGAGGAGGAGCACTGGGCCCGGCTT
CAGAGGTACTCCTGTGCTCCAGAGAGAGAACTTCAAGAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001284395

Insert Size: 1977 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001284395.1](#)

RefSeq Size: 2361 bp

RefSeq ORF: 1977 bp

Locus ID: 56964

UniProt ID: [Q6P2C0](#)

Cytogenetics: 15q26.1

MW: 74.3 kDa