

Product datasheet for **SC337135**

OSBP2 (NM_001282741) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OSBP2 (NM_001282741) Human Untagged Clone
Tag:	Tag Free
Symbol:	OSBP2
Synonyms:	HLM; ORP-4; ORP4; OSBPL1; OSBPL4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337135 representing NM_001282741.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGGACGCGGCTGCTGGGACACGGCTCCGTGCGCTCTTCTGCAGCTCCGGTGCCTGAATAAA
TTTAAGAAGGATGACTCTGGGGACGACGACGAGGCTACCACCCAGCCGACAAGAGCGAGCTGCACCAC
ACCCTGAAGAATCTTCCCTGAAGTTAGATGACCTCAGCACGTGCAATGACCTCATCGCCAAGCACGGC
GCTGCACTCCAGCGCTCCCTGACAGAGCTGGACGGCCTCAAGATCCCATCTGAGAGTGGGGAGAAGCTG
AAGGTGGTGAATGAGCGGGCCACCCTCTCCGCATCACATCCAATGCTATGATCAACGCCTGCAGGGAC
TTCTTGGAACTAGCAGAGATACACAGTCGGAATGGCAGCGGGCACTGCAGTATGAGCAGGAGCAGCGC
GTGCACTTGGAGGAAACCATTGAGCAGCTGGCGAAGCAGCACAACAGCCTCGAGCGGGCCTTCCACAGT
GCCCTGGCCGGCCGCAACCCCTCCAAGAGCTTCATTGAGGGAAGCCTCTTGACTCCCAAGGAGAG
GACAGTGAGGAAGATGAAGATACCGAGTACTTTGATGCCATGGAAGACTCCACATCCTTCATCACCGTG
ATCACCGAGGCCAAGGAAGACAGCAGAAAAGCTGAAGGTAGCACCGGGACAAGTCCGTGGACTGGAGC
TCAGCAGACAATGTACTAGATGGTGCCTCGCTCGTCCCAAGGGTTCATCCAAAGTCAAGAGGCGAGTC
CGCATTCCCAACAAGCCCAACTACAGCCTTAACCTCTGGAGCATCATGAAGAACTGCATCGGCCGGGAG
CTCTCCAGGATCCCCATGCCGGTGAACCTCAATGAGCCCCTGTCCATGCTCCAGCGGCTGACAGAGGAC
CTGGAGTACCACCACCTGCTGGACAAGGCACTGCAGTGCACCAGCTCAGTGGAGCAGATGTGCCTGGTG
GCCGCTTCTCTGTGCTCCTACTCCACCACAGTGCACCGCATCGCCAAGCCCTCAACCCCATGCTG
GGGAGACCTTCGAGCTGGACCGCCTCGACGACATGGGCTGCGCTCCCTCTGTGAGCAGGTGAGCCAC
CACCCCCCTCAGCTGCGCACTACGTGTTCTCCAAGCATGGCTGGAGCCTCTGGCAGGAGATCACCATC
TCCAGCAAGTCCGGGGAAAATACATCTCCATCATGCCGCTAGGTGCCATCCACTTAGAATCCAGGCC
AGTGGGAATCACTACGTGTGGAGGAAGAGCACCTCAACTGTTCAACAACATCATCGTGGGCAAGCTCTGG
ATCGACCAGTCAGGGGACATCGAGATTGTGAACCATAAGACCAATGACCGGTGCCAGCTGAAGTTCTG
CCCTACAGCTACTTCTCCAAAGAGGCAGCCCGAAGGTGACAGGAGTGGTGAAGTACAGCCAGGGCAAG
GCCATTACGTGCTGTCCGGCTCGTGGGATGAACAAATGGAGTGTCCAAGGTCATGCATAGCAGTCCC
AGCAGCCCCAGCTCTGACGGGAAGCAGAAGACAGTGTACCAGACCCTGTCCAGCAAGCTGCTGTGGAAG
AAGTACCCGCTGCCGGAGAACGCGGAGAACATGTACTACTTCTCAGAGCTGGCCCTGACCCTCAACGAG
CACGAGGAGGGCGTAGCCCAACCGACAGCCGCTGCGGCCGACCAGCGGCTGATGGAGAAGGGCCGT
TGGGACGAGGCCAATACCGAGAAGCAGCGGCTGGAGGAGAAGCAGCGCCTGTCCGGCGCCGGCGGCTG
GAGGCCTGCGGGCCGGGACGAGCTGCAGCTCGGAGGAAGCAGAGAAGGAGGCGGATGCCTACACGCCA
CTGTGTTTGAAGAGGCTGGATCCGCTGACCGGGAGATGGCCTGTGTGTACAAGGGCGGCTACTGG
GAGGCCAAGGAGAAGCAAGACTGGCATATGTGCCCAACATCTTCTGA
ACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001282741

Insert Size: 1980 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282741.1](#)

RefSeq Size: 3589 bp

RefSeq ORF: 1980 bp

Locus ID: 23762

UniProt ID: [Q969R2](#)

Cytogenetics: 22q12.2

MW: 74.7 kDa

Gene Summary:

The protein encoded by this gene contains a pleckstrin homology (PH) domain and an oxysterol-binding region. It binds oxysterols such as 7-ketocholesterol and may inhibit their cytotoxicity. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2013]

Transcript Variant: This variant (6) differs in the 5' UTR, lacks a portion of the 5' coding region, and uses an alternate in-frame splice site, compared to variant 1. It represents use of an internal promoter and initiates translation at an alternate start codon. The encoded isoform (f) is shorter and has a distinct N-terminus, compared to isoform a.