

Product datasheet for **SC337140**

ZNF133 (NM_001283006) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF133 (NM_001283006) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF133
Synonyms:	pHZ-13; pHZ-66; ZNF150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337140 representing NM_001283006.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGAGCTGCCGGCCTTCAGGGCGCCCTGGTGGCCCGGAGCATCCTGGGAGGTGTGGTCCGGGCGCT
GTGCCGAGGATTCGGGGTAGTGTAGTCCTGGCGCCCGCTGGAGGAGCTCCCAGCTGGTGGCTGGAGC
GACCCCTTGTCTTTGGTGGCGCTGAACCTCATACCCAGCTGGAGCAAGGAAAGAGACCTGGAGAGAG
GAAAAAATGTTACCCGCAACCTGTCCAGCAGATCCAGAGCCAGAGCTCTACCTCGATCCTTTCTGC
CCTCCGGTCTCCAGTCAGAAATTCATGCAGCATGTGCTGTGTAATCATCCCCCTGGATCTTC
ACATGCTGTGTGCAGAAGTAACATCCAGCCTGGGGATCCGGGCCAGGGGACCAGGAGAAGCAGCAA
CAAGCCTCTGAGGGGAGACCTGGAGTGATCAAGCAGAAGGTCTGAGGGAGAAGGTGCCATGCCTTTG
TTTGAAGAACCAAGAAAGGACTCTGGGAGCGTTCTCCAGGCCACCCAGAGCCAGCCAGTCAAGCTCT
CGGAACGGCTCAGAGGGTGGAGTTAGAAGCCAGCCAGCTCAGACAGGGAACCCTGAGGAAACAGAC
AAATTGTTGAAGAGGATAGAAGTCTTAGGATTTGGAACAGTCAACTGTGGAGAGTGTGGACTGAGCTTC
AGCAAGATGACAAACCTGCTCAGTCACAGCGGATACACTCAGGGGAGAAGCCCTACGTGTGTGGGTA
TGTGAGAAGGGCTTCAGCCTAAAGAAGAGCCTCGCCAGACACCAGAAGGCACACTCGGGGAGAAGCCA
ATTGTGTGCAGGGAGTGTGGACGAGGCTTAAACCGAAGTCAACGCTAATCATACACGAACGGACACAC
TCCGGTGAGAAACCTTACATGTGCAGTGTGGGCGAGGCTTCAGCCAGAAGTCAAACCTCATCATA
CACCAGAGGACACACTCAGGGGAAAAGCCTTACGTGTGCCGGAAATGTGGCAAAGGCTTCAGCCAGAAG
TCAGCTGTGTCGAGACACCAGAGGACACACTTGGAGGAGAAGACCATCGTGTGCAGTGTGGCTG
GGCTTCAGCGACAGGTCAAACCTCATCTCCACCAGAGGACGCACTCTGGGGAAGCCCTACGCCTGC
AAGGAGTGTGGCGATGCTTCAGGCAGAGGACCCTTGTCAACCACCAGAGGACACACTCAAAGGAG
AAGCCCTATGTGTGCGGGTGTGTGGGCACAGCTTCAGCCAGAATTCAACCCTCATCTCTCACAGGCGG
ACACACACTGGGGAAGCCGTATGTTTGTGGGTGTGTGGGCGAGGCTTGTCTCAAGTCACACCTC
AACAGACACCAGAACATACTCAGGAGAGAAGCCATTGTGTGCAAGGACTGTGGCCGGGGCTTCAGC
CAGCAATCCAACCTCATCAGACACCAGAGGACGCACTCAGGCGAGAAGCCATGGTGTGTGGGAGTGC
GGGCGAGGCTTCAGCCAGAAGTCAAACCTTGTGTCACACCAGAGGACGCACTCAGGGGAGAGGCCGTAT
GTGTGCCGAGAGTGCAGGCGAGGCTTGTAGCCACCAGGCGGCTCATCAGGCACAAGCGGAAGCACTCG
AGGGAAGCCCTACATGTGCAGGAGTGTGGACTGGCTTTGGCAATAAGTCAAGTCTAATTACACAC
AAGCGGGCTCACTCGGAAGAGAAGCCTTGTGTGTGCAGAGAGTGTGGCCAAGGCTTCTCCAAAAGTCA
CACCTCACCTTACATCAATGACACATACGGGGGAGAAGCCATATGTGTGCAAGAGCTGTGGCGGGGC
TTCAGCCTCAAGTCTCACCTCAGCAGACACAGGAAGACCAGTCTGTCCACCACAGACTGCCAGTGCAG
CCCAGCCTGAGCCGTGTGCAGGGCAACCTTCGGATTCTTATACTCTCTCTGA
ACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001283006

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001283006.1](#)

RefSeq Size: 2342 bp

RefSeq ORF: 1986 bp

Locus ID: 7692

Cytogenetics: 20p11.23

Protein Families: Transcription Factors

MW: 73.2 kDa

Gene Summary:

May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (14) lacks four internal exons, initiates translation at an alternate start codon, and uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. The encoded isoform (g) has a distinct N-terminus and is longer than