

Product datasheet for **SC337142**

PDE4 (PDE4B) (NM_001297441) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE4 (PDE4B) (NM_001297441) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDE4B
Synonyms:	DPDE4; PDEIVB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC337142 representing NM_001297441.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGAAAAGCAGGAGTGTGATGACGGTGTGCTGATGTAATTTTAAAGAACTTTGATGTGGAA
AATGGCCCTTCCCAGGTCGGAGTCCACTGGATCCCCAGGCCAGCTCTTCCGCTGGGCTGGTACTTAC
GCCACCTTTCTGGGCACAGCCAGCGCAGAGAGTCATTTCTCTACAGATCAGACAGCGACTATGACTTG
TCACCAAAGGCGATGTCGAGAACTCTTCTCTTCCAAGCGAGCAACACGGCGATGACTTGATTGTAAC
CCTTTTGGCCAGGTCCTTGGCAGCTTGGCAAGTGTGAGAAACAACCTTCACTATACTGACAAAACCTTCA
GGTACATCTAACAAAGAGGTCCTCCAGCTGCTAGTCAGCCTCCTGTCTCCAGAGTCAACCCACAAGAAGAA
TCTTATCAAAAATTAGCAATGGAACCGTGGAGGAATTAGACTGGTGTAGACCAGCTAGAGACCATA
CAGACCTACCGGTCTGTCACTGAGATGGCTTCAACAAGTTCAAAGAATGCTGAACCGGGAGCTGACA
CACCTCTCAGAGATGAGCCGATCAGGGAACCGGTGTCTGAATACATTTCAAATACTTTCTTAGACAAG
CAGAATGATGTGGAGATCCATCTCTACCCAGAAAGACAGGGAGAAAAAGAAAAAGCAGCAGCTCATG
ACCCAGATAAGTGGAGTGAAGAAATTAATGCATAGTTCAAGCCTAAACAATACAAGCATCTCACGCTT
GGAGTCAACACTGAAAATGAAGATCACCTGGCCAAGGAGCTGGAAGACCTGAACAATGGGGTCTTAAC
ATCTTTAATGTGGCTGGATATTCTCACAATAGACCCCTAACATGCATCATGTATGCTATATTCCAGGAA
AGAGACCTCCTAAAGACATTCAGAACTCTCATCTGACACATTTATAACCTACATGACTTTAGAAGAC
CATTACCATTCTGACGTGGCATATCACACAGCCTGCACGCTGCTGATGTAGCCAGTCGACCCATGTT
CTCCTTTCTACACCAGCATTAGACGCTGTCTTACAGATTTGGAGATCCTGGCTGCCATTTTGCAGCT
GCCATCCATGACGTTGATCATCTGGAGTCTCCAATCAGTTTCTCATCAACAAAATTCAGAACTTGCT
TTGATGTATAATGATGAATCTGTGTTGGAAAATCATCACCTTGCTGTGGGTTTCAAACCTGTCAGAA
GAACACTGTGACATCTTCATGAATCTACCAAGAAGCAGCGTCAGACACTCAGGAAGATGGTTATTGAC
ATGGTGTAGCAACTGATATGTCTAAACATATGAGCCTGCTGGCAGACCTGAAGACAATGGTAGAAACG
AAGAAAGTTACAAGTTCAGGCGTCTTCTCCTAGACAACATAACCGATCGCATTAGGTCCTTCGCAAC
ATGGTACACTGTGACAGCCTGAGCAACCCACCAAGTCTTGGAAATGATCGGCAATGGACAGACCCGC
ATCATGGAGGAATTTTCCAGCAGGGAGACAAAGAGCGGGAGAGGGGAATGGAAATTAGCCCAATGTGT
GATAAACACACAGCTTCTGTGAAAAATCCCAGGTTGGTTTCATCGACTACATTGTCATCCATTGTGG
GAGACATGGGCAGATTTGGTACAGCCTGATGCTCAGGACATTCTCGATACCTTAGAAGATAACAGGAAC
TGGTATCAGAGCATGATACCTCAAAGTCCCTCACCACCCTGGACGAGCAGAACAGGGACTGCCAGGGT
CTGATGGAGAAGTTTCAGTTTGAAGTACTCTCGATGAGGAAGATTCTGAAGACCTGAGAAGGAGGGA
GAGGGACACAGCTATTTACAGCAGCAAAAGACGCTTTGTGTGATTGATCCAGAAAACAGAGATCCCTG
GGAGAGACTGACATAGACATTGCAACAGAAGACAAGTCCCCCGTGGATACATAA
ACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001297441

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297441.1](#)

RefSeq Size: 3963 bp

RefSeq ORF: 1986 bp

Locus ID: 5142

UniProt ID: [Q07343](#)

Cytogenetics: 1p31.3

Protein Families: Druggable Genome

Protein Pathways: Progesterone-mediated oocyte maturation, Purine metabolism

MW: 75.2 kDa

Gene Summary: This gene is a member of the type IV, cyclic AMP (cAMP)-specific, cyclic nucleotide phosphodiesterase (PDE) family. The encoded protein regulates the cellular concentrations of cyclic nucleotides and thereby play a role in signal transduction. Altered activity of this protein has been associated with schizophrenia and bipolar affective disorder. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014]
Transcript Variant: This variant (f) differs in the 5' UTR and contains multiple differences in the 5' coding region, compared to variant 1. The encoded isoform (5) is shorter than isoform 1.